



Small-RNA sequencing libraries with greatly reduced adaptor-dimer background

Conventional methods for preparing small-RNA-seq libraries by adaptor ligation generate a significant amount of adaptor dimer, thereby resulting in wasted sequencing reads. These methods also do not capture small 5'-capped and 5'-triphosphorylated RNAs. The ScriptMiner™ small-RNA-seq library preparation technology overcomes these limitations. ScriptMiner™ small-RNA libraries contain greatly reduced amounts of adaptor-dimer when compared to conventional methods and also generate coverage that is characteristic of the entire small-RNA transcriptome.

The small-RNA transcriptome contains a diverse array of RNAs, such as microRNA (miRNA), small nucleolar RNA (snoRNA), small nuclear RNA (snRNA), small interfering RNA (siRNA) and piwi-interacting RNA (piRNA). The biogenesis of these classes of RNAs results in a variety of different modifications at the 5' end: 5'-monophosphate (5' pN), 5'-triphosphate (5'pppN) or 5'-cap (GpppN). In addition, the 3' end of the RNA may be modified to include a 2'-O-methyl, 3'-OH instead of the normal 2',3'-OH.

Current methods of preparing small-RNA-seq libraries involve ligation-tagging (adaptor-ligation) of the 3' and 5' ends of the RNA, reverse transcription of the di-tagged RNA into cDNA, and PCR amplification. However, these methods suffer from two major drawbacks. First, they amplify a significant amount of adaptor-dimer that reduces the efficiency of the 5' ligation reaction and contaminates the sequencing library, leading to a large number of nonproductive sequencing reads. Second, the conventional methods do not capture small 5'-capped and 5'-triphosphorylated RNAs.

The ScriptMiner™ small-RNA-seq library preparation method significantly reduces the amount of adaptor-dimers in the library and enables the user to capture the entire small-RNA transcriptome, including small 5'-capped and 5'-triphosphorylated RNAs. The result is a more sensitive and comprehensive representation of the small transcriptome in the library.

Overview

Figure 1 presents an overview of the ScriptMiner™ process. Briefly, total RNA or size-selected RNA is tagged at its 3' end with a preadenylated 3'-adaptor oligonucleotide. A significant proportion of the excess

3'-adaptor oligonucleotide, which can form undesired adaptor-dimers, is then enzymatically removed. The default procedure provides the option to tag only 5' monophosphorylated RNAs, such as miRNA. An alternative step—treatment with tobacco acid pyrophosphatase




Figure 1 | Overview of the ScriptMiner™ small-RNA-seq library preparation procedure. The ScriptMiner™ procedure includes a novel enzymatic step to greatly reduce the amount of excess 3' adaptor and, thus, adaptor-dimers in the library. The procedure also enables the user to selectively prepare libraries from only 5'-monophosphorylated RNA or from the entire small-RNA transcriptome.

Jim Pease

EPICENTRE Biotechnologies, Madison, Wisconsin, USA.
Correspondence should be addressed to: Jim Pease. E-mail: jim.pease@epibio.com.

APPLICATION NOTES

Table 1 | Summary of sequencing data from ScriptMiner libraries with and without TAP treatment. The ScriptMiner library was prepared including TAP treatment to permit capture of the entire small-RNA transcriptome. The untreated ('No-TAP') library captured only 5'-monophosphorylated small RNA. BrRR, human brain reference RNA; HeLa, HeLa total RNA.

RNA sample	Enzyme treatment	Reads passing filter	Reads aligned to hg19	% reads aligned
BrRR	TAP	1,384,314	1,380,799	99.75
BrRR	No TAP	712,841	711,926	99.87
HeLa	TAP	1,560,105	1,555,406	99.70
HeLa	No TAP	662,365	661,362	99.85

(TAP: supplied in the kits)—enables capture of the entire small-RNA transcriptome. The di-tagged RNA is reverse-transcribed into cDNA, and the cDNA amplified by PCR. In addition to amplifying the library, the PCR also incorporates the necessary platform-specific adaptor sequences into the library and adds a barcode (index read) to the library, if desired. The PCR-amplified library is gel-purified, and the extracted small-RNA library is ready for cluster generation before sequencing.

Reduced adaptor-dimer


The ScriptMiner™ procedure uses an optimized strategy for degrading excess 3' adaptor oligonucleotide in order to suppress the formation of undesired adaptor-dimers. **Figure 2** shows the greatly reduced level of adaptor-dimer formed by the ScriptMiner™ process compared to a conventional small-RNA-seq library method. By significantly reducing the amount of adaptor-dimer, more of the desired small-RNA transcriptome is amplified by PCR, and less adaptor-dimer contaminates the final sequencing library.

The ScriptMiner™ method captures the entire small-RNA transcriptome

Following degradation of excess 3' adaptor oligonucleotide, the user has the option of preparing libraries either from small RNA with a 5'-monophosphate (such as miRNA) or from the entire small-RNA transcriptome, depending on the enzymes chosen in the degradation steps (**Fig. 1**). By treating the sample with TAP, small RNA with a 5' monophosphate, a 5' triphosphate or a 5' cap will all be included in the library. **Table 1** shows libraries prepared from only 5'-monophosphorylated RNA ('No TAP') and from TAP-treated samples. The TAP-treated library produced twice as many aligned reads as the 'No-TAP' library, indicating that the TAP-treated library captured more of the small-RNA transcriptome than conventionally produced ('No-TAP') libraries.

Conclusions

The ScriptMiner™ procedure significantly reduces the amount of adaptor-dimers in small-RNA libraries compared to conventional methods, thereby reducing the number of wasted sequencing reads. ScriptMiner™ libraries can be prepared from the whole small-RNA transcriptome, providing a more detailed picture of the regulatory processes that are mediated by small RNAs within a cell. Current



ScriptMiner™ kits permit the preparation of both non-barcoded (singleplex) and barcoded (multiplex) Illumina®-compatible libraries.

ACKNOWLEDGMENTS

We are grateful to Scott Kuersten and Ramesh Vaidyanathan for the data used in this application note.

This article was submitted to *Nature Methods* by a commercial organization and has not been peer reviewed. *Nature Methods* takes no responsibility for the accuracy or otherwise of the information provided.