

AppliChem

DNA decontamination: DNA-ExitusPlus in comparison with conventional reagents

Here we present a completely new DNA decontamination reagent DNA-ExitusPlus. In comparison with conventional products, DNA-ExitusPlus guarantees fast and efficient destruction of nucleic acids without corrosive or toxic properties. DNA-ExitusPlus was developed by AppliChem GmbH, Darmstadt, in cooperation with multiBIND biotec GmbH, Dortmund.

The polymerase chain reaction (PCR) has set the standard for DNA amplification technology¹, and the rapid extension of PCR technologies has resulted in multiple new applications^{2,3}. The latest protocols and enzyme variants for PCR reactions allow the detection of even single DNA molecules¹. One important consequence of the improved sensitivity in DNA amplification is the necessity to avoid any contamination from unwanted external DNA molecules.

The latest research in genetic technology has demonstrated that the presence of free DNA molecules can pose serious problems^{4,5}. Infections, biological transformations or recombination of genomes can be generated by free DNA plasmids or fragments³⁻⁷. Changes in viral and bacterial infectivity are observed and are related to the phenomenon of resistance against multiple antibiotics⁸. As a result, the detection of DNA contamination or prevention of amplification artifacts in PCR experiments is essential for all applications as well as biological containment and safety.

Investigations of the properties of conventional DNA decontamination reagents revealed two major problems. First, none of the reagents studied destroyed DNA molecules efficiently, and second, existing reagents typically contain components with corrosive or even toxic properties. As a consequence, we saw the necessity to develop new solutions for effective DNA decontamination.


DNA decontamination reagents use three different molecular principles for destruction or inactivation of genetic material: modification, denaturation and degradation. Safe DNA decontamination depends on the degradation of DNA into very small fragments. We developed a DNA degradation test to compare conventional decontamination reagents with the new DNA-ExitusPlus. This test allows sensitive quantification of the fragmentation process (**Fig. 1**).

Karl-Heinz Esser¹, Wolfram H Marx² & Thomas Lisowsky¹

¹multiBIND biotec GmbH, Otto-Hahn-Str. 15, D-44227, Dortmund, Germany.

²AppliChem GmbH, Ottoweg 4, D-64291 Darmstadt, Germany. Correspondence should be addressed to T.L. (thomas.lisowsky@gmx.de) or W.H.M. (w.marx@appliChem.de).

Figure 1 | Comparison of DNA degradation for selected conventional DNA decontamination reagents and DNA-ExitusPlus. For each sample, aliquots of 1 µg covalent, closed, circular (CCC) plasmid DNA dissolved in 5 µl sterile Tris buffer (1 mM; pH 8.0) were treated with 5 µl of the listed reagents, respectively, for 2 min at 20°C. The control (C) contains intact CCC plasmid DNA after treatment with sterile water. The products X1, X2 and X4 cause no detectable degradation of the test DNA. For the product X3 and D (conventional DNA-Exitus) only partial degradation is observed. Only DNA-ExitusPlus (D+) causes very rapid and nearly complete DNA degradation.



One surprising result was that some of the conventional reagents are based only on the principles of modification or denaturation. This is revealed by the complete absence of any degraded DNA molecules (**Fig. 1**). Thus the genetic information is only masked but not destroyed. By chemical demasking reactions the DNA sequences of these molecules would be available and active again. Based on our knowledge about gene technology and the principles of recombination, we concluded that these reagents are no longer appropriate.

However, even reagents that degrade DNA cause only limited and partial destruction. Hence, very large DNA fragments containing the complete genetic information still survive treatment. Only DNA-ExitusPlus achieves rapid and efficient degradation.

Another severe disadvantage of conventional reagents is revealed in a test of their corrosive potential. For this purpose different metal plates were incubated with aliquots of the reagents (**Fig. 2**). This test demonstrates that all currently available products contain aggressive chemicals with corrosive, harmful or even toxic effects. Known ingredients of conventional reagents are azides, mineral acids like

APPLICATION NOTES




Figure 2 | The corrosive potential of selected conventional DNA decontamination reagents in comparison with DNA-ExitusPlus. Aliquots of 10 μ l from each listed reagent were applied to metal plates that are typically used for laboratory materials and equipment, and then incubated for 20 min. Sterile water was used as a control (C). The reagents X2, X3 and D cause irreversible corrosion and damage to many of the metal surfaces. For DNA-ExitusPlus (D+) no damage is observed for any of the surfaces. In some cases, one observes a polishing effect owing to the removal of dirt or oxide layers.

phosphoric acid or hydrochloric acid, aggressive peroxides or strong alkaline substances like sodium hydroxide. Therefore, after only 20 min of incubation, irreversible damage of metal surfaces is observed (Fig. 2). The newly developed reagent DNA-ExitusPlus exhibits its unique characteristics especially in this test. For all metal surfaces no corrosion is observed; DNA-ExitusPlus was also tested on many different plastic surfaces without any indication of damage (data not shown). DNA-ExitusPlus offers an efficient, gentle and environmentally safe alternative and proves its superiority to other available decontamination reagents. DNA-ExitusPlus not only degrades and removes all DNA molecules with high efficiency but is neither toxic nor corrosive.

In summary, one observes the following new and unique characteristics: (i) catalytic and cooperative effects guarantee rapid nonenzymatic degradation of nucleic acids, (ii) all components of DNA-ExitusPlus are biodegradable and not harmful or toxic, (iii) no aggressive mineral acids or alkaline substances are used. Equipment and materials are not damaged or corroded even after prolonged incubation.

Currently, the most effective method for decontamination appears to be autoclaving. Under standard conditions for autoclaving, DNA molecules are degraded into fragments of 20 to 30 base pairs. However, recent investigations with highly sensitive PCR analysis demonstrate that even after autoclaving larger DNA fragments can persist⁹. Furthermore, autoclaving can only be used for decontamination of heat-resistant materials and equipment that fit into the autoclave. Decontamination of laboratory benches or larger equipment is impossible.

Efficient degradation of DNA molecules by DNA-ExitusPlus was monitored by PCR analysis (Fig. 3) proving that no amplifiable DNA templates are present. Today, only very different nonstandardized PCR tests are used as controls for successful DNA decontamination. In the case of large DNA control templates, low DNA concentrations and high dilutions in the washing steps, evidence for successful DNA decontamination is very limited. Therefore one has to be very cautious




Figure 3 | PCR test for the complete removal of DNA contamination by DNA-ExitusPlus. (a,b) Selected amounts of test DNA were lyophilized on the inner surface of PCR tubes, after which the samples were incubated for 20 s with either sterile water (a) or DNA-ExitusPlus (b). For the PCR test reaction, mixtures of 50 μ l were added to each tube. This reaction mixture contains primers for the amplification of control and test DNA. The control DNA is added to each sample and shows that the PCR reaction is not inhibited. The negative control (with sterile water) exhibits DNA bands for the test and control templates. The PCR reaction after treatment with DNA-ExitusPlus amplifies the control DNA only.

about using a single PCR test as evidence for complete decontamination because such a PCR test would also be negative in the case where DNA is only modified or masked. For complete evaluation of the potential of a DNA decontamination reagent one has to use PCR analysis in combination with a sensitive DNA degradation test.

The tests described here reveal the unique characteristics of DNA-ExitusPlus. These properties offer new opportunities for potential applications in the health sector, the life sciences, medical hygiene, food production and the household. We are convinced that this product defines a new standard for efficient, rapid and gentle DNA decontamination. According to the latest results on biological activities of free DNA molecules, such a product is critical for the new tasks concerning biological containment and safety.

- Innis, M.A. *et al.* (eds.) *PCR Protocols: A Guide to Methods and Applications* (Academic Press, Inc., San Diego, 1990).
- Cavalli-Sforza, L.L. The Human Genome Diversity Project: past, present and future. *Nat. Rev. Genet.* **6**, 333–340 (2005).
- Tumpey, T.M. *et al.* Characterization of the reconstructed 1918 Spanish influenza pandemic virus. *Science* **310**, 77–80 (2005).
- Burns, P.A. *et al.* Transformation of mouse skin endothelial cells *in vivo* by direct application of plasmid DNA encoding the human T24 H-ras oncogene. *Oncogene* **6**, 1973–1978 (1991).
- Moniz, M. *et al.* HPV DNA vaccines. *Front. Biosci.* **8**, d55–d68 (2003).
- Gibbs, M.J. *et al.* Recombination in the hemagglutinin gene of the 1918 “Spanish flu.” *Science* **293**, 1842–1845 (2001).
- Kaiser, J. Biocontainment. 1918 flu experiments spark concerns about biosafety. *Science* **306**, 591 (2004).
- Guyot, A. *et al.* Molecular epidemiology of multi-resistant *Escherichia coli*. *J. Hosp. Infect.* **43**, 39–48 (1999).
- Elhafi, G. *et al.* Microwave or autoclave treatments destroy the infectivity of infectious bronchitis virus and avian pneumovirus but allow detection by reverse transcriptase-polymerase chain reaction. *Avian Pathol.* **33**, 303–306 (2004).

This article was submitted to *Nature Methods* by a commercial organization and has not been peer reviewed. *Nature Methods* takes no responsibility for the accuracy or otherwise of the information provided.