2 © 1995 Nature Publishing Group http://www.nature.com/naturebiotechnology

/NEW PRODUCTS



Fermentor

Bioengineering AG (Wald, Switzerland) introduces a fermentor with downstream processing to allow sterile and continuous ultrafiltration of fermentor broth. The modular unit is integrated in a standard fermentor. Filtration capacity is flexible, as the number of filtration disks is variable. Installing blade stirrers in between the filtration disks minimizes fouling on the membranes.

Write in 923 on Reader Service Card.

Autosampler



The Model 708 fixed loop, positive displacement HPLC autosampler from Alcott Chromatography (Nor-cross, GA) includes a 128-vial tray, integrated high-speed injector, and "onthe-fly" programming functions. Optional refrigeration and system interfacing (IEEE, R-232, and BCD) are available.

Write in 920 on Reader Service Card.

Blood analyzer

AVL Medical Instruments (Stone, U.K.) introduces the AVL Omni analyzer, capable of measuring pH/blood gas, electrolytes, and Co-Oximeter profile from one blood sample. It is capable of analyzing up to 40 samples per hour. All user operations are performed by an interactive color-coded touch screen with icon-driven software.

Write in 924 on Reader Service Card.

INTERNATIONAL ADVERTISING MANAGER (NEW YORK) Sande T. Giaccone 1 (800) 989-1112

INTERNATIONAL ADVERTISING MANAGER (LONDON) Kathryn Wayman 44 (171) 843-4960

> EAST COAST Eric Gordon 1 (800) 989-1112 345 Park Avenue South New York, NY 10012 FAX: 1 (212) 726-9200

MIDWEST John D. Carr/Steve Loerch 1 (708) 498-4520 255 Revere Drive, Suite 106 Northbrook, IL60062 FAX: 1 (708) 498-5911 WEST COAST & SOUTHWEST Lynne Stickrod/Heidi Y. Gregory 1 (800) 229-9758 582 Market Street, Suite 1407 San Francisco, CA04104 FAX: 1 (415) 781-3805

CANADA Peter Drake 1 (416) 690-2423 17 Pine Crescent Toronto, Ontario M4E IL1 FAX: 1 (416) 868-1213

UNITED KINGDOM Kathryn Wayman 44 (171) 843-4960 Porter's South, Crinan Street London, NI 9XW FAX: 44 (171) 843-4996

SOUTHERN GERMANY & AUSTRIA Sabine Fürst 49 (89) 52-40-81 Sandstrasse 41, D80335 München FAX: 49 (89) 523-2222

Image analysis software

The BioMax ID image analysis software line from Eastman Kodak (Rochester, NY) is designed to analyze nucleic acid and protein electrophoresis gels and provide molecular weight and mass data. Acquired images can be saved as TIFF or PICT files. Images for analysis may be acquired using a TWAIN compliant scanner, a digital camera, or may be opened from any TIFF image file.

Write in 921 on Reader Service Card.

Microplate Iuminometer

Anthos Labtec (Salzburg, Austria) introduces the Anthos Lucy 1 microplate luminometer, designed for flexibility with its highly sensitive detection system, reagent injectors, random well access, and accurate thermoelectric temperature control. The luminometer incorporates a photomeric measurement unit that may be used for colormetric applications.

Write in 925 on Reader Service Card.

Mycoplasma detection

American Type Culture Collection (Rockville, MD) introduces a PCRbased mycoplasma detection kit for detecting the eight commonly encountered mycoplasma contaminants. The kit contains a first- and second-stage primer mixture and two positive DNA controls providing enough material for 50 PCR reactions in each kit.


Write in 922 on Reader Service Card.


NORTHERN GERMANY Andrea Küchler 49 (40) 3800-137 Prof-Brix-Weg 8, D-22767 Germany FAX: 49 (40) 3800-166

SCANDINAVIA Alexandra Hewardt 46 (8) 644-0005 Andrew Karnig & Associates Finnbodavagen I, 13131 NACKASweden FAX: 46 (8) 642-3150

FRANCE/HOLLAND/BELGIUM Patricia Hummel 33 (14) 887-4217 5 Rue Joseph Sansboeuf 75008 Paris France FAX: 33 (14) 387-4215

SWITZERLAND Verena Loewenthal-Hugi (41) 54 61 3030 Interpress Internationale Mediabenatung Postfach CH-8267 Berlingen FAX: (41) 54 61 3200







JAPAN Phillip Hamill 81 (3) 267-8751 Shin-Mitsuke Building (4F) 3-61chigaya Tamachi Shinjuku-ku Tokyo 162 FAX: 81 (3) 267-8746

AUSTRALIA Anne Murrey 61 (2) 299-5677 World Media Network Pty, Level 1 89 York Street, Sydney NSW 2000 FAX: 61 (2) 299-6178

> ISRAEL Heiko Sieger, Intermedia 972 (3) 528-7815 76 King George St. PO Box11820 Tel Aviv, Israel 61116 FAX: 972 (3) 299259