

- 729–730 (1997).
- Alexander, M.Y. & Akhurst, R.J. Liposome-mediated gene transfer and expression via the skin. *Hum. Mol. Genet.* **4**, 2279–2285 (1995).
 - Li, L. & Hoffman, R.M. The feasibility of targeted selective gene therapy of the hair follicle. *Nature Med.* **1**, 705–706 (1995).
 - Domashenko, A., Gupta, S. & Cotsarelis, G. Efficient delivery of transgenes to human hair follicle progenitor cells using topical lipoplex. *Nature Biotechnol.* **18**, 420–423 (2000).
 - McAllister, D.V., Allen, M.G. & Prauznitz, M.R. Microfabricated microneedles for gene and drug delivery. *Ann. Rev. Biomed. Engineering* **2**, 289–313 (2000).
 - US Patent 5928207. Microneedles with isotropically etched tips, and method of fabrication (July 27, 1999).
 - Grove, G., Grove, M., Zerweck, C. & Pierce, E. Computerized evaporimetry using the DermLab TEWL probe. *Skin Res. Technol.* **4**, 1–5 (1998).
 - Patrick, E. & Maibach, H. in *Principles and Methods of Toxicology*, 3rd edn. (ed. Hayes, W.E.) 767–803 (Raven, New York, 1994).
 - Gracely, R.H., McGrath, P.A. & Dubner, R. Validity and sensitivity of ratio scales of sensory and affective verbal pain descriptors: Manipulation of affect by diazepam. *Pain* **5**, 19–29 (1978).
 - Fan, H., Lin, Q., Morrissey, G.R. & Khavari, P.A. Immunization via hair follicles by topical application of naked DNA to normal skin. *Nature Biotechnol.* **17**, 870–872 (1999).
 - Schneck, J.P. Monitoring antigen-specific T cells using MHC-Ig dimers. *Immunol. Invest.* **29**, 163–169 (2000).
 - Raz, E. *et al.* Intradermal gene immunization: The possible role of DNA uptake in the induction of cellular immunity to viruses. *Proc. Natl. Acad. Sci. USA* **91**, 9519–9523 (1994).
 - Paul, A., Cevc, G. & Bachhawat, B.K. Transdermal immunization with an integral membrane component, gap junction protein, by means of ultradeformable drug carriers, transferrinosomes. *Vaccine* **16**, 188–195 (1998).
 - Watabe, S. *et al.* Protection against influenza virus challenge by topical application of influenza DNA vaccine. *Vaccine* **19**, 4434–4444 (2001).
 - W.H. *et al.* Topical gene delivery to murine skin. *J. Invest. Dermatol.* **112**, 370–375 (1999).
 - McCluskie, M.J. *et al.* Route and method of delivery of DNA vaccine influence immune responses in mice and non-human primates. *Mol. Med.* **5**, 287–300 (1999).

ON THE MARKET

ASSAYS



Circulating immune complexes assay.

ALPCO Diagnostics has announced the release of an **immunoassay for the quantitative determination of circulating immune complexes (CIC)** in human serum or plasma. The CIC-C1q EIA is intended to aid in diagnosing and monitoring patients with various autoimmune and other CIC-related diseases. This assay, developed by Buhlmann Laboratories of Switzerland, is a solid-phase enzyme immunoassay in which CIC present in the patient sample is first immobilized to the pre-coated microplate wells, and then incubated with an alkaline phosphatase-labeled antibody. After a second washing step, substrate is added to generate a chromophoric reaction that is subsequently stopped, and the optical density then measured with a microplate reader at 405 nm. The kit is offered in a 96-well microplate format (12 x 8-well strips) and comes with ready-to-use reagents, including six standards and a bi-level control set. This product recently received premarket FDA clearance for in vitro diagnostic use.

Tel. (+1) 800-592-5726

Fax (+1) 603-898-6854

www.alpco.com

Step-saving ELISA kits.



Requiring only a few working steps, Instant **ELISA kits** from Bender MedSystems are designed to greatly reduce the potential for experimental errors when detecting biological molecules and provide more precise results than standard ELISAs. The kits are available for several TNF receptors, interleukins and CD markers, as well as ICAMs and selectins. The company says that the Instant ELISAs eliminate many of the processing steps needed to perform an assay, including repeated washings and incubations. The components, which include standard protein, detection antibody, streptavidin HRP and sample diluent, come freeze-dried on the microtiter plate. To begin an assay, the researcher simply adds distilled water for reconstitution and then places samples in the respective wells.

Tel. (+43) 1-796-4040-0

Fax (+43) 1-796-40-40-400

www.bendermedsystems.com

IMMUNOPRECIPITATION/
ISOLATION KITS

Vector Laboratories offers a wide selection of Fusion-Aid kits for the **immunoprecipitation/isolation of tagged fusion proteins**. Each kit contains an affinity-purified antibody to the specific 'tag'. The antibody is coupled through a stable, hydrophilic spacer arm to

crosslinked 4% agarose beads. A total of 0.5 ml of settled gel is provided per kit. In addition, five spin columns and collection tubes, as well as a protocol are included. The Fusion-Aid kits are designed for the following specific tags: green fluorescent protein (GFP) recognizing wild-type, recombinant or any variants of *Aequorea victoria* GFP; maltose binding protein, c-Myc; human influenza virus hemagglutinin. Compared to traditional immunoprecipitation methods, elution of the fusion protein from the Fusion-Aid affinity matrix is not accompanied by contaminating antibody to the tag, says Vector Labs. The spin columns provided allow rapid and convenient processing with maximal protein recovery.

Tel. (+1) 650-697-3600

Fax (+1) 650-697-0339

www.vetorlabs.com

PROCESSING ENZYMES

Alexis Biochemicals offers a new technical bulletin describing the **endoproteolytic processing enzymes called proprotein convertases**. These enzymes have a specificity for basic motifs (such as RXX/R, RXXR and K/RR) in proteins in the secretory pathway, and are capable of generating bioactive proteins by limited endoproteolysis of the corresponding precursor proteins. Furin is the prototype of this novel family of processing enzymes. The bulletin contains information about the complete range of monoclonal antibodies directed against all major domains of human and mouse furin convertase. A broad selection of antisera to other mammalian proprotein convertases is also available.

Tel. (+1) 858-658-0065

Fax (+1) 858-658-9224

www.alexis-corp.com