

SCIENTIFIC REPORTS

OPEN

Accurate quantification of creatinine in serum by coupling a measurement standard to extractive electrospray ionization mass spectrometry

Received: 14 August 2015
Accepted: 07 December 2015
Published: 13 January 2016

Keke Huang¹, Ming Li², Hongmei Li², Mengwan Li², You Jiang² & Xiang Fang²

Ambient ionization (AI) techniques have been widely used in chemistry, medicine, material science, environmental science, forensic science. AI takes advantage of direct desorption/ionization of chemicals in raw samples under ambient environmental conditions with minimal or no sample preparation. However, its quantitative accuracy is restricted by matrix effects during the ionization process. To improve the quantitative accuracy of AI, a matrix reference material, which is a particular form of measurement standard, was coupled to an AI technique in this study. Consequently the analyte concentration in a complex matrix can be easily quantified with high accuracy. As a demonstration, this novel method was applied for the accurate quantification of creatinine in serum by using extractive electrospray ionization (EESI) mass spectrometry. Over the concentration range investigated (0.166 ~ 1.617 $\mu\text{g/mL}$), a calibration curve was obtained with a satisfactory linearity ($R^2 = 0.994$), and acceptable relative standard deviations (RSD) of 4.6 ~ 8.0% ($n = 6$). Finally, the creatinine concentration value of a serum sample was determined to be $36.18 \pm 1.08 \mu\text{g/mL}$, which is in excellent agreement with the certified value of $35.16 \pm 0.39 \mu\text{g/mL}$.

Renal failure (known as kidney failure) has received significant attention since even moderate reduction in kidney function is associated with high rate of mortality¹. Early diagnosis of kidney disease, followed by appropriate medical treatments, can prevent or postpone kidney failure. Creatinine, a byproduct of muscle metabolism, is removed from the blood chiefly by the kidneys. The serum creatinine level is the most widely used and commonly accepted measure of renal function in clinical medicine^{2,3}. Many techniques including electrochemical sensor⁴, Raman spectroscopy⁵, spectrophotometry^{6,7}, capillary electrophoresis^{8,9}, high performance liquid chromatography¹⁰, are available for the determination of creatinine in serum. Alternatively, hyphenated mass spectrometric methods including gas chromatography-mass spectrometry (GC/MS)^{11,12}, liquid chromatography-mass spectrometry (LC/MS)¹²⁻¹⁵ have been developed for the accurate quantification of creatinine in serum. Among these techniques, liquid chromatography-isotope dilution mass spectrometry (LC-IDMS) was recommended as the reference technique by National Kidney Disease Education Program¹⁶. However, traditional MS method suffers from the need of tedious, time- and cost-consuming sample pretreatment steps (derivatization, extraction, chromatographic separation, etc.) before complex raw samples are analysed. Recently introduced ambient ionization techniques, including desorption electrospray ionization (DESI)¹⁷⁻¹⁹, direct analysis in real time (DART)²⁰, desorption atmospheric pressure chemical ionization (DAPCI)²¹, dielectric barrier discharge ionization (DBDI)²², low temperature plasma probe (LTP)²³, easy ambient sonic-spray ionization (EASI)²⁴, and extractive electrospray ionization (EESI)²⁵⁻³⁰, allow direct desorption/ionization of chemicals in raw samples under ambient conditions. Ambient MS methods require no or minimal sample preparation and tolerate chemical contamination of an ion source. It is notable that DESI was adopted for the determination of urine creatinine level¹⁹. Nonetheless, current obstacle on the quantitative analysis of ambient mass spectrometry is the matrix effect³¹, especially in the case

¹State Key Laboratory of Inorganic Synthesis and Preparative Chemistry, College of Chemistry, Jilin University, Changchun, 130012, P. R. China. ²Chemistry Department, National Institute of Metrology, Beijing, 100013, P. R. China. Correspondence and requests for materials should be addressed to Ming Li. (email: mingutah@hotmail.com)




Figure 1. A typical EESI mass spectrum of a serum sample with 50 times dilution (inset shows an MS² spectrum of m/z 114).

that samples have high complex matrix and high quantitative accuracy is required³². For example, the signal of protonated analytes will dramatically decrease in the case: i) the proton affinities (PA) of the species in matrix are larger than PA of the analyte, ii) alkali-metals (e.g. Na⁺, K⁺) in the matrix cause the formation of alkali-metal adducts of analyte. Therefore, a standard addition method or/and isotope dilution MS were coupled to ambient ionization technique to improve the analytical accuracy^{33,34}. In this case, samples were prepared by spiking a series of stock solutions before being analyzed. For each sample analysis, at least five spiked samples had to be prepared, and consequently at least five MS analysis had to be performed^{33,34}. It is still tedious, time- and cost-consuming to prepare and to analyse the spiked samples, especially for the high throughput analysis.

In this study, a simple, rapid ambient MS method was developed for accurate analysis of creatinine in serum by coupling a measurement standard to EESI.

Results

EESI-MS spectra of creatinine.

To optimize the parameters of the instrument, a creatinine standard solution (0.1 μg/mL) was introduced into EESI source and a typical mass spectrum of creatinine was recorded (see Figure S1). The signal at m/z 114 corresponds to protonated creatinine molecules. The product ion spectrum of the mass selected ions of m/z 114 (inset of Figure S1) generated ionic fragments of m/z 86 by the neutral loss of CO³⁵. A serum sample (GBW 09170) diluted 50 times with water was directly analysed by EESI-MS, and the corresponding mass spectrum was recorded (see Fig. 1). The signal at m/z 114 corresponds to protonated creatinine molecules, which was confirmed by tandem mass spectrum (inset of Fig. 1). The signals at m/z 61, 83, 136 correspond to protonated urea, sodiated urea, and sodiated creatinine, respectively. Actually, sodium ion adducts are very common in MS³⁶, and sodiated creatinine was also observed in ambient ionization mass spectrometry such as DESI¹⁹, EESI³⁰, and contactless atmospheric pressure ionization³⁷.

Demonstration of the matrix effect.

A series of creatinine standard solutions prepared with pure water were used to establish a calibration curve (see Figure S2). Afterwards, 50 fold dilution of GBW09171 was analysed by using EESI-MS. The creatinine level in the original serum sample was determined to be 13.96 ± 0.42 μg/mL, which is significantly lower than the certified value of 35.16 ± 0.39 μg/mL¹⁵. This is because matrix effect still exists during the ionization process even the serum sample was diluted 50 times with water. In addition, distribution of ion current between protonated creatinine (m/z 114) and sodiated creatinine (m/z 136) can also be a source of error.

Establishment of a calibration curve with a matrix reference material.

For reliable, accurate determination of creatinine level in serum sample, a matrix reference material, which is a particular form of measurement standard, was coupled to EESI technique. GBW09170 was chosen for constructing a calibration curve as it has a lower creatinine level. Five standard solutions, prepared with the matrix reference material, were analysed by using EESI-MS in MS² mode. The intensity of quantitative ions of m/z 86, derived from the parent ions of m/z 114, was plotted as a function of creatinine concentration (see Fig. 2). Over the concentration range investigated (0.166 ~ 1.617 μg/mL), a calibration curve was obtained with a satisfactory linearity ($R^2 = 0.994$), and acceptable relative standard deviations (RSD) of 4.6 ~ 8.0% ($n = 6$). The regression equation of the calibration curve was given as $y = 95.5x + 2.93$. The limit of quantification ($S/N = 10$) was calculated to be 0.03 μg/mL (i.e. 1.5 μg/mL in original serum sample). So the method reported here can be adopted for analysing the serum samples with creatinine level ranging 1.5 ~ 80.85 μg/mL, which is a suitable range for clinical application³⁸. To assess the inter-day variation of the response, five calibration curves were established in five consecutive days with the same series of standard solutions (see Figure S3). The relative standard deviation of five slopes was calculated to be 0.49%, however it is suggested that a calibration curve should be re-established each time the EESI source is re-assembled. Here, the standard solutions were prepared with a matrix reference material, which has similar chemical characteristics to the samples being tested¹⁵. Therefore, this calibration curve can be used to accurately interpolate the creatinine level in serum samples.




Figure 2. Calibration curve for determination of creatinine in serum with a transition ion pair m/z 114-86 as quantitative ions (Error bars designate the standard deviation, $n = 6$).

Sample number	conc. by LC-IDMS ($\mu\text{g/mL}$)	conc. by EESI-MS ($\mu\text{g/mL}$)	Relative deviation (%)
1	7.67	7.42	3.26
2	9.03	8.85	2.00
3	22.52	23.19	2.98
4	33.34	32.01	0.99

Table 1. Comparison of creatinine level in serum by using LC-IDMS and EESI-MS.

Quantification of creatinine in serum samples. A serum sample (GBW 09171) diluted 50 times with water was directly analysed by using EESI. The signal intensity in MS^2 mode was recorded and the creatinine level in the original sample was calculated to be $36.18 \pm 1.08 \mu\text{g/mL}$. This value is in excellent agreement with the certified value of $35.16 \pm 0.39 \mu\text{g/mL}$, which was confirmed by an international comparison¹⁵. To assess the inter-assays variation, the same sample was analysed after one week, and the result was calculated to be $35.83 \pm 1.01 \mu\text{g/mL}$. The relative deviation between the results obtained from two different days is less than 0.5%. To further confirm the accuracy of this novel method, four serum samples were quantified via the proposed method and a reference method (i.e. LC-IDMS) developed previously in our lab¹⁵. The relative deviations of the results from these two methods were calculated to be 0.99–3.26%, as listed in Table 1.

Discussion

Matrix effects are a current obstacle on the quantitative analysis of ambient mass spectrometry because of ion suppression, and distribution of ion current between protonated and sodiated molecules. To resolve this problem, a matrix reference material, which is considered as a “micro-ruler” to check the accuracy of tests and analytical procedures, was coupled to EESI technique. As a demonstration, this method was applied for the accurate quantification of creatinine in serum. The results show that the accurate analyte concentration in a complex matrix can be easily obtained with only one analysis for each sample by coupling a measurement standard to EESI.

Methods

Instrumentation and Working Conditions. All the experiments were performed using a Thermo Finnigan LTQ mass spectrometer (San Jose, CA, USA) equipped with a homemade EESI source, which is described elsewhere^{24–29}. A schematic of EESI source is shown in Fig. 3. A methanol/water/formic acid (v/v/v, 50:50:0.1) solution was injected into channel 1 using a syringe pump with a rate of $5 \mu\text{L}/\text{min}$ and nebulized with a sheath gas (N_2). A high voltage of 4 kV was applied to channel 1 to generate the primary ions. Samples or standard solutions were injected into channel 2 using another syringe pump with a rate of $5 \mu\text{L}/\text{min}$ and also nebulized with a sheath gas (N_2). The neutral analytes were ionized when intersecting with the primary ions. The angle (α) formed between the channel 1 and the channel 2 was 60° . The angle (β) formed between the channel 1 and the MS inlet was 150° . The distance (a) between the emitters of these two channels was 2 mm. The distance (b) between the EESI source and the MS inlet was 2.5 cm. The temperature of the ion entrance capillary was maintained at 390°C during the experiment. For tandem mass spectrometry, the precursor ions of interest were isolated with a window width of 1.5 mass/charge units (full-width) and helium was used as a collisional gas. MS^2 experiments were performed with an activation time of 40 ms and collision energy of 24% (arbitrary unit). All the mass spectra were recorded using Xcalibur[®] software. The mean of the mass spectra accumulated for 1.5 min was recorded as one measurement result. An average of 6 measurements for one sample was used as the final result.

Statement. All experimental protocols were approved by the ethics committee of National Institute of Metrology, China and adhered to the tenets of the Declaration of Helsinki. Additionally, the benefits and risks of




Figure 3. Schematic of the EESI source coupling with a measurement standard.

this study were clearly presented to the volunteers, and thereafter the written consent was obtained if they agreed to join the study.

Samples Preparation. Methanol (HPLC grade), and formic acid were purchased from Chinese Chemical Reagent Co., Ltd. (Shanghai, China). Ultrapure water ($18.2 \text{ M}\Omega \cdot \text{cm}^{-1}$) was prepared with Milli-Q Direct water purification system (Millipore, USA). Pure creatinine (SRM914a) with a purity of 99.7% was purchased from national institute of standards and technology (NIST), USA. Serum matrix reference materials GBW 09170 and GBW 09171 with certified creatinine values of $8.30 \pm 0.10 \mu\text{g/mL}$ and $35.16 \pm 0.39 \mu\text{g/mL}$, were developed in our lab (i.e. national institute of metrology, China). Four serum samples, collected from volunteers including two healthy persons and two patients with kidney disease, were provided by the academy of military medical sciences (Beijing, China). All the chemicals were used directly without any further treatment.

References

1. United States Renal Data System, 2014 Annual Data Report: Epidemiology of kidney disease in the United States. National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases (Bethesda, MD 2014).
2. Narayanan, S. & Appleton, H. D. Creatinine: a review. *Clin. Chem.* **26**, 1119–1126 (1980).
3. Perrone, R. D., Madias, N. E. & Levey, A. S. Serum creatinine as an index of renal function: new insights into old concepts. *Clin. Chem.* **38**, 1933–1953 (1992).
4. Hassan, S. S. M., Elhemma, E. M. & Mohamed, A. H. K. Novel biomedical sensors for flow injection potentiometric determination of creatinine in human serum. *Electroanal.* **17**, 2246–2253 (2005).
5. Stosch, R., Henrion, A., Schiel, D. & Guttler, B. Surface enhanced Raman scattering based approach for quantitative determination of creatinine in human serum. *Anal. Chem.* **77**, 7386–7392 (2005).
6. Chromy, V., Rozkosna, K. & Sedlak, P. Determination of serum creatinine by Jaffe method and how to calibrate to eliminate matrix interference problems. *Clin. Chem. Lab. Med.* **46**, 1127–1133 (2008).
7. Guterres, M. V., Volpe, P. O. L. & Ferreira, M. M. C. Multiway calibration for creatinine determination in human serum using the Jaffe reaction. *Appl. Spectrosc.* **58**, 54–60 (2004).
8. Zinellu, A., Carru, C., Usai, M. F., Sotgia, S. & Deiana, L. Determination of creatinine in human serum by short-end injection capillary zone electrophoresis. *Electrophoresis* **25**, 1096–1101 (2004).
9. Paroni, R. *et al.* Creatinine determination in serum by capillary electrophoresis. *Electrophoresis* **25**, 463–468 (2004).
10. Krcmova, L. *et al.* Determination of neopterin, kynurenine, tryptophan and creatinine in human serum by high throughput HPLC. *Talanta* **85**, 1466–1471 (2011).
11. Fernández-Fernández, M. *et al.* Development of an isotope dilution GC-MS procedure for the routine determination of creatinine in complex serum samples. *Clin. Chim. Acta* **431**, 96–102 (2014).
12. Fernández-Fernández, M. *et al.* Simultaneous determination of creatinine and creatine in human serum by double-spike isotope dilution liquid chromatography-tandem mass spectrometry (LC-MS/MS) and gas chromatography-mass spectrometry (GC-MS). *Anal. Chem.* **87**, 3755–3763 (2015).
13. Rasoul, A. K. *et al.* Determination of serum creatinine by liquid chromatography isotope dilution tandem mass spectrometry (LC-IDMS/MS): a candidate reference method. *Clin. Biochem.* **41**, 1273–1273 (2008).
14. Stokes, P. & O'Connor, G. Development of a liquid chromatography-mass spectrometry method for the high-accuracy determination of creatinine in serum. *J. Chromatogr. B* **794**, 125–136 (2003).
15. Dai, X. *et al.* Certified reference materials (GBW09170 and 09171) of creatinine in human serum. *J. Chromatogr. B* **879**, 429–435 (2011).
16. Randviir, E. P. & Banks, C. E. Analytical methods for quantifying creatinine within biological media. *Sensor. Actuat. B-Chem.* **183**, 239–252 (2013).
17. Takáts, Z., Wiseman, J. M., Gologan, B. & Cooks, R. G. Mass spectrometry sampling under ambient conditions with desorption electrospray ionization. *Science* **306**, 471–473 (2004).
18. Takáts, Z., Wiseman, J. M. & Cooks, R. G. Ambient mass spectrometry using desorption electrospray ionization (DESI): instrumentation, mechanisms and applications in forensics, chemistry, and biology. *J. Mass Spectrom.* **40**, 1261–1275 (2005).
19. Devenport, N. A. *et al.* Direct determination of urinary creatinine by reactive-thermal desorption-extractive electrospray-ion mobility-tandem mass spectrometry. *Anal. Chem.* **86**, 357–361 (2014).
20. Cody, R. B., Laramée, J. A. & Durst, H. D. Versatile new ion source for the analysis of materials in open air under ambient conditions. *Anal. Chem.* **77**, 2297–2302 (2005).
21. Li, M. *et al.* Document authentication at molecular level using desorption atmospheric pressure chemical ionization mass spectrometry imaging. *J. Mass Spectrom.* **48**, 1042–1049 (2013).

22. Na, N., Zhao, M., Zhang, S., Yang, C. & Zhang, X. Development of a dielectric barrier discharge ion source for ambient mass spectrometry. *J. Am. Soc. Mass Spectrom.* **18**, 1859–1862 (2007).
23. Harper, J. D., Charipar, N. A., Mulligan, C. C., Cooks, R. G. & Ouyang, Z. Low temperature plasma probe for ambient desorption ionization. *Anal. Chem.* **80**, 9097–9104 (2008).
24. Haddad, R., Milagre, H. M. S., Catharino, R. R. & Eberlin, M. N. Easy ambient sonic-spray ionization mass spectrometry combined with thin-layer chromatography. *Anal. Chem.* **80**, 2744–2750 (2008).
25. Chen, H. *et al.* Sensitive detection of native proteins using extractive electrospray ionization mass spectrometry. *Angew. Chem. Int. Ed.* **49**, 3053–3056 (2010).
26. Chen, H. & Zenobi, R. Neutral desorption sampling of biological surfaces for rapid chemical characterization by extractive electrospray ionization mass spectrometry. *Nat. Protoc.* **3**, 1467–1476 (2008).
27. Li, M. *et al.* Extractive electrospray ionization mass spectrometry towards in situ analysis without sample pretreatment. *Anal. Chem.* **81**, 7724–7731 (2009).
28. Pan, S. *et al.* Quantitative detection of nitric oxide in exhaled human breath by extractive electrospray ionization mass spectrometry. *Sci. Rep.* **5**, 8725 (2015).
29. Li, M. *et al.* Facilitated diffusion of acetonitrile revealed by quantitative breath analysis using extractive electrospray ionization mass spectrometry. *Sci. Rep.* **3**, 1205 (2013).
30. Liang, D., Fang, X., Li, M., Chinglin, K. & Li, H. Direct determination of creatinine in urine and analysis of pure aniline by extractive electrospray ionization mass spectrometry. *Anal. Lett.* **48**, 2002–2010 (2015).
31. Shelley, J. T. & Hieftje, G. M. Ionization matrix effects in plasma-based ambient mass spectrometry sources. *J. Anal. At. Spectrom.* **25**, 345–350 (2010).
32. Albert, A., Shelley, J. T. & Engelhard, C. Plasma-based ambient desorption/ionization mass spectrometry: state-of-the-art in qualitative and quantitative analysis. *Anal. Bioanal. Chem.* **406**, 6111–6127 (2014).
33. Li, X. *et al.* Direct quantification of creatinine in human urine by using isotope dilution extractive electrospray ionization tandem mass spectrometry. *Anal. Chim. Acta* **748**, 53–57 (2012).
34. Cai, J., Li, M., Xiong, X., Fang, X. & Xu, R. Detection of histamine in beer by nano extractive electrospray ionization mass spectrometry. *J. Mass Spectrom.* **49**, 9–12 (2014).
35. Hušková, R., Chrastina, P., Adam, T. & Schneiderka, P. Determination of creatinine in urine by tandem mass spectrometry. *Clin. Chim. Acta* **350**, 99–106 (2004).
36. Perry, R. H. *et al.* Capturing fleeting intermediates in a catalytic C-H amination reaction cycle. *Proc. Natl. Acad. Sci.* **109**, 18295–18299 (2012).
37. Hsieh, C., Meher, A. K. & Chen, Y. Automatic sampling and analysis of organics and biomolecules by capillary action-supported contactless atmospheric pressure ionization mass spectrometry. *PLoS ONE* **8**, e66292 (2013).
38. Thongprayoon, C. *et al.* Optimum methodology for estimating baseline serum creatinine for the acute kidney injury classification. *Nephrology* **20**, 881–886 (2015).

Acknowledgements

This work was jointly supported by National Natural Science Foundation of China (No. 21275134) and Ministry of Science and Technology of China (No. 2014DFA31810 and 2011YQ06008403)

Author Contributions

M.L. conceived the experiments, K.H. and H.L. conducted the experiments, Y.J. and X.F. analysed the results, M.L., K.H. and M.W.L. wrote the manuscript. All authors reviewed the manuscript.

Additional Information

Supplementary information accompanies this paper at <http://www.nature.com/srep>

Competing financial interests: The authors declare no competing financial interests.

How to cite this article: Huang, K. *et al.* Accurate quantification of creatinine in serum by coupling a measurement standard to extractive electrospray ionization mass spectrometry. *Sci. Rep.* **6**, 19283; doi: 10.1038/srep19283 (2016).



This work is licensed under a Creative Commons Attribution 4.0 International License. The images or other third party material in this article are included in the article's Creative Commons license, unless indicated otherwise in the credit line; if the material is not included under the Creative Commons license, users will need to obtain permission from the license holder to reproduce the material. To view a copy of this license, visit <http://creativecommons.org/licenses/by/4.0/>