

Array-based comparative genomic hybridization: clinical contexts for targeted and whole-genome designs

Swaroop Aradhya, PhD, and Athena M. Cherry, PhD

Array-based comparative genomic hybridization is ushering in a new standard for analyzing the genome, overcoming the limits of resolution associated with conventional G-banded karyotyping. The first genomic arrays were based on bacterial artificial chromosome clones mapped during the initial phases of the Human Genome Project. These arrays essentially represented multiple fluorescence in situ hybridization assays performed simultaneously. The first arrays featured a targeted design, consisting of hundreds of bacterial artificial chromosome clones limited mostly to genomic regions of known medical significance. Then came whole-genome arrays, which contained bacterial artificial chromosome clones from across the entire genome. More recently, alternative designs based on oligonucleotide probes have been developed, and all these are high-density whole-genome arrays with resolutions between 3 and 35 kb. Certain clinical circumstances are well suited for investigation by targeted arrays, and there are others in which high-resolution whole-genome arrays are necessary. Here we review the differences between the two types of arrays and the clinical contexts for which they are best suited. As array-based comparative genomic hybridization is integrated into diagnostic laboratories and different array designs are used in appropriate clinical contexts, this novel technology will invariably alter the testing paradigm in medical genetics and will lead to the discovery of novel genetic conditions caused by chromosomal anomalies. **Genet Med 2007;9(9):553–559.**

Key Words: array comparative genomic hybridization, bacterial artificial chromosome, oligonucleotide, copy number variant, whole genome array, targeted array, fluorescence in situ hybridization

Since 1956, when Tjio and Levan made the landmark observation that human cells contain 46 chromosomes, the practice of cytogenetics has witnessed many innovations, including chromosome banding techniques, fluorescence in situ hybridization (FISH), spectral karyotyping, and comparative genomic hybridization (CGH).¹ Of these advances, G-banded karyotyping is the most useful method for whole-genome analysis and has hence been the gold standard for analyzing chromosomes for the past several decades. As these advances helped illuminate the causes of a range of genetic disorders and cancers, experts in the field formulated criteria for their use in clinical diagnostics. The most recent and perhaps the most consequential advancement yet, microarray-based CGH (array CGH), finally offers a method for whole-genome analysis at resolutions much greater than that possible with conventional karyotyping.² This technology holds considerable promise for the field of medical genetics, but several stipula-

tions must be appreciated as it begins to take its place in cytogenetics laboratories.

BACTERIAL ARTIFICIAL CHROMOSOME (BAC) ARRAYS

Completion of the Human Genome Project and improvements in microarray technology led to the development of the first “genomic arrays,” which differed from previous arrays that consisted mainly of expressed sequence probes. The first genomic arrays, developed by various academic centers and commercial entities (Table 1), incorporated large clones derived from the initial physical mapping stages of the Human Genome Project.^{3–5} Most of these were BAC clones that were larger and more stable than alternative types of clones. The first BAC arrays were “targeted” formats, containing proprietary collections of clones from specific genomic regions, designed to detect deletions or duplications associated with known genetic disorders.^{3,4} These targeted arrays were followed by “whole-genome” arrays, which contained clones selected at regular intervals across the entire genome. For example, one commercial whole-genome BAC array contains about 2600 clones spread across the genome with an average distance of 1 Mb between clones.^{6,7} Some academic centers have also developed whole-genome arrays with similar resolution, as well as tiling arrays that consist of overlapping clones spanning the entire genome.^{8–10} Irrespective of the design of BAC arrays, their underlying concept is the same. These

From the Department of Pathology, Stanford University School of Medicine, Palo Alto, California.

Swaroop Aradhya, PhD, GeneDx, 207 Perry Parkway, Gaithersburg, MD 20877. Tel: 301-519-2100. Email: swaroop@genedx.com

Disclosure: The authors are members of the Scientific Advisory Committee for Cytogenetics at Affymetrix, Inc. (Santa Clara, CA).

Submitted for publication March 27, 2007.

Accepted for publication June 20, 2007.

DOI: 10.1097/GIM.0b013e318149e354

Table 1
Examples of current array CGH products and services

Vendor	Type of array	Resolution	Feature type
Baylor College of Medicine, Houston, TX	Targeted array	Tiling path in regions of known medical significance. Resolution in these regions = BAC length (70–200 kb). Oligo array under development	BACs
Wellcome Trust Sanger Institute, Hinxton, Cambridge, UK	Targeted array	Tiling path in regions of known medical significance. Resolution in these regions = BAC length (70–200 kb)	BACs
	Whole-genome array	Tiling path across genome: resolution approximately 70–200 kb	
Signature Genomics, Inc., Spokane, WA	Targeted arrays: SignatureChip, MarkerChip	Tiling path in regions of known medical significance. Resolution in these regions = BAC length (70–200 kb)	BACs
Spectral Genomics, Inc., Houston, TX	Targeted array: Constitutional Chip	Tiling path in regions of known medical significance. Resolution in these regions = BAC length (70–200 kb)	BACs
	Whole-genome array: Spectral Chip 2600	Averaged resolution: 1-Mb resolution across entire genome	
Agilent Technologies, Inc., Santa Clara, CA	All are whole-genome arrays:	Averaged resolution (gene centric):	60-mer oligo
	44B (43,000+ probes)	35 kb	
	105K (99,000+ probes)	15 kb	
	244K (236,000+ probes)	6.4 kb	
Affymetrix, Inc., Santa Clara, CA	All are whole-genome and SNP-based ^a :	Averaged resolution:	25-mer oligo
	100K mapping array	26 kb	
	500K mapping array	5 kb	
	Genome-wide human SNP array 5.0	4 kb	
Illumina, Inc., San Diego, CA	All are whole-genome and SNP based ^a :	Averaged resolution:	50-mer oligos on beads
	Infinium HumanHap300 (318,000 probes)	5 kb	
	Infinium HumanHap550 (555,000 probes)	2.8 kb	
	Infinium HumanHap650 (655,000 probes)	2 kb	
Nimblegen, Inc., Madison, WI	All are whole-genome arrays (385,000 probes):	Averaged resolution (repeats masked):	50–75-mer oligo
	Single array	6.27 kb	
	4-set arrays	1.57 kb	
	8-set arrays	713 bp	

CGH, comparative genomic hybridization; BAC, bacterial artificial chromosome; SNP, single nucleotide polymorphism.

^aRequire whole genome amplification before DNA labeling.

arrays essentially represent FISH experiments en masse and can be used to screen for many genetic disorders simultaneously.

OLIGONUCLEOTIDE ARRAYS

Soon after the first BAC arrays were introduced, various companies developed arrays containing single-stranded oligonucleotide probes (Table 1). Unlike BAC arrays, oligonucleotide arrays do not have the same underlying design concept. Variations include manufacturing method, sample throughput capacity, probe length (25–80 bases), and presence or ab-

sence of single nucleotide polymorphisms (SNPs) within the probes.¹¹ Although a discussion of these differences is important and has relevance to diagnostic testing, it is beyond the scope of this review. However, it is worthwhile to emphasize that oligonucleotide probes are constructed in situ using different patented technologies that allow the arrays to be quickly customized. All available oligonucleotide platforms are whole-genome arrays, but they can be customized either to function as targeted arrays (similar to available targeted BAC arrays) or to substantially increase resolution in a specific genomic region of interest. This versatility will eventually prove useful as array-CGH technology becomes widespread.

RESOLUTION OF TARGETED VERSUS WHOLE-GENOME ARRAYS

The resolution of any genomic array is defined by the length of the individual probes and the distance between adjacent probes on that array. Classifying the resolution on an array as macroresolution or microresolution may also be useful. We define macroresolution as the average space between probes across the entire genome, and it applies only to whole-genome arrays. Microresolution would be the density of probes within a specific region of interest and is relevant to both targeted and whole-genome arrays. For example, a targeted BAC array may have a 100-kb microresolution in the 1p36 region, whereas a whole-genome oligonucleotide array may have a 150-kb macroresolution across the entire genome, but a 10-kb microresolution within the critical region for Williams syndrome. Higher microresolution gives users greater confidence in observed data if multiple probes within a region indicate the same result. Ultimately, the macro- and microresolution on an array have direct bearing on the use of array CGH for surveying regions associated with known genomic disorders, for identifying novel chromosomal anomalies, and for refining the boundaries of observed chromosomal imbalances.

CLINICAL CONTEXTS FOR TARGETED ARRAYS

Targeted BAC arrays were originally designed to test for microdeletion and microduplication syndromes (genomic disorders) and for gains or losses of subtelomeric DNA with a single assay (Table 2). The medical genetics community has, however, debated the utility of such an assay. Some downplay the utility of targeted BAC arrays, stating that most conditions that they test for can be clinically diagnosed under most circumstances and that a simple FISH test can usually verify those diagnoses. Others argue that a clinical diagnosis cannot be made with certainty in every situation because phenotypic variability in many disorders is not fully appreciated, even within well-characterized syndromes. For instance, a targeted array may help clarify the underlying cause if a patient is too young to have manifested the full phenotype for a disorder, if the disorder shows significant variable expressivity, or if a single patient shows characteristics of multiple disorders. For these reasons and because targeted BAC arrays can also detect subtelomeric abnormalities, they have been in considerable demand since their inception.

The American College of Medical Genetics recommends high-resolution karyotyping and subtelomere-targeted FISH as part of the diagnostic workup for children with developmental delay.¹³ These patients often demonstrate mental retardation as they grow, and many suffer from congenital anomalies. Conventional subtelomere-targeted FISH detects cytogenetic imbalances in about 2.5% of these patients.¹⁴ Targeted BAC arrays have already uncovered subtelomeric anomalies in many patients whose results from subtelomere-targeted FISH were normal.^{4,15} This is because pathogenic rearrangements can occur along greater lengths of subtelomeric regions than

Table 2

Clinical contexts for targeted versus high-resolution whole-genome array CGH		
Clinical context	Targeted array	High-resolution whole-genome array
Whole chromosome aneuploidy	Yes	Yes
Subtelomere anomalies	Yes (see comment) ^a	Yes
Characterized microdeletion syndromes	Yes	Yes
Unknown microdeletions/duplications	Lower chance of detection ^a	Higher chance of detection ^a
Marker chromosomes	Yes	Yes
Uniparental disomy	No	Yes ^b
Unbalanced translocations	Yes (with subtelomeric coverage) ^a	Yes
Defining unbalanced translocation boundaries	Yes, less precise (see comment) ^a	Yes, with better boundary definition ^a
Translocation breakpoint anomalies (del/dup)	Lower chance of detection ^a	Higher chance of detection ^a
Insertions	Lower chance of detection ^a	Higher chance of detection ^a
Balanced translocations ^c	No ^c	No ^c
Inversions	No	No

This table applies only to “over-the-counter” arrays, not to customized arrays. CGH, comparative genomic hybridization.

^aDepends on probe coverage. With subtelomere anomalies, some targeted arrays may not detect the proximal boundary, for example, if a large deletion extends beyond the local subtelomeric coverage. A targeted array may define translocation breakpoint anomalies (i.e., local deletions or duplications at the breakpoint) if it has coverage at the breakpoint, although anomalies at translocation breakpoints are typically small (few base pairs to few kilobases) and may be detected only by an oligonucleotide array. A targeted array may also define the boundary of a translocation on a derivative chromosome if it has sufficient coverage at the breakpoint, but with less precision than an oligonucleotide array.

^bOnly on single nucleotide polymorphism–based oligonucleotide arrays.

^cDetection of balanced translocations would be possible if derivative chromosomes were flow-sorted before array CGH (see Gribble et al.¹²).

those covered by the current set of FISH probes used to analyze subtelomeric rearrangements. For example, 1p36 deletion syndrome (OMIM #607872) is caused by large deletions that vary in size and location within 1p36. Targeted BAC arrays with high microresolution within 1p36 have defined the boundaries for these deletions and provided insight into genotype-phenotype correlations.^{16,17}

Targeted BAC arrays are also useful for classifying supernumerary marker chromosomes (Table 2).^{18–20} Conventional methods for analyzing marker chromosomes are time-consuming and labor-intensive and require complex cocktails of FISH probes. With these methods, identifying the exact origin of a marker chromosome is difficult, and determining the amount of euchromatin present is harder still. This has direct relevance to predicting phenotypic outcome in those individ-

uals carrying marker chromosomes. Provided that the marker chromosome is present in a reasonable percentage of cells ($\geq 25\%$ mosaicism to 100% of cells), targeted BAC arrays with adequate microresolution in the pericentromeric regions of chromosomes can more easily identify the origin of a marker chromosome and also define its euchromatic segments. However, it is important to recognize that large BAC clones in the pericentromeric regions can contain repeats, and, unless each clone on a targeted array is carefully selected to ensure that it maps to the appropriate chromosome and produces signals with good specificity, results may be equivocal and require follow-up.¹⁸ Targeted BAC arrays are also limited in their ability to define euchromatic segments to intervals smaller than the BACs themselves, which are 80 to 200 kb.

Targeted BAC arrays also have some more general advantages. First, they can be used to quickly examine regions of known clinical relevance in a multiplex fashion. Second, when a deletion or duplication is identified, confirmation is straightforward because the clones that initially detected the abnormality can themselves be used as FISH probes. Third, well-designed targeted arrays generally detect only abnormalities that have been previously described (with the exception of some subtelomeric rearrangements), hence making data interpretation straightforward. On some targeted arrays, sparse coverage outside regions of known medical relevance can identify novel rearrangements, but this is generally infrequent.

A major limitation of targeted arrays, of course, is that when clinical signs and symptoms do not fit any of the disease-associated chromosomal anomalies for which the targeted array is intended to screen, then the arrays cannot be used to survey other regions of the genome in detail. For instance, children with global developmental delay are typically evaluated by high-resolution karyotyping and subtelomere-targeted FISH analysis.¹³ In this regard, targeted BAC arrays are useful for detecting some anomalies in subtelomeric regions, but they lack even coverage in the rest of the genome and therefore cannot address the high-resolution karyotyping component. This is a context in which whole-genome arrays prove indispensable.

CLINICAL CONTEXTS FOR WHOLE-GENOME ARRAYS

At the simplest level, whole-genome arrays offer everything that targeted arrays offer and more. However, the benefits of that "more" in a diagnostics setting are under debate, as the advantages of whole-genome arrays are accompanied by specific stipulations. Table 2 summarizes the appropriate clinical contexts for using the two types of arrays. The major advantage to using whole-genome arrays is clear. Geneticists must often evaluate conditions that do not fit the criteria for known genetic disorders but that are strongly suspected to have a chromosomal component. In these cases, whole-genome inspection allows an in-depth survey that could lead to the discovery of novel chromosomal imbalances. Numerous novel pathogenic deletions and duplications have already been identified using a whole-genome approach.^{5,8,10,21–24} An excellent exam-

ple is the recently discovered 17q21 microdeletion syndrome (OMIM #610443). The deletion is < 1 Mb and was identified by four separate groups using whole-genome array CGH.^{21,25,26} This disorder is thought to account for as many as 1% of all cases of mental retardation.^{25,27} Another example of a cytogenetic anomaly discovered with whole-genome array CGH is the deletion responsible for Nablus masklike facial syndrome (OMIM #608156), a multiple malformation syndrome.²⁸ Several more novel disorders will be discovered in the near future as whole-genome array CGH is more widely applied.

High-resolution whole-genome arrays are of significant interest to the cytogenetics community because of their genome-wide assaying capacity as well as versatility. Commercially available whole-genome BAC arrays have moderate resolution (1 Mb), and no high-resolution BAC arrays (such as tiling arrays) are available on the market. Thus, high-resolution whole-genome arrays will be defined as high-density oligonucleotide arrays for the purposes of this review.

One major advantage of high-resolution whole-genome arrays is that they can potentially define deletion or duplication boundaries with greater precision than BAC arrays. Knowing the extent of some cytogenetic alterations may provide insight into genotype-phenotype correlations. One illustrative example is the multiple malformation syndrome, Greig cephalopolysyndactyly syndrome (GCPS; OMIM #175700), which is caused by haploinsufficiency for *GLI3* and adjacent genes in 7p14.1. The prognosis in GCPS patients depends on the size of the deletion, which is difficult to determine by FISH or DNA sequencing. Focused array CGH with very high microresolution (730 bp) in the *GLI3* region has recently been applied to a cohort of GCPS patients to accurately define different deletion boundaries.²⁹

Most contiguous gene deletion or duplication syndromes occur due to segmental duplications that predispose a specific region of the genome to rearrangement.³⁰ Hence, patients with one of these disorders usually carry rearrangements with the same breakpoints. Only a few disorders, such as GCPS, have been characterized by deletions or duplications of varying sizes, but more likely exist, particularly in complex repeat-rich genomic regions that are vulnerable to rearrangements. Indeed, part of the reason that more such disorders have not been identified may be because sensitive technologies like high-resolution whole-genome array CGH were not available to precisely define the deletion or duplication intervals. Moreover, even some well-known genomic disorders may have subtle differences in deletion sizes that will come to light as array CGH is more widely used. Hence, disorders mapped to gene-dense regions or segmental duplication-rich regions²⁷ and those linked to large genes would be amenable to focused high-resolution array CGH analysis. Some of these disorders include Duchenne and Becker muscular dystrophies (OMIM #310200, 300376), α -thalassemia (OMIM #141800), 1p36 deletion syndrome (OMIM #607872), Prader-Willi and Angelman syndromes (OMIM #176270, 105830), DiGeorge/velocardiofacial syndrome (OMIM #188400, 192430).

Another advantage of high-resolution whole-genome array CGH is that it may be useful for evaluating the breakpoints in apparently “balanced” translocations. A proportion of translocations that appear to be balanced by conventional karyotyping are actually unbalanced with small losses or gains of DNA at the breakpoints. While it may prove useful for evaluating a patient’s phenotype, investigating such breakpoints with conventional methods can be difficult and labor-intensive. High-resolution whole-genome array CGH can potentially identify such alterations, and customized oligonucleotide arrays with high microresolution in the regions of interest can quickly define deletions or duplications at the translocation breakpoints. Numerous studies have already emphasized the value of this approach.^{31–34}

The versatility of high-density oligonucleotide arrays is a function of in situ probe synthesis chemistry, the digital workflow that controls that synthesis, and creative feature and array formats. Users can generate customized arrays quickly for two possible clinical uses. First, users may prepare focused high-resolution arrays to evaluate findings on a high-resolution whole-genome array. This may serve as a confirmation method [similar to multiplex ligation-dependent amplification or quantitative polymerase chain reaction (qPCR)] and provide more detail about the cytogenetic imbalance in question. Second, it is only a matter of time before users begin adding probes in various submicroscopic regions, such as individual exons, promoter regions, and imprinting centers, to detect small deletions or duplications that may be pathogenic. This will be particularly true for large genes. One study has already confirmed the feasibility of this approach by detecting deletions in the *DMD* gene, which is mutated in Duchenne and Becker muscular dystrophies (OMIM #310200, 300376), and in other disease-associated genes.³⁵ Genomic regions that contain imprinting centers, such as that associated with Prader-Willi and Angelman syndromes in 15q11-q13 (OMIM #176270, 105830), are also good candidates for high-resolution whole-genome array CGH since copy imbalance in these regions can cause imprinting disorders.³⁶

Some oligonucleotide array platforms are based on SNP probes. Although initially produced for genotyping to facilitate disease-association studies, these arrays have been successfully used to detect genomic copy number.^{37,38} SNP arrays cannot only detect gains or losses of DNA but can also leverage their genotyping capability to detect uniparental disomy and loss of heterozygosity. Moreover, these arrays could eventually include probes to screen for disease-related point mutations, such as those that cause cystic fibrosis (OMIM #602421), β -thalassemia (OMIM #141900) and a variety of other single-gene disorders.

There are clear advantages to high-resolution whole-genome arrays, but, paradoxically, the complexity of data generated from these arrays can itself be a major constraint in their use. As the macroresolution of an array increases, more sequence variants can be detected, and the proportion of small variants among the total number of variants detected also increases. This is the single most challenging aspect of using oli-

gonucleotide arrays designed with high-density probes at regular intervals. Therefore, it is important to keep background noise to a minimum so that the data from an array are highly reproducible and accurate and can identify all potential alterations. This can be addressed by using dye-reversal hybridization, a process in which two arrays (non-SNP arrays only) are separately hybridized with differentially labeled patient and control DNA samples. The two samples hybridized to the first array are identical to the two hybridized to the second array, but the labeling patterns are reversed (e.g., patient-Cy5/control-Cy3 on the first array and patient-Cy3/control-Cy5 on the second array). Data from the two arrays are then superimposed to ensure that results are consistent. This approach is especially necessary as the size of the smallest observed sequence variant in a given region approaches the microresolution of that region. In this regard, high-density probe coverage can be very useful because multiple adjacent probes would likely detect a single-sequence variant, thereby improving confidence in that observation.

Two difficulties arise due to the high rate of detecting sequence variants on high-resolution whole-genome arrays. First, the smaller a sequence variant, the more difficult it is to find reagents (e.g., BAC clones) to verify that variant by FISH. If the variant is smaller than the length of an average BAC, then molecular methods such as qPCR or multiplex ligation-dependent amplification must be used. However, these methods cannot distinguish between different types of chromosome rearrangements, each of which can have very different recurrence risks.

Another serious issue related to the detection of sequence variants on high-resolution arrays is recognizing pathogenic sequence variants from benign ones. It has become evident that the human genome contains a much greater degree of normal variation than previously expected in sequence lengths between 1 kb and 3 Mb. Accounting for as much as 5% of the total human genome, these large-scale variants are referred to as copy number variants (CNVs) to distinguish them from smaller variants, such as SNPs, short tandem repeats, short interspersed repeats, and long interspersed repeats.^{37,39–41} CNVs are generally considered benign in terms of disease association, and some are known to be polymorphic (i.e., present in at least 1% of the population). In the context of array CGH, more CNVs will be detected as the macroresolution increases, thereby making data interpretation more complicated.

Ongoing efforts to catalog human CNVs have been successful, but the record is still incomplete.^{7,37,41} Two public databases that are attempting to capture information on all human CNVs are proving to be particularly valuable for clarifying array CGH data.^{42,43} Within the next 3 to 5 years, these databases will become more comprehensive and thus of greater assistance in interpreting array CGH data in a diagnostic setting. Hence, as a first step in analyzing array CGH data generated from high-resolution whole-genome array, potential chromosomal alterations should be matched against the CNV databases to exclude benign variants. Any remaining alterations would represent potentially pathogenic changes that can be verified by comparison with array CGH data or follow-up

FISH data from parental samples. If the parents do not carry the same rearrangement as their child (and nonpaternity is not an issue), clinicians can be a little more confident that the de novo rearrangement in the child is pathogenic. Finally, considering the gene content in that region and surveying the literature for similar alterations (either pathogenic or benign) will further clarify the results from a proband sample. Eventually, the CNV databases, along with other genomics databases, may integrate into the Online Mendelian Inheritance in Man (OMIM) database to provide quick and clear explanations for all observed cytogenetic alterations.

Another approach to aid the interpretation of high-resolution array CGH data are to use stringent algorithms to evaluate the data. Algorithms that measure signals from oligonucleotide arrays are versatile in that they can be designed to evaluate signals only from regions of interest. Hence, regions known to contain benign CNVs can be excluded from analysis; alternatively, until all CNVs are carefully studied, they can be flagged and given scores for pathogenicity based on gene content, repeat content, and correlations with information in CNV databases. In addition, it is possible to survey only regions of known clinical significance, hence allowing a whole-genome array to perform more like a targeted array. If necessary, at an added service cost, the initial masked results can later be unmasked to survey the rest of the genome for potential rearrangements. Notwithstanding potential medical and/or ethical dilemmas, this approach would be comparable with routine versus high-resolution karyotyping and may be a logical first step as high-resolution whole-genome array CGH becomes part of clinical diagnostics.

AN IDEAL ARRAY?

Most cytogenetics laboratories will likely incorporate array CGH into their testing services in the next few years. Many may initially opt for targeted arrays because they can adequately recognize known microdeletions and microduplications and detect unbalanced subtelomeric alterations. Importantly, prenatal testing will rely exclusively on targeted arrays, at least in the foreseeable future, because data interpretation and FISH confirmation are straightforward. However, most facilities will likely recognize that high-resolution whole-genome arrays are more versatile, can perform as targeted arrays with appropriate algorithms if desired, and are essential for analyzing patients with disorders of unknown cause.

As array CGH enters the mainstream and users gain more confidence in whole-genome arrays, an ideal resolution for whole-genome analysis needs to be defined. A moderate resolution of 1 Mb has proven insufficient for detecting many pathogenic anomalies using a whole-genome BAC array.²¹ Clearly, greater macroresolution is needed, but the optimal resolution is not easy to define and will vary for different regions of the genome. This has to be empirically determined. Given the outcomes of recent array CGH studies, it is clear that a macroresolution <100 kb is essential because novel disorders have been described with rearrangements as small as 35 to 100 kb.^{21,44} Eventually, the American College of Medical Genetics

should issue guidelines for the minimum components of a “consensus” array that will bring some agreement on the content of a diagnostic genomic array.

We predict that the general trend will be toward carefully designed, high-resolution whole-genome (i.e., high-density oligonucleotide) arrays with high overall macroresolution (perhaps 50–100 kb) and higher microresolution in gene-rich intervals and regions of known medical significance. However, a couple of issues might slow this trend. First, high-resolution whole-genome arrays are relatively expensive, although the costs are expected to decrease as the technology gains wider acceptance. Second, different users may use different array CGH platforms and even different designs within a single platform, making comparison of findings on two separate arrays difficult. However, enforcing nomenclature standards that include information on sequence coordinates for all genomic alterations reported can alleviate this problem.⁴⁵

A third and more pressing hurdle in the increased use of high-resolution whole-genome arrays is the difficulty of confirming de novo alterations that they detect. For this reason, as explained earlier, parental DNA samples must be analyzed when a potentially new alteration is found. Typically, novel findings on an array should be followed up by FISH to distinguish between chromosomal imbalances, such as simple deletions and unbalanced translocations or simple duplications, and marker chromosomes. In cases of novel rearrangements in complex genomic regions, such as 15q11,^{21,46,47} parental samples may have to be examined by array CGH as well. Where FISH proves unreliable because the deletion interval is too small or rich with repeats, follow-up would require molecular methods such as qPCR. This prolonged follow-up could cause considerable anxiety not only for patients and their families but also for laboratory directors, clinicians, and genetic counselors involved in the cases. To address these circumstances, we recommend an initiative within the cytogenetics community to build a library of genome-wide reagents that are quickly accessible for confirming array CGH results. These reagents could include FISH-tested BAC clones or pooled oligonucleotides, multiplex ligation-dependent amplification probes, and isothermal qPCR primers. Ultimately, this targeted versus whole-genome array issue will not exist because all clinically relevant cytogenetic imbalances discovered on whole-genome arrays will eventually be incorporated into the testing menu of a comprehensive targeted array that will have representative confirmation probes for every locus and will also exclude regions of known benign variation.

Cytogenetics will soon witness a seismic shift in testing paradigm as high-resolution whole-genome array CGH is introduced into a growing number of laboratories. This will tremendously benefit the field of medical genetics as novel genomic disorders are discovered and subtleties in both new and previously characterized disorders are revealed. With appropriate regulatory oversight and careful integration into cytogenetics laboratories, this sophisticated technology should help diagnose many more patients with genetic conditions than previously thought possible.

ACKNOWLEDGMENTS

The authors thank Dr. Melanie A. Manning for critically reviewing and Kerry W. Aradhya for editing this manuscript.

References

- Trask BJ. Human cytogenetics: 46 chromosomes, 46 years and counting. *Nat Rev Genet* 2002;3:769–778.
- Speicher MR, Carter NP. The new cytogenetics: blurring the boundaries with molecular biology. *Nat Rev Genet* 2005;6:782–792.
- Cheung SW, Shaw CA, Yu W, Li J, et al. Development and validation of a CGH microarray for clinical cytogenetic diagnosis. *Genet Med* 2005;7:422–432.
- Shaffer LG, Kashork CD, Saleki R, Rorem E, et al. Targeted genomic microarray analysis for identification of chromosome abnormalities in 1500 consecutive clinical cases. *J Pediatr* 2006;149:98–102.
- Shaw-Smith C, Redon R, Rickman L, Rio M, et al. Microarray based comparative genomic hybridisation (array-CGH) detects submicroscopic chromosomal deletions and duplications in patients with learning disability/mental retardation and dysmorphic features. *J Med Genet* 2004;41:241–248.
- Tyson C, McGillivray B, Chijiwa C, Rajcan-Separovic E. Elucidation of a cryptic interstitial 7q31.3 deletion in a patient with a language disorder and mild mental retardation by array-CGH. *Am J Med Genet A* 2004;129:254–260.
- Iafate AJ, Feuk L, Rivera MN, Listewnik ML, et al. Detection of large-scale variation in the human genome. *Nat Genet* 2004;36:949–951.
- Erdogan F, Ullmann R, Chen W, Schubert M, et al. Characterization of a 5.3 Mb deletion in 15q14 by comparative genomic hybridization using a whole genome “tiling path” BAC array in a girl with heart defect, cleft palate, and developmental delay. *Am J Med Genet A* 2007;143:172–178.
- Fiegler H, Redon R, Andrews D, Scott C, et al. Accurate and reliable high-throughput detection of copy number variation in the human genome. *Genome Res* 2006;16:1566–1574.
- Tzschach A, Menzel C, Erdogan F, Schubert M, et al. Characterization of a 16 Mb interstitial chromosome 7q21 deletion by tiling path array CGH. *Am J Med Genet A* 2007;143:333–337.
- Gershon D. DNA microarrays: more than gene expression. *Nature* 2005;437:1195–1198.
- Gribble SM, Prigmore E, Burford DC, Porter KM, et al. The complex nature of constitutional de novo apparently balanced translocations in patients presenting with abnormal phenotypes. *J Med Genet* 2005;42:8–16.
- Shaffer LG. American College of Medical Genetics guideline on the cytogenetic evaluation of the individual with developmental delay or mental retardation. *Genet Med* 2005;7:650–654.
- Ravnán JB, Tepperberg JH, Papenhausen P, Lamb AN, et al. Subtelomere FISH analysis of 11 688 cases: an evaluation of the frequency and pattern of subtelomere rearrangements in individuals with developmental disabilities. *J Med Genet* 2006;43:478–489.
- Bejjani BA, Saleki R, Ballif BC, Rorem EA, et al. Use of targeted array-based CGH for the clinical diagnosis of chromosomal imbalance: is less more? *Am J Med Genet A* 2005;134:259–267.
- Heilstedt HA, Ballif BC, Howard LA, Lewis RA, et al. Physical map of 1p36, placement of breakpoints in monosomy 1p36, and clinical characterization of the syndrome. *Am J Hum Genet* 2003;72:1200–1212.
- Yu W, Ballif BC, Kashork CD, Heilstedt HA, et al. Development of a comparative genomic hybridization microarray and demonstration of its utility with 25 well-characterized 1p36 deletions. *Hum Mol Genet* 2003;12:2145–2152.
- Ballif BC, Rorem EA, Sundin K, Lincicum M, et al. Detection of low-level mosaicism by array CGH in routine diagnostic specimens. *Am J Med Genet A* 2006;140:2757–2767.
- Yu S, Barbooth D, Benke PJ, Warburton PE, et al. Characterization of a neocentric supernumerary marker chromosome originating from the Xp distal region by FISH, CENP-C staining, and array CGH. *Cytogenet Genome Res* 2007;116:141–145.
- Vermeesch JR, Melotte C, Salden I, Riegel M, et al. Tetrasomy 12pter-12p13.31 in a girl with partial Pallister-Killian syndrome phenotype. *Eur J Med Genet* 2005;48:319–327.
- Aradhya SManningMA CherryAM. High-resolution whole-genome aCGH identifies novel contiguous gene deletions and duplications associated with developmental delay, mental retardation, and dysmorphic features. *Am J Med Genet A* 2007;143:1431–1441.
- Miyake N, Shimokawa O, Harada N, Sosonkina N, et al. BAC array CGH reveals genomic aberrations in idiopathic mental retardation. *Am J Med Genet A* 2006;140:205–211.
- Schoumans J, Ruivenkamp C, Holmberg E, Kyllerman M, et al. Detection of chromosomal imbalances in children with idiopathic mental retardation by array based comparative genomic hybridisation (array-CGH). *J Med Genet* 2005;42:699–705.
- Redon R, Baujat G, Sanlaville D, Le Merrer M, et al. Interstitial 9q22.3 microdeletion: clinical and molecular characterisation of a newly recognised overgrowth syndrome. *Eur J Hum Genet* 2006;14:759–767.
- Koolen DA, Vissers LE, Pfundt R, de Leeuw N, et al. A new chromosome 17q21.31 microdeletion syndrome associated with a common inversion polymorphism. *Nat Genet* 2006;38:999–1001.
- Shaw-Smith C, Pittman AM, Willatt L, Martin H, et al. Microdeletion encompassing MAPT at chromosome 17q21.3 is associated with developmental delay and learning disability. *Nat Genet* 2006;38:1032–1037.
- Sharp AJ, Hansen S, Selzer RR, Cheng Z, et al. Discovery of previously unidentified genomic disorders from the duplication architecture of the human genome. *Nat Genet* 2006;38:1038–1042.
- Shieh JT, Aradhya S, Novelli A, Manning MA, et al. Nablus mask-like facial syndrome is caused by a microdeletion of 8q detected by array-based comparative genomic hybridization. *Am J Med Genet A* 2006;140:1267–1273.
- Johnston JJ, Walker RL, Davis S, Facio F, et al. Zoom-in comparative genomic hybridisation arrays for the characterisation of variable breakpoint contiguous gene syndromes. *J Med Genet* 2007;44:e59.
- Shaw CJ, Lupski JR. Implications of human genome architecture for rearrangement-based disorders: the genomic basis of disease. *Hum Mol Genet* 2004;13 Spec No 1:R57–64.
- Lennon PA, Cooper ML, Curtis MA, Lim C, et al. Array-based comparative genomic hybridization facilitates identification of breakpoints of a novel der(1)(t(1;18)(p36.3;q23)dn in a child presenting with mental retardation. *Am J Med Genet A* 2006;140:1156–1163.
- Hayashi S, Kurosawa K, Imoto I, Mizutani S, et al. Detection of cryptic chromosome aberrations in a patient with a balanced t(1;9)(p34.2;p24) by array-based comparative genomic hybridization. *Am J Med Genet A* 2005;139:32–36.
- Gunn SR, Mohammed M, Reveles XT, Viskochil DH, et al. Molecular characterization of a patient with central nervous system dysmyelination and cryptic unbalanced translocation between chromosomes 4q and 18q. *Am J Med Genet A* 2003;120:127–135.
- Hwang KS, Pearson MA, Stankiewicz P, Lennon PA, et al. Cryptic unbalanced translocation t(17;18)(p13.2;q22.3) identified by subtelomeric FISH and defined by array-based comparative genomic hybridization in a patient with mental retardation and dysmorphic features. *Am J Med Genet A* 2005;137:88–93.
- Dhumi P, Coffey AJ, Abbs S, Vermeesch JR, et al. Exon array CGH: detection of copy-number changes at the resolution of individual exons in the human genome. *Am J Hum Genet* 2005;76:750–762.
- Buiting K, Gross S, Lich C, Gillissen-Kaesbach G, et al. Epimutations in Prader-Willi and Angelman syndromes: a molecular study of 136 patients with an imprinting defect. *Am J Hum Genet* 2003;72:571–577.
- Redon R, Ishikawa S, Fitch KR, Feuk L, et al. Global variation in copy number in the human genome. *Nature* 2006;444:444–454.
- Peiffer DA, Le JM, Steemers FJ, Chang W, et al. High-resolution genomic profiling of chromosomal aberrations using Infinium whole-genome genotyping. *Genome Res* 2006;16:1136–1148.
- Feuk L, Carson AR, Scherer SW. Structural variation in the human genome. *Nat Rev Genet* 2006;7:85–97.
- Freeman JL, Perry GH, Feuk L, Redon R, et al. Copy number variation: new insights in genome diversity. *Genome Res* 2006;16:949–961.
- Wong KK, deLeeuw RJ, Dosanjh NS, Kimm LR, et al. A comprehensive analysis of common copy-number variations in the human genome. *Am J Hum Genet* 2007;80:91–104.
- Database of Genomic Variants [database online]. Toronto, Canada: The Centre for Applied Genomics. Available at: <http://projects.tcag.ca/variation/>. Updated March 19, 2007.
- Database of Chromosomal Imbalance and Phenotype in Humans using Ensembl Resources (DECIPHER) [database online]. Cambridge, UK: The Sanger Institute. Available at: <http://www.sanger.ac.uk/PostGenomics/decipher/>. Updated August 10, 2007.
- Friedman JM, Baross A, Delaney AD, Ally A, et al. Oligonucleotide microarray analysis of genomic imbalance in children with mental retardation. *Am J Hum Genet* 2006;79:500–513.
- Shaffer LGTommerupN, editors. ISCN 2005. Basel, Switzerland: Karger, 2005.
- Chai JH, Locke DP, Grealley JM, Knoll JH, et al. Identification of four highly conserved genes between breakpoint hotspots BP1 and BP2 of the Prader-Willi/Angelman syndromes deletion region that have undergone evolutionary transposition mediated by flanking duplicons. *Am J Hum Genet* 2003;73:898–925.
- Fantes JA, Mewborn SK, Lese CM, Hedrick J, et al. Organisation of the pericentromeric region of chromosome 15: at least four partial gene copies are amplified in patients with a proximal duplication of 15q. *J Med Genet* 2002;39:170–177.