



OPEN

Direct PCR amplification from saliva sample using non-direct multiplex STR kits for forensic DNA typing

Pankaj Shrivastava¹✉, Toshi Jain¹ & R. K. Kumawat²

Due to its proficiency to provide the most discriminating results for forensic applications, medical research and anthropological studies, multiplex PCR based STR analysis has been established as the most efficient technique in the forensic DNA analysis. Several multiplex amplification kits based on 4, 5 and 6 dyes chemistry are commercially available and used in forensic DNA typing across the globe. These multiplex PCR systems are routinely used for amplification of multiple STR loci (Autosomal, Y and/or X STR's) in the DNA extracted from various biological samples. In the routine forensic DNA testing, DNA profile obtained is compared with the DNA profile of the reference sample, which takes a certain turnaround time and employs costly lab resources. Successive development in forensic DNA typing have resulted in advent of improved multiplex kits which have reduced the effective analysis time, cost and minimized the number of steps required in comparison to conventional forensic DNA typing. Specialized direct amplification compatible multiplex kits are also available nowadays. These kits are relatively costlier but still require few pre-processing steps, which does not make them worth the hefty cost. Herein, this study, we have used non-direct multiplex STR kits to assess their efficacy for direct amplification. In the present study, 103 saliva samples were directly amplified without any pre-treatment of the samples using thirteen non-direct multiplex kits (4 dyes, 5 dyes and 6 dyes chemistry based) for forensic DNA typing. Here, we report a validated direct PCR amplification protocol from the reference saliva samples by omitting DNA extraction and quantification steps, which resulted in 80% reduction of the turnaround time. The developed protocol is cost effective, time efficient and it does not compromise with the quality of DNA profiles. To the best of our knowledge, this is the first report for direct amplification of DNA with the most commonly used non-direct multiplex STR kits without any pre-treatment of the sample. Complete DNA profiles matching all the essential quality parameters were obtained successfully from all the tested samples.

In early 1990s, the advent of the polymerase chain reaction (PCR) technology and the use of short tandem repeat (STR) polymorphism¹ was a major breakthrough in forensic DNA technology. In last two decades, PCR-based STR typing has become a routine technique in the forensic investigations due to associated features such as specificity, sensitivity, feasibility to simultaneously amplify several loci (multiplexing) and automation^{2,3}. Despite several improvements in the methodology, time required for the conventional forensic STR analysis could not be reduced and the process is still lengthy needing hours to several days for processing of most of the samples⁴. Conventional forensic DNA typing is a multistep process involving several steps of DNA extraction, quantification, amplification, genotyping and analysis (Fig. 1). These steps are not just time consuming but also result in loss of the valuable DNA^{5,6}, which is already in minute quantities in the forensic samples. Also multiple processing steps need more human intervention which further increases the possibility of contamination and error. If these steps are reduced and the amplification could be achieved without the extraction, quantification, and concentration processes, higher yield of DNA could be achieved. Also, the chances of manual error and contamination could be reduced. Approximately 20–76% of DNA is lost from the swab samples during the DNA extraction step⁷. Various studies have reported that column-based DNA extraction techniques also result in the loss of DNA, thereby affecting the genotyping using this method^{8–14}. Most of the time in routine forensic DNA typing is consumed in the pre-processing steps of DNA extraction and Quantitation⁴. The development of semi-automated DNA

¹DNA Fingerprinting Unit, State Forensic Science Laboratory, Sagar, MP 470001, India. ²DNA Division, State Forensic Science Laboratory, Rajasthan, Jaipur 302016, India. ✉email: pankaj.shrivastava@rediffmail.com




Figure 1. Workflow of conventional forensic DNA typing.

extraction kits and automated DNA extraction systems using pre-formulated kits has been a major breakthrough in this area so far. Some short non-automated extraction protocols like Chelex or alkaline lysis method and/or magnetic, paramagnetic beads based rapid DNA isolation are also in use. Since last decade, most of the forensic laboratories have encountered a considerable rise in the number of forensic DNA casework^{15–17}. Towards the effort to avoid the delay in justice, a lot of research work has been conducted to improve the analysis speed and to develop faster processing methods to save the reporting time. Development of rapid, compact and portable devices capable of producing real-time results from forensic samples at the crime scene has been emphasised upon^{15–17}. However, these efforts were not of much help to reduce the cost and time efficiency of the forensic DNA analysis, therefore were not adopted in most of the laboratories¹⁸. The analysis of mandatory reference samples in the forensic cases, which may vary from one to many, also require a similar lengthy processing time, as that of the case samples. This drawback is associated with both the methods viz. automated or specific direct amplification compatible chemistries¹⁹ which adds to the further delay in reporting the case. To curb down the analysis time and fasten the process of case reporting, there has been a growing interest towards direct amplification of case and reference samples such as blood and buccal swabs, single hair follicle and tape lifts of clothing, swab on FTA card^{20–24}, blood²⁵, fabrics⁶, hairs^{26,27}, touch DNA^{28,29}, blood stain³⁰, fingernails³¹, tissues^{32,33} and fibers³⁴. Direct amplification protocol eliminates the DNA extraction and quantitation steps, and accommodates the sample directly to the PCR step. The streamlined process of direct amplification for processing the DNA samples has obvious benefits such as simplify the process, low risk of contamination due to less handling steps, cost efficiency and has low chances of loss of forensic DNA, which already is in minute quantities^{32,33}. Methods to successfully genotype different samples from varied origins in a single shot by direct amplification sounds promising and challenging too. Forensic cases have various sorts of circumstantial exhibits which could include anything found at the scene of crime. One such forensic case exhibits are buccal swabs and in the case, the subjects had a habit of smoking, chewing tobacco etc., it could cause inhibition in the process of PCR. Although, the buccal cells are stable for weeks or months if stored in the laboratory at $-20\text{ }^{\circ}\text{C}$ till final processing, they are quite vulnerable to the bacterial growth, nucleic acid degradation. Reduced amplification ability after a span of 4 days or so has been observed, particularly under warm and moist conditions³⁵. Keeping these points into consideration, saliva sample was chosen for the direct analysis. The fact that saliva or any cell suspension is expected to have some cell free DNA formed the ground for this study. Also, another important aspect considered was the nature and texture of the saliva cells as these cells are not strong enough to stand high temperatures like $95\text{ }^{\circ}\text{C}$, which is a prerequisite in most of the amplification kits for activation of hot start Taq polymerase³⁶.

The process of genotyping nowadays does not require the high yield of DNA. This could be attributed to the improvement in the amplification step by using the upgraded and highly sensitive STR kits which require very low inputs of DNA (approx less than 500pg). The advantages associated with these advanced kits not only include the increase in the number of loci to be examined, as per the recommendation of the Scientific Working Group on DNA Analysis Methods (SWGDM)³⁷, the master mix has also been improved to reduce the amplification time and to cope with the PCR inhibitors. Rapid amplification of multiplex STRs for human identification has also been demonstrated by changing the master mix and/or enzyme^{22,23,38,39}. There are certain multiplex kits specially designed for direct amplification of single source samples for the population data purpose. However, still most of the commonly used multiplex STR amplification kits need the pre-processing of samples before DNA amplification and genotyping⁴⁰. Also, these specially designed kits vary in their markers. Therefore, if the standard kits are used for direct and rapid amplification, resultant profiles will not be comparable for all the tested markers with the profiles from routine casework samples and also from the existing DNA databases. Most of the widely used commercially available kits (non-direct) are not meant or reported to amplify the samples without DNA extraction. Working on a research project on X-STRs, we collected saliva samples of unrelated individuals and used Investigator Argus X-12 multiplex STR kit (Qiagen, Germany) for DNA amplification as

per manufacturer's recommendation. DNA was isolated and extracted from 188 saliva samples using automated DNA extraction system 12 GC (Precision System Science Co., Ltd., Matsudo, Japan) following manufacturer's instructions and pre formulated kits⁴¹. Few samples were tried for direct amplification of saliva samples with the same multiplex kit, using the same protocol for amplification. Promising results in terms of good quality DNA profile from the pilot study resulted into this detailed study on the direct amplification of saliva samples using commercially available non-direct multiplex kits. Kits used in the study were from leading brands including Thermo Fisher Scientific (AMPFLSTR IDENTIFILER, AMPFLSTR IDENTIFILER PLUS, GLOBALFILER, AMPFLSTR Y FILER and AMPFLSTR YFILER PLUS), Promega (POWERPLEX 16HS SYSTEM, POWERPLEX 21 SYSTEM, POWERPLEX FUSION 6C SYSTEM and POWERPLEX Y 23 SYSTEM) and Qiagen (INVESTIGATOR IDPLEX PLUS, and INVESTIGATOR ARGUS X-12 MULTIPLEX KIT) were tested for direct amplification protocol. This novel protocol was also tested using VERIFILER PLUS PCR AMPLIFICATION KIT (Thermo Fisher Scientific) and SURE ID PANGLOBAL HUMAN DNA IDENTIFICATION KIT (Health Gene Technologies) as well newly launched but still not used in routine forensic DNA typing kits. The details of multiplex kits used in this study are mentioned in Table 1. This study was designed with the aim to test and validate the direct amplification protocol using commercially available 4, 5 and 6 dyes chemistry based non-direct multiplex kits of leading brands in forensics.

Result and discussion

Optimizing parameters for direct PCR amplification. For the pilot study, initially 10 saliva samples were directly amplified using both the recommended full reaction volume of 25 μ l and a reaction mixture with reduced volume of 10 μ l with the help of PCR using all the tested multiplex kits. The DNA profile obtained from both the reaction volumes was evaluated, considering amplification of all the markers to showcase a full DNA profile, peak heights including inter and intra marker balance, observation of stutter and/or other artifacts. In this preliminary test, all the tested samples produced quality DNA profiles with both the reaction volumes and showed concordance. Our few previous studies have also reported use of 10 μ l PCR reaction volume for an efficient reaction product, which further supported the reduction of final PCR reaction volume^{42,43}. Direct amplification protocol was further validated for 10 μ l reaction volume PCR reaction as per recommended PCR conditions of the particular multiplex kits. Using the above said direct amplification protocol, all the tested samples were directly amplified using non-direct 4 dyes chemistry based autosomal STR marker multiplex kit POWERPLEX 16HS SYSTEM (now not commonly used in forensics), 5 dyes chemistry based autosomal STR marker multiplex kits AMPFLSTR IDENTIFILER, AMPFLSTR IDENTIFILER PLUS, POWERPLEX 21 SYSTEM and INVESTIGATOR IDPLEX PLUS, 6 dyes chemistry based autosomal STR marker multiplex kits GLOBALFILER, POWERPLEX FUSION 6C SYSTEM, VERIFILER PLUS and SURE ID PANGLOBAL HUMAN DNA IDENTIFICATION KIT; Y-STR kits AMPFLSTR YFILER, POWERPLEX Y 23 SYSTEM and AMPFLSTR YFILER PLUS MULTIPLEX KITS; X-STR kit INVESTIGATOR ARGUS X-12 MULTIPLEX KIT. For all the used kits, standard recommended conditions except final reaction volume were followed. Complete and quality DNA profile could be obtained with all the tested multiplex kits using this protocol.

All the tested samples were directly amplified at all the loci with the all the above mentioned multiplex PCR kits used in this study. However, it should be noted that the performance evaluation by the peak heights between the kits is only an approximation. It is because the same loci are labeled with different fluorescent dyes and have different amplicon sizes in different kits²² and for the proper assessment of DNA profile, the overall quality of the profile using profile quality measures is more important⁴⁴.

In conventional PCR method, Tris EDTA (TE) buffer is used for the final volume adjustment, which has also been reported to inhibit PCR amplification to some extent^{19,45}. In this study, amplification grade water was used for final volume adjustment, which reduced the chances of reaction inhibition due to TE. In the case of DNA degradation or inhibition, larger/high molecular weight loci (usually more than 300 bp loci) reflect low peak height, higher peak imbalance or allele drop outs. This effect leads to development of a DNA profile with a slope peak height, also termed as SKI slope effect^{19,46}. Noticeably, in our study almost all the DNA profiles obtained were having balanced peak heights without SKI slope effect. Rarely, very few DNA profiles were observed with the SKI slope effect which was due to the degradation of the samples. This was confirmed by the usual automated extraction followed by RT PCR of those few samples (data not shown).

Conventional DNA typing versus direct amplification methods. To determine the difference between the conventional DNA typing (including DNA extraction, quality and quantity check, amplification and Genotyping), and the present method of direct amplification, and to estimate the DNA quantity likely to be present in the direct amplification reaction, 10 saliva samples were used for amplification with both the methods. The mean quantity of DNA was found to be 0.45 ng for all the 10 tested samples and most of the multiplex kits claim to be more sensitive for more than 0.25 ng input DNA^{47,48}. With conventional and direct amplification methods, all the alleles were observed at all the loci of respective multiplex kits and balanced high quality DNA profiles were obtained, thus establishing concordance. Further, all the tested samples of this study were directly amplified using commonly used non-direct multiplex kits used in forensic DNA typing (Table 1). Direct PCR amplification protocol from the reference saliva samples by omitting DNA extraction and quantification steps, which resulted in 80% reduction of the turnaround time (Table 2).

Assessment of DNA profile quality. To assure the quality of DNA profile, Total Peak Height (TPH), Peak Height Ratio (PHR), inter and intra locus balance (local and global balance) of the obtained DNA profiles were analyzed.

PCR multiplex kit	Make/manufacturer of the multiplex kit	Dye set used in the kit	Launching year	Turnaround time for PCR amplification as per recommended protocol by the manufacturer	STR markers included in the multiplex kit							Total number of markers	
					Markers in the multiplex	Autosomal STRs	Sex determination	SNP	Y-STRs	X-STRs	Other		
AMPFISTR IDENTIFILER PCR AMPLIFICATION KIT	Thermo Fisher Scientific	5 Dye	2001	Approx 150 min	D8S1179, D21S11, D7S820, CSF1PO, D3S1358, TH01, D13S317, D16S539, D2S1338, D19S433, vWA, TPOX, D18S51, Amelogenin, D5S818 and FGA	15	Amelogenin	-	-	-	-	-	16
AMPFISTR IDENTIFILER PLUS PCR AMPLIFICATION KIT	Thermo Fisher Scientific	5 Dye	2010	Approx 150 min	D8S1179, D21S11, D7S820, CSF1PO, D3S1358, TH01, D13S317, D16S539, D2S1338, D19S433, vWA, TPOX, D18S51, Amelogenin, D5S818 and FGA	15	Amelogenin	-	-	-	-	-	16
POWERPLEX 16HS SYSTEM	Promega	4 Dye	2009	Approx 60 min	D18S51, D21S11, TH01, D3S1358, Penta E, FGA, TPOX, D8S1179, vWA, CSF1PO, D16S539, D7S820, D13S317, D5S818, Penta D and Amelogenin	15	Amelogenin	-	-	-	-	-	16
POWERPLEX 21 SYSTEM	Promega	5 Dye	2012	Approx 90 min	D1S1656, D2S1338, D3S1358, D5S818, D6S1043, D7S820, D8S1179, D12S391, D13S317, D16S539, D18S51, D19S433, D21S11, Amelogenin, CSF1PO, FGA, Penta D, Penta E, TH01, TPOX and vWA	20	Amelogenin	-	-	-	-	-	21
INVESTIGATOR IDPLEX PLUS AMPLIFICATION KIT	Qiagen	5 Dye	2014	Approx 90 min	D3S1358, TH01, D21S11, D18S51, D5S818, D13S317, D7S820, D16S539, CSF1PO, vWA, D8S1179, TPOX, FGA, D2S1338, D19S433 and Amelogenin	15	Amelogenin	-	-	-	-	-	16
INVESTIGATOR ARGUS X-12 KIT	Qiagen	5 Dye	2015	Approx 60 min	DXS10103, DXS8378, DXS7132, DXS10134, DXS10074, DXS10101, DXS10135, DXS7423, DXS10146, DXS10079, HPRTB and DXS10148	-	Amelogenin	-	-	12	-	-	12
GLOBALFILER PCR AMPLIFICATION KIT	Thermo Fisher Scientific	6 Dye	2013	Approx 90 min	D3S1358, vWA, D16S539, CSF1PO, TPOX, D8S1179, D21S11, D18S51, D2S441, D19S433, TH01, FGA, D22S1045, D5S818, D13S317, D7S820, SE33, D10S1248, D1S1656, D12S391, D2S1338, DYS391, Y indel and Amelogenin	21	Amelogenin	Y Indel	1	-	-	-	24
Continued													

PCR multiplex kit	Make/manufacturer of the multiplex kit	Dye set used in the kit	Launching year	Turnaround time for PCR amplification as per recommended protocol by the manufacturer	STR markers included in the multiplex kit							Total number of markers
					Markers in the multiplex	Autosomal STRs	Sex determination	SNP	Y-STRs	X-STRs	Other	
SURE ID PANGLOBAL HUMAN DNA IDENTIFICATION KIT	Health Gene Technologies	6 Dye	2018	Approx 60 min	D3S1358, TH01, D21S11, D18S51, Penta E, D12S391, D6S1043, D2S1338, D1S1656, D5S818, D13S317, D7S820, D19S433, CSF1PO, Penta D, D2S441, vWA, D8S1179, TPOX, FGA, D16S539, D22S1045, SE33, D10S1248, Alelogenin, Y Indel and DYS391	24	Amelogenin	Y Indel	1	-	-	27
POWERPLEX FUSION 6C SYSTEM	Promega	6 Dye	2015	Approx 90 min	CSF1PO, FGA, TH01, vWA, D1S1656, D2S1338, D2S441, D3S1358, D5S818, D7S820, D8S1179, D10S1248, D12S391, D13S317, D16S539, D18S51, D19S433, D21S11, Amelogenin, DYS391, Penta D, Penta E, D22S1045, TPOX, SE33, DYS570 and DYS576	23	Amelogenin	-	3	-	-	27
VERIFILER PLUS PCR AMPLIFICATION KIT	Thermo Fisher Scientific	6 Dye	2018	Approx 90 min	D3S1358, vWA, D16S539, CSF1PO, D6S1043, D8S1179, D21S11, D18S51, D5S818, D2S441, D19S433, FGA, D10S1248, D22S1045, D1S1656, D13S317, D7S820, Penta E, Penta D, TH01, D12S391, D2S1338, TPOX, IQCS and IQCL, Y indel and Amelogenin	23	Amelogenin	Y Indel	-	-	2 Quality sensor Markers	27
POWERPLEX Y 23 SYSTEM	Promega	5 Dye	2012	Approx 90 min	DYS576, DYS389I, DYS448, DYS389II, DYS19, DYS391, DYS481, DYS549, DYS533, DYS438, DYS437, DYS570, DYS635, DYS390, DYS439, DYS392, DYS643, DYS393, DYS458, DYS385a/b, DYS456 and Y-GATA-H4	-	-	-	23	-	-	23
AMPFLSTR YFILER PCR AMPLIFICATION KIT	Thermo Fisher Scientific	5 Dye	2004	Approx 150 min	DYS456, DYS389 I, DYS390, DYS389 II, DYS458, DYS19, DYS385 a/b, DYS393, DYS391, DYS439, DYS635, DYS392, Y GATA H4, DYS437, DYS438 and DYS448	-	-	-	17	-	-	17

Continued

PCR multiplex kit	Make/manufacturer of the multiplex kit	Dye set used in the kit	Launching year	Turnaround time for PCR amplification as per recommended protocol by the manufacturer	STR markers included in the multiplex kit							Total number of markers
					Markers in the multiplex	Autosomal STRs	Sex determination	SNP	Y-STRs	X-STRs	Other	
YFILER PLUS PCR AMPLIFICATION KIT	Thermo Fisher Scientific	6 Dye	2014	Approx 90 min	DYS576, DYS389I, DYS635, DYS389II, DYS627, DYS460, DYS458, DYS19, YGATAH4, DYS448, DYS391, DYS456, DYS390, DYS438, DYS392, DYS518, DYS570, DYS437, DYS385, DYS449, DYS393, DYS439, DYS481, DYF387S1 and DYS533	-	-	-	27	-	-	27

Table 1. Comparative list of various multiplex kits.

Turnaround time for DNA isolation				DNA quantitation		Total time required from DNA isolation to DNA quantitation (approx in minutes)		Turnaround time for PCR and genotyping
Manual	Automatic DNA extraction system			Quantitation using quantifiler/quantifiler duo kit	Quatitation using quantifiler trio kit	Manual	Automatic	
Organic extraction Method (phenol chloroform isoamyl alcohol-PCIA)	Pre-processing time	Time on machine	Total time					
Approx 300 min	Approx 30–70 min	Approx 18–30 min	Approx 100 min	Approx 100 min	Approx 50 min	Approx 350/400 min	Approx 200/250 min	Approx 120–200 min
Method in practice		Effective turnaround time of complete DNA Processing up to pre-PCR (including multistep handling time) approx 80%						Around 20%
Proposed method		Not required						Remains same

Table 2. Turnaround time for DNA typing process.

TPH is the sum of peak heights of both the heterozygous alleles in the DNA profile. PHR for heterozygous alleles defines the ratio of lower peak height and higher peak height at the specific marker. This value ranges from 0 to 1, where 1 represents the two identical peaks with equal heights and 0 represents the situation where one of the peaks is not observed which might be because of allele drop out due to PCR inhibition or degraded DNA⁴⁴.

$$PHR = \frac{\text{Height of lower peak}}{\text{Height of higher peak}} \text{---for a heterozygous allele at particular locus}$$

PHR = 1 for true homozygous allele at particular locus

Inter and Intra locus balance was measured and represented in the local balance and global balance respectively, which is mean of TPH for all the tested loci in the multiplex kits.

Overall TPH, PHR and global balance quality parameters of DNA profile were evaluated statistically using one way analysis of variance (ANOVA) non parametric applying Friedman test at 5 percent of significant level and Dunn’s Multiple Comparison test using Prism GraphPad v5 software⁴⁹.

Autosomal STR multiplex kits. *AMPFLSTR IDENTIFILER AND AMPFLSTR IDENTIFILER PLUS PCR AMPLIFICATION KIT.* DNA profiles obtained with these multiplex kits were found to be complete and balanced. The average TPH and PHR ranged from 6553 to 7833 and from 0.875 to 0.984, respectively for both the kits. Out of 103 DNA profiles, only five DNA profiles were observed with ski slope effect. In Dunn’s Multiple Comparison test, at 5 percent significant level, all the DNA profiles showed no significant variation in the global balance (Figs. 2, 3), (Supplementary Tables S1, S2).

POWERPLEX 16HS SYSTEM. DNA profiles obtained with this multiplex kit was also found to be balanced with an average TPH and PHR ranging from 6593 to 7889 and 0.863 to 0.993 respectively. Seven DNA profiles were observed with ski slope effect. In Dunn’s Multiple Comparison test, at 5 percent significant level, all the DNA profiles exhibited no significant variation in the global balance (Fig. 4), (Supplementary Table S3).




Figure 2. Quality assessment of DNA profile by Dunn's multiple comparison test for the AMPFLSTR IDENTIFILER PCR AMPLIFICATION KIT (A) Total peak height (B) Peak height ratio (C) Global balance.

POWERPLEX 21 SYSTEM. DNA profiles obtained with this multiplex kit were found to have average TPH and PHR ranging from 6548 to 7925 and from 0.865 to 0.993 respectively. Out of the total studied DNA profiles, ten DNA profiles were observed to have ski slope effect. In Dunn's Multiple Comparison test, at 5 percent significant level, all the DNA profiles accounted no significant variation in the global balance (Fig. 5), (Supplementary Table S4).

INVESTIGATOR IDPLEX PLUS KIT. DNA profiles obtained with this multiplex kit were found to have an average TPH and PHR ranging from 6593 to 7883 and from 0.866 to 0.993, respectively. Out of the total studied DNA profiles, fourteen DNA profiles were observed with ski slope effect. In Dunn's Multiple Comparison test, at 5 percent significant level, all the DNA profiles exhibited no significant variation in the global balance (Fig. 6), (Supplementary Table S5).

POWERPLEX FUSION 6C SYSTEM. DNA profiles obtained with this multiplex kit were found to have an average TPH and PHR values ranging from 6548 to 7927 and from 0.864 to 0.994, respectively. Out of the total studied DNA profiles, ten DNA profiles were observed with ski slope effect. In Dunn's Multiple Comparison test, at 5 percent significant level, all the DNA profiles accounted no significant variation in the global balance (Fig. 7), (Supplementary Table S6).

GLOBALFILER PCR AMPLIFICATION KIT. DNA profiles obtained with this multiplex kit showed average TPH and PHR values ranging from 6546 to 7924 and from 0.861 to 0.994, respectively. Out of the total studied DNA profiles, ten DNA profiles showed ski slope effect. In Dunn's Multiple Comparison test, at 5 percent significant level, all the DNA profiles exhibited no significant variation in the global balance (Fig. 8), (Supplementary Table S1).

VERIFILER PLUS PCR AMPLIFICATION KIT. DNA profiles obtained with this multiplex kit had an average TPH and PHR values ranging from 6549 to 7928 and from 0.867 to 0.995 respectively. Out of the total studied DNA profiles, nine DNA profiles were with ski slope effect. In Dunn's Multiple Comparison test, at 5 percent significant level, all the DNA profiles exhibited no significant variation in the global balance (Fig. 9), (Supplementary Table S8).




Figure 3. Quality assessment of DNA profile by Dunn's multiple comparison test for the AMPFLSTR IDENTIFILER PLUS PCR AMPLIFICATION KIT (A) Total peak height (B) Peak height ratio (C) Global balance.




Figure 4. Quality assessment of DNA profile by Dunn's multiple comparison test for the POWERPLEX 16HS SYSTEM (A) Total peak height (B) Peak height ratio (C) Global balance.




Figure 5. Quality assessment of DNA profile by Dunn's multiple comparison test for the POWERPLEX 21 SYSTEM (A) Total peak height (B) Peak height ratio (C) Global balance.

SURE ID PANGLOBAL HUMAN DNA IDENTIFICATION KIT. DNA profiles obtained with this multiplex kit were found to have an average TPH and PHR values ranging from 6556 to 7922 and from 0.866 to 0.996 respectively. Out of the total studied DNA profiles, fifteen DNA profiles showed ski slope effect. In Dunn's Multiple Comparison test, at 5 percent significant level, all the DNA profiles did not exhibit any significant variation in the global balance (Fig. 10), (Supplementary Table S9).

Y-STR multiplex kits. All the male samples in this study were directly amplified using Y STR multiplex kits. The haplotype data of these samples were statistically evaluated on TPH and global balance parameters. The obtained profiles matched the standard criterion of profile quality index for most of the tested sample.

POWERPLEX Y 23 SYSTEM. DNA profiles obtained with this multiplex kit showed an average TPH range from 3224 to 4019. Out of the total studied DNA profiles, six DNA profiles were observed with ski slope effect. In Dunn's Multiple Comparison test, at 5 percent significant level, all the DNA profiles accounted no significant variation in the global balance (Fig. 11), (Supplementary Table S10).

AMPFLSTR YFILER PCR AMPLIFICATION KIT. DNA profiles obtained with this multiplex kit were found to have an average TPH values ranging from 3211 to 4015. Out of the total studied DNA profiles, eleven DNA profiles were observed with ski slope effect. In Dunn's Multiple Comparison test, at 5 percent significant level, all the DNA profiles did not show any significant variation in the global balance (Fig. 12), (Supplementary Table S11).

AMPFLSTR YFILER PLUS PCR AMPLIFICATION KIT. DNA profiles obtained with this multiplex kit were found to have an average TPH ranging from 3213 to 4051. Out of the total studied DNA profiles, seven DNA profiles had ski slope effect. In Dunn's Multiple Comparison test, at 5 percent significant level, all the DNA profiles showed no significant variation in the global balance (Fig. 13), (Supplementary Table S12).




Figure 6. Quality assessment of DNA profile by Dunn's multiple comparison test for the INVESTIGATOR IDPLEX PLUS KIT (A) Total peak height (B) Peak height ratio (C) Global balance.

X-STR multiplex kit. Since male sample possesses only one X chromosome, it resulted in single peak at particular loci of X-STR marker, and female samples showed two peaks at particular loci due to the two chromosomes. Thus, the DNA profile obtained from the male samples were evaluated for TPH and global balance and the DNA profile obtained from the female samples were evaluated for TPH, PHR and global balance as well. Here we used the mean value of THP, PHR for this evaluation. Conclusive outcomes have been represented in Fig. 14 and Supplementary Table S13 for both the male and female DNA profiles accordingly.

INVESTIGATOR ARGUS X-12 MULTIPLEX KIT. DNA profiles obtained with this multiplex kit were found to have an average TPH value ranging from 6593 to 7949 (for both the male and female DNA profiles) and the mean PHR ranged from 0.847 to 0.999 (for female DNA profiles). Out of the total studied DNA profiles, twenty three DNA profiles were observed with ski slope effect. In Dunn's Multiple Comparison test, at 5 percent significant level, all the DNA profiles showed no significant variation in the global balance (Fig. 14), (Supplementary Table S13).

Overall data suggests that this novel protocol of direct amplification worked well with all the multiplex kits used in this study. The proposed method will not only save the effective cost, but will also curb down the turnaround time. Moreover, this novel direct amplification protocol will supposedly be useful for the speedy analysis of forensic DNA cases and ultimately will lead to the justice dissemination at the earliest.

Quality control. All the authors have passed proficiency test and quality control exercise for DNA fingerprinting from GITAD, Spain (<http://gitad.ugr.es/principal.htm>). Also internal laboratory control standards were followed and controls provided with multiplex kits were used.




Figure 7. Quality assessment of DNA profile by Dunn's multiple comparison test for the POWERPLEX FUSION 6C SYSTEM (A) Total peak height (B) Peak height ratio (C) Global balance.




Figure 8. Quality assessment of DNA profile by Dunn's multiple comparison test for the GLOBALFILER PCR AMPLIFICATION KIT (A) Total peak height (B) Peak height ratio (C) Global balance.




Figure 9. Quality assessment of DNA profile by Dunn's multiple comparison test for the VERIFILER PLUS PCR AMPLIFICATION KIT (A) Peak height ratio (B) Total peak height (C) Global balance.

Methods

Sample collection. 103 saliva samples from anonymous donors, (50 males and 53 females) who hailed from the Ratlam district of Madhya Pradesh were collected in pre-sterile collection tubes. Written informed consent was taken from the subjects in compliance with the Declaration of Helsinki⁵⁰. The above said population has already been depicted in our earlier publication, in which the authors typed the same samples for 12 XSTRs⁴¹ and 15 autosomal STRs⁵¹. All methods were carried out in accordance with the relevant guidelines and regulations. All the experimental protocols were approved by the Institutional committee of DNA Fingerprinting Unit, State Forensic Science Laboratory, Sagar, MP, India. The samples were brought to the laboratory at 4 °C and were stored at – 20 °C. Saliva samples processed with direct amplification protocol were kept in the laboratory at – 20 °C until final processing. A maximum of one month old samples were processed for the direct amplification.

Polymerase chain reaction (PCR) and capillary electrophoresis (CE) conditions. 103 saliva samples were used directly for PCR using the following multiplex systems to test the efficacy of commonly used multiplex systems viz. AMPFLSTR IDENTIFILER, AMPFLSTR IDENTIFILER PLUS, POWERPLEX 21 SYSTEM, INVESTIGATOR IDPLEX PLUS, GLOBALFILER, POWERPLEX FUSION 6C SYSTEM, VERIFILER PLUS, SURE ID PANGLOBAL HUMAN DNA IDENTIFICATION KIT, AMPFLSTR YFILER, POWERPLEX Y 23 SYSTEM, AMPFLSTR YFILER PLUS and INVESTIGATOR ARGUS X-12 MULTIPLEX KITS available for forensic and population genetic purpose. The modifications in the reaction volume, cycle number, and input saliva samples for the respective kits have been shown in Tables 3 and 4.




Figure 10. Quality assessment of DNA profile by Dunn's multiple comparison test for the SURE ID PANGLOBAL HUMAN DNA IDENTIFICATION KIT (A) Total peak height (B) Peak height ratio (C) Global balance.

Capillary electrophoresis. In all the cases, 0.3 μ l of the total PCR product was run on an Applied Biosystems 3100/3500XL Genetic Analyzer using standard parameters and with the size standards/allelic ladder provided with the respective non-direct multiplex kit and a 10 seconds injection time. The obtained data was analyzed using GeneMapper V3.5/GeneMapper ID version 1.4 softwares (Applied Biosystems, Foster City, CA). All the amplifications were performed using ABI thermal cycler 9700 (Applied Biosystems, Foster City, CA), with the modifications (as mentioned with the respective multiplex kits). Alleles were designated on the basis of number of allele repeats with the help of allelic ladder provided along with the respective multiplex kits. Peak detection threshold was set to 50 RFUs for the allele designation.

Statistical evaluation. To evaluate the quality of DNA profile, the Total Peak Height (TPH), Peak Height Ratio (PHR) and global balance parameters were considered. The overall DNA profile quality of direct amplification method was statistically examined applying ANOVA non parametric by Friedman test at 5 percent of significant level Dunn's Multiple Comparison test using Prism GraphPad v5 software⁴⁹.




Figure 11. Quality assessment of DNA profile by Dunn's multiple comparison test for the POWERPLEX Y 23 SYSTEM (A) Total peak height (B) Global balance.




Figure 12. Quality assessment of DNA profile by Dunn's multiple comparison test for the AMPFLSTR YFILER PCR AMPLIFICATION KIT (A) Total peak height (B) Global balance.




Figure 13. Quality assessment of DNA profile by Dunn's multiple comparison test for the AMPF/STR YFILER PLUS PCR AMPLIFICATION KIT (A) Total peak height (B) Global balance.




Figure 14. Quality assessment of DNA profile by Dunn's multiple comparison test for the INVESTIGATOR ARGUS X-12 MULTIPLEX KIT (A) Total peak height (B) Peak height ratio (C) Global balance.

S. no.	PCR MULTIPLEX KIT	Number of samples used in the study	Final PCR volume used in the study (in μL)	Input volume of saliva sample (in μL)	Modification in PCR cycle number form the recommended by the respective multiplex system	Complete DNA profile obtained	No. of DNA profile with SKI slope	% of balanced DNA profile obtained from direct amplification protocol	Time after sample collection
1	AMPFISTR IDENTIFILER PCR AMPLIFICATION KIT	103	10	1	Two cycle increase	103	5	95	Immediately to more than one month storage of sample
2	AMPFISTR IDENTIFILER PLUS PCR AMPLIFICATION KIT	103	10	1	Two cycle increase	103	5	95	
3	POWERPLEX 16HS SYSTEM	103	10	1	One cycle decrease	103	7	93	
4	POWERPLEX 21 SYSTEM	103	10	0.5	Two cycle decrease	103	10	90	
5	INVESTIGATOR IDPLEX PLUS AMPLIFICATION KIT	103	10	1	Two cycle decrease	103	14	86	
6	INVESTIGATOR ARGUS X-12 KIT	103	10	0.5	One cycle decrease	103	23	78	
7	GLOBALFILER PCR AMPLIFICATION KIT	103	10	1	One cycle decrease	103	10	90	
8	SURE ID PANGLOBAL HUMAN DNA IDENTIFICATION KIT	103	10	1	No change in PCR cycle	103	15	85	
9	POWERPLEX FUSION 6C SYSTEM	103	10	0.5	Two cycle decrease	103	10	90	
10	VERIFILE PLUS PCR AMPLIFICATION KIT	103	10	1	No change in PCR cycle	103	9	91	
11	POWERPLEX Y 23 SYSTEM	50	10	0.5	Two cycle decrease	50	6	88	
12	AMPFLSTR YFILER PCR AMPLIFICATION KIT	50	10	1	One cycle increase	50	11	78	
13	YFILER PLUS PCR AMPLIFICATION KIT	50	10	1	No change in PCR cycle	50	7	86	

Table 3. Validated parameters of novel direct amplification protocol.

Multiplex kit	PCR reaction volume recommended by the manufacturer					Reduced PCR reaction volume for the direct amplification				
	Master mix (in μ L)	Primer set (volume in μ L)	Sample input (containing DNA 500 pg to 1 ng) (in μ L)	Amplification grade water to achieve final reaction volume (in μ L)	Final reaction volume (in μ L)	Master mix (in μ L)	Primer set (in μ L)	Sample input (in μ L)	Amplification grade water to achieve final reaction volume (in μ L)	Final reaction volume (in μ L)
AMPFISTR IDENTIFILER PCR AMPLIFICATION KIT	10.5 + (0.5 AmpliTaq Gold DNA Polymerase)	5.5	1	7.5	25	5 + (0.2 AmpliTaq Gold DNA Polymerase)	2	1	1.8	10
AMPFISTR IDENTIFILER PLUS PCR AMPLIFICATION KIT	10	5	1	9	25	5	2	1	2	10
POWERPLEX 16HS SYSTEM	5	2.5	1	16.5	25	2.25	1.15	1	5.6	10
POWERPLEX 21 SYSTEM	5	5	1	14	25	2	2	0.5	5.5	10
INVESTIGATOR IDPLEX PLUS AMPLIFICATION KIT	7.5	2.5	1	14	25	3	1	1	5	10
INVESTIGATOR ARGUS X-12 KIT	5 + (0.6 Multi Taq2 polymerase)	2.5	1	15.9	25	2 + (0.2 Multi Taq2 polymerase)	1	0.5	5.3	10
GLOBALFILER PCR AMPLIFICATION KIT	7.5	2.5	1	14	25	3	1	1	5	10
SURE ID PANGLOBAL HUMAN DNA IDENTIFICATION KIT	7.5	2.5	1	14	25	3	1	1	5	10
POWERPLEX FUSION 6C SYSTEM	5	5	1	14	25	2	2	0.5	5.5	10
VERIFILE PLUS PCR AMPLIFICATION KIT	5	2.5	1	16.5	25	2	1	1	6	10
POWERPLEX Y 23 SYSTEM	5	2.5	1	16.5	25	2	1	0.5	5.5	10
AMPFLSTR YFILER PCR AMPLIFICATION KIT	9.2 + (0.8 AmpliTaq Gold DNA Polymerase)	5	1	9	25	4 + (0.2 AmpliTaq Gold DNA Polymerase)	1.25	1	5.55	10
YFILER PLUS PCR AMPLIFICATION KIT	10	5	1	9	25	3	1	1	5	10

Table 4. Details of recommended and reduced PCR reaction volume.

Received: 28 July 2016; Accepted: 15 March 2021

Published online: 29 March 2021

References

- Edwards, A. L., Hammond, H. A., Jin, L., Caskey, C. T. & Chakraborty, R. Genetic variation at five trimeric and tetrameric tandem repeat loci in four human population groups. *Genomics* **12**, 241–253 (1992).
- Schoske, R., Vallone, P. M., Ruitberg, C. M. & Butler, J. M. Multiplex PCR design strategy used for the simultaneous amplification of 10 y chromosome short tandem repeat (STR) loci. *Anal. Bioanal. Chem.* **375**, 333–343 (2003).
- Laurin, N. & Frégeau, C. Optimization and validation of a fast amplification protocol for AmpFISTR[®] Profiler Plus[®] for rapid forensic human identification. *Forensic Sci. Int. Genet.* **6**, 47–57 (2012).
- Butler, J. M. *Advanced Topics in Forensic DNA Typing: Methodology*. (Academic Press, 2011).
- Van Oorschot, R. A. H., Ballantyne, K. N. & Mitchell, R. J. Forensic trace DNA: A review. *Investig. Genet.* **1**, 1–17 (2010).
- Linacre, A., Pekarek, V., Swaran, Y. C. & Tobe, S. S. Generation of DNA profiles from fabrics without DNA extraction. *Forensic Sci. Int. Genet.* **4**, 137–141 (2010).
- Van Oorschot, R. A. H. *et al.* Are you collecting all the available DNA from touched objects? In *International Congress Series* vol. 1239 803–807 (Elsevier, 2003).
- Barta, J. L. *et al.* One of the key characteristics of ancient DNA, low copy number, may be a product of its extraction. *J. Archaeol. Sci.* **46**, 281–289 (2014).
- Mumy, K. L. & Findlay, R. H. Convenient determination of DNA extraction efficiency using an external DNA recovery standard and quantitative-competitive PCR. *J. Microbiol. Methods* **57**, 259–268 (2004).
- Kemp, B. M., Winters, M., Monroe, C. & Barta, J. L. How much DNA is lost? Measuring DNA loss of short-tandem-repeat length fragments targeted by the PowerPlex 16[®] system using the Qiagen MinElute purification kit. *Hum. Biol.* **86**, 313–329 (2014).

11. Doran, A. E. & Foran, D. R. Assessment and mitigation of DNA loss utilizing centrifugal filtration devices. *Forensic Sci. Int. Genet.* **13**, 187–190 (2014).
12. Garvin, A. M. & Fritsch, A. Purifying and concentrating genomic DNA from mock forensic samples using Millipore Amicon filters. *J. Forensic Sci.* **58**, S173–S175 (2013).
13. Norén, L., Hedell, R., Ansell, R. & Hedman, J. Purification of crime scene DNA extracts using centrifugal filter devices. *Investig. Genet.* **4**, 1–8 (2013).
14. Barbaro, A., Staiti, N., Cormaci, P. & Saravo, L. DNA profiling by different extraction methods. In *International Congress Series* vol. 1261 562–564 (Elsevier, 2004).
15. Liu, P. *et al.* Integrated portable polymerase chain reaction-capillary electrophoresis microsystem for rapid forensic short tandem repeat typing. *Anal. Chem.* **79**, 1881–1889 (2007).
16. Gangano, S. *et al.* DNA investigative lead development from blood and saliva samples in less than two hours using the RapidHIT™ Human DNA Identification System. *Forensic Sci. Int. Genet. Suppl. Ser.* **4**, e43–e44 (2013).
17. Dawnay, N. *et al.* Developmental validation of the ParaDNA® screening system-A presumptive test for the detection of DNA on forensic evidence items. *Forensic Sci. Int. Genet.* **11**, 73–79 (2014).
18. Shrivastava, P., Jain, T. & Trivedi, V. B. DNA fingerprinting: A substantial and imperative aid to forensic investigation. *Eur. J. Forensic Sci.* **3**, 23 (2016).
19. Butler, J. M. *Forensic DNA Typing: Biology, Technology, and Genetics of STR Markers*. (Elsevier, 2005).
20. Weispenning, R. *et al.* Doing more with less: Implementing direct amplification with the PowerPlex® 18D System. *Forensic Sci. Int. Genet. Suppl. Ser.* **3**, e409–e410 (2011).
21. Verheij, S., Hartevelde, J. & Sijen, T. A protocol for direct and rapid multiplex PCR amplification on forensically relevant samples. *Forensic Sci. Int. Genet.* **6**, 167–175 (2012).
22. Myers, B. A., King, J. L. & Budowle, B. Evaluation and comparative analysis of direct amplification of STRs using PowerPlex® 18D and Identifier® Direct systems. *Forensic Sci. Int. Genet.* **6**, 640–645 (2012).
23. Vallone, P. M., Hill, C. R. & Butler, J. M. Demonstration of rapid multiplex PCR amplification involving 16 genetic loci. *Forensic Sci. Int. Genet.* **3**, 42–45 (2008).
24. Wang, D. Y., Chang, C.-W., Oldroyd, N. J. & Hennessy, L. K. Direct amplification of STRs from blood or buccal cell samples. *Forensic Sci. Int. Genet. Suppl. Ser.* **2**, 113–114 (2009).
25. Mercier, B., Gaucher, C., Feugeas, O. & Mazurier, C. Direct PCR from whole blood, without DNA extraction. *Nucleic Acids Res.* **18**, 5908 (1990).
26. Ottens, R., Taylor, D., Abarno, D. & Linacre, A. Optimising direct PCR from anagen hair samples. *Forensic Sci. Int. Genet. Suppl. Ser.* **4**, e109–e110 (2013).
27. Ottens, R., Taylor, D., Abarno, D. & Linacre, A. Successful direct amplification of nuclear markers from a single hair follicle. *Forensic Sci. Med. Pathol.* **9**, 238–243 (2013).
28. Ottens, R. *et al.* Application of direct PCR in forensic casework. *Forensic Sci. Int. Genet. Suppl. Ser.* **4**, e47–e48 (2013).
29. Templeton, J. *et al.* Genetic profiling from challenging samples: Direct PCR of touch DNA. *Forensic Sci. Int. Genet. Suppl. Ser.* **4**, e224–e225 (2013).
30. Gray, K., Crowle, D. & Scott, P. Direct amplification of casework bloodstains using the Promega PowerPlex® 21 PCR Amplification System. *Forensic Sci. Int. Genet.* **12**, 86–92 (2014).
31. Ottens, R., Taylor, D. & Linacre, A. DNA profiles from fingernails using direct PCR. *Forensic Sci. Med. Pathol.* **11**, 99–103 (2015).
32. Shokralla, S., Singer, G. A. C. & Hajibabaei, M. Direct PCR amplification and sequencing of specimens' DNA from preservative ethanol. *Biotechniques* **48**, 305–306 (2010).
33. Sorensen, A. *et al.* Direct-to-PCR tissue preservation for DNA profiling. *Int. J. Legal Med.* **130**, 607–613 (2016).
34. Blackie, R., Taylor, D. & Linacre, A. DNA profiles from clothing fibers using direct PCR. *Forensic Sci. Med. Pathol.* **12**, 331–335 (2016).
35. Harty, L. C. *et al.* Collection of buccal cell DNA using treated cards. *Cancer Epidemiol. Prev. Biomark.* **9**, 501–506 (2000).
36. Lebedev, A. V. *et al.* Hot start PCR with heat-activatable primers: A novel approach for improved PCR performance. *Nucleic Acids Res.* **36**, e131–e131 (2008).
37. Hares, D. R. Selection and implementation of expanded CODIS core loci in the United States. *Forensic Sci. Int. Genet.* **17**, 33–34 (2015).
38. Giese, H., Lam, R., Selden, R. & Tan, E. Fast multiplexed polymerase chain reaction for conventional and microfluidic short tandem repeat analysis. *J. Forensic Sci.* **54**, 1287–1296 (2009).
39. Wang, D. Y., Chang, C.-W. & Hennessy, L. K. Rapid STR analysis of single source DNA samples in 2 h. *Forensic Sci. Int. Genet. Suppl. Ser.* **2**, 115–116 (2009).
40. Sahajpal, V., Rajput, S., Sharma, T., Sharma, A. & Thakar, M. K. Development and evaluation of a novel DNA purification buffer and protocol for blood samples on FTA cards. *Forensic Sci. Int. Rep.* **1**, 100014 (2019).
41. Shrivastava, P., Jain, T., Gupta, U. & Trivedi, V. B. Genetic polymorphism study on 12 X STR loci of investigator Argus X STR kit in Bhil tribal population of Madhya Pradesh, India. *Leg. Med.* **17**, 214–217 (2015).
42. Kumawat, R. K. *et al.* Peopling of Rajasthan, India: Evaluating the gene flow from east and west. *Gene Rep.* **22**, 100990 (2020).
43. Kumawat, R. K., Shrivastava, P., Shrivastava, D., Mathur, G. K. & Dixit, S. Genomic blueprint of population of Rajasthan based on autosomal STR markers. *Ann. Hum. Biol.* 1–6 (2020).
44. Hedman, J., Ansell, R. & Nordgaard, A. A ranking index for quality assessment of forensic DNA profiles. *BMC Res. Notes* **3**, 290 (2010).
45. Huggett, J. F. *et al.* Differential susceptibility of PCR reactions to inhibitors: An important and unrecognised phenomenon. *BMC Res. Notes* **1**, 1–9 (2008).
46. Govindarajan, N., Lemalu, A. & Patel, J. Forensic casework methodology for direct PCR amplification of blood swabs. *Forensic Sci. Int. Genet.* **42**, 125–134 (2019).
47. Ludeman, M. J. *et al.* Developmental validation of GlobalFiler™ PCR amplification kit: A 6-dye multiplex assay designed for amplification of casework samples. *Int. J. Legal Med.* **132**, 1555–1573 (2018).
48. Ensenberger, M. G. *et al.* Developmental validation of the PowerPlex® fusion 6C system. *Forensic Sci. Int. Genet.* **21**, 134–144 (2016).
49. Swift, M. L. GraphPad prism, data analysis, and scientific graphing. *J. Chem. Inf. Comput. Sci.* **37**, 411–412 (1997).
50. Rickham, P. P. Human experimentation. Code of ethics of the world medical association. Declaration of Helsinki. *Br. Med. J.* **2**, 177 (1964).
51. Shrivastava, P., Jain, T., Gupta, U. & Trivedi, V. B. Genetic variation at 15 autosomal STR loci in Bhil tribal population of Central India. *Ann. Hum. Biol.* **43**, 81–84 (2016).

Acknowledgements

Authors thankfully acknowledge Zeal India, Promega (India), Qiagen and Health Gene Technologies for providing the multiplex STR kits used in the study. We are also thankful to Dr. M.K. Verma, Director and Dr. Veena Ben Trivedi, Joint Director, State Forensic Science laboratory, Sagar (MP) India for providing the laboratory

facilities and Madhya Pradesh Council of Science and Technology, Bhopal, for sanctioning the financial grant for the research project.

Author contributions

P.S. conceptualized the study, did genotyping and wrote the manuscript, T.J. helped in genotyping. R.K. did the statistical analysis and wrote the manuscript. All authors reviewed the manuscript.

Competing interests

The authors declare no competing interests.

Additional information

Supplementary Information The online version contains supplementary material available at <https://doi.org/10.1038/s41598-021-86633-0>.

Correspondence and requests for materials should be addressed to P.S.

Reprints and permissions information is available at www.nature.com/reprints.

Publisher's note Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.



Open Access This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons licence, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons licence, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons licence and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this licence, visit <http://creativecommons.org/licenses/by/4.0/>.

© The Author(s) 2021