# scientific reports

Check for updates

# **OPEN** Efficacy of a new refractive multifocal contact lens for presbyopia

Do Young Kim<sup>1,2</sup>, Hyunmin Ahn<sup>3</sup>, Sukyung Lee<sup>4</sup>, Ikhyun Jun<sup>1</sup>, Kyoung Yul Seo<sup>1</sup>, Sang Yeop Kim<sup>5</sup> & Tae-Im Kim<sup>1⊠</sup>

This prospective single-arm study aimed to evaluate the clinical efficacy and safety of a refractive multifocal contact lens for the correction of presbyopia in 22 patients. The participants underwent clinical examinations before and 1 week after wearing a refractive multifocal contact lens (OptaCon ZOOM). The primary endpoints were the corrected distance visual acuity (CDVA) and distancecorrected near visual acuity (DCNVA). Defocus curve, contrast sensitivity, and ocular surface disease index (OSDI) were analyzed. A slit-lamp examination was performed for safety analysis. Contact lens comfort and patient satisfaction were assessed using a questionnaire. No significant difference in CDVA was observed before and 1 week after refractive multifocal contact lens use (p = 0.127), whereas DCNVA was significantly improved after 1 week (p < 0.001). The contrast sensitivity was not significantly affected at any spatial frequency under mesopic or photopic conditions. OSDI was significantly increased (p = 0.023). The patient-reported satisfaction scores were 96.2, 91.9, and 85.0 out of 100 at far, intermediate, and near distances, respectively. No significant adverse events were observed. Refractive multifocal contact lenses improved near vision while maintaining distance vision in presbyopic patients, without compromising contrast sensitivity. The study results suggest that OptaCon ZOOM can be considered safe and effective for the correction of presbyopia.

Keywords Multifocal lens, Contact lens, Presbyopia

Presbyopia is the age-related decrease in human accommodative amplitude that begins in the second decade of life and becomes clinically significant between 40 and 45 years of age<sup>1,2</sup>. Given the expected increase in life expectancy and proportion of aged population, presbyopia is a significant challenge for eye care practitioners due to its increasing prevalence with age<sup>3</sup>. The global prevalence of presbyopia is estimated to be 1.8 billion, or approximately 22% of the world's population. Uncorrected or under-corrected presbyopia is estimated to be present in around 1.3 billion people and is significantly associated with reduced quality of life and decreased productivity in daily activities<sup>4-6</sup>.

Various methods including contact lenses, spectacles, and surgical interventions that modify the optics of the cornea or replace the crystalline lens have been proposed to correct the refractive error associated with presbyopia<sup>4,7</sup>. There has been constant interest in using contact lens for presbyopia because of their optical stability on the visual axis with the eye movement and cosmetic appeal<sup>8</sup>. Currently, various forms of multifocal contact lenses are available for presbyopia correction. The alternating-vision design has separate areas for near and distance vision. However, this design requires a translational change to relocate the most pupil to near zones and has the issue of image jumps, resulting in decreased depth perception and longer adaptation time<sup>4,9</sup>. In addition, simultaneous vision using multifocal contact lenses has been increasingly used for presbyopia correction. Simultaneous vision refers to foveal viewing of the images located at different distances along the optical axis for presbyopia correction. The pupil size and power distribution across the lens are essential for providing optical zones with different optical powers<sup>1,10</sup>. The concentric design features a central distance or near power zone surrounded by one or more annulus of opposite power in repeating pattern. The sharp transition between the near and distance zones can compromise the visual performance<sup>4,11</sup>.

The aspheric multifocal contact lens provides a gradual transition of the lens power between the distance and near powers, which increases the depth of focus by inducing spherical aberration. The induced spherical

<sup>1</sup>Department of Ophthalmology, Severance Hospital, Yonsei University College of Medicine, 50 Yonsei-ro, Seodaemun-gu, Seoul 03722, Korea. <sup>2</sup>Siloam Eye Hospital, Seoul, Republic of Korea. <sup>3</sup>Eyejun Ophthalmic Clinic, Seoul, Republic of Korea. <sup>4</sup>Yonsei Plus Eye Center, Seongnam, Republic of Korea. <sup>5</sup>Clear Eye Clinic, Pyeongtaek, Republic of Korea. <sup>⊠</sup>email: TIKIM@yuhs.ac

aberration may compromise the retinal image quality but it allows for increase in the vergence range at the same time<sup>3</sup>. Aspheric lens designs include center-near and center-distance designs. The center-near design has the highest plus power in the center, with a gradual decrease in power toward the periphery, whereas the center-distance design has the lowest plus power in the center, with a gradual increase in power toward the periphery<sup>3,12,13</sup>.

This study aimed to evaluate the efficacy of the new center-near aspheric refractive multifocal contact lenses that use simultaneous-image principle for the correction of presbyopia and report any side effects or complications after its application.

## Methods

### Ethics

This single-arm prospective study was performed in Severance Hospital, Yonsei University College of medicine, Seoul, Korea from June 2022 to August 2022. The study was approved by the Severance Hospital Institutional Review Board (IRB approval no. 2021-3385-001). It adhered to the principles of the Declaration of Helsinki. Informed consent was obtained from all subjects after the explanation of the nature and possible consequences of the study.

#### **Study participants**

The current study was conducted on 24 patients with presbyopia aged 40 years and older. Two participants dropped out as they were lost to follow-up and 22 participants completed the study. The inclusion criteria were as follows: (1) presbyopia, addition of 0.75 D or more in both eyes; (2) astigmatism: up to -1.00 D of cylinder power in both eyes (according to subjective refraction); (3) distance visual acuity of 0.3 logarithm of the minimal angle of resolution (logMAR) or better with contact lens fitting; and (4) no clinical history of ocular surface disease or dry eye symptoms. The exclusion criteria were as follows: (1) patients unsuitable for contact lens wear in their activities or occupations; (2) ocular disease or abnormalities that required regular topical medication or might have interfered with the ophthalmic examinations required for the current study; (3) the presence of chronic or moderate-to-severe ocular surface disease, including dry eye disease; and (4) a history of ocular surgery within 3 months. The dry eye symptoms were identified based on the report of the diagnostic methodology subcommittee of the international dry eye workshop (2007)<sup>14</sup>.

#### Measurements

All participants underwent a complete ophthalmic examination prior to the study. Screening for ocular and systemic diseases, visual acuity, subjective refraction, near add power, eye dominance, and comprehensive slit lamp examination were performed. All the measurements were conducted by the same researcher (Tae-im, Kim). Keratometry was performed using a Topcon KR 8800 auto-kerato-refractometer (Topcon Corporation, Tokyo, Japan). The primary outcomes were corrected distance visual acuity (CDVA) and distance-corrected near visual acuity (DCNVA). Subjective distance refraction was performed using standard optometric measures and all measurements were taken with distance corrections. The VA measurement was cut at 0.00 logMAR. Distance visual acuity was measured at 4 m using an early treatment diabetic retinopathy study (ETDRS) logMAR chart (Precision Vision, La Salle, IL, USA) and near visual acuity was measured with the Logarithmic Visual Acuity Chart 2000 New ETDRS chart 1 for testing at 40 cm (Precision Vision, La Salle, IL, USA). Near add was defined as the minimum plus power over full subjective distance refraction required to read 0.0 logMAR print on the ETDRS chart. For defocus curve measurement, the eyes were defocused to + 3.00 D spherical from a subjective refraction value. Then, the trial lenses with negative spherical power were added in increments of -0.50 D until a power of -5.00 D was reached. The binocular visual acuity was measured using the ETDRS logMAR chart at 4 m with every increment. Defocus curves were obtained before and after 1 week of multifocal contact lens wear and were compared in each participant. We intended to measure each subject's binocular visual acuity to simulate the subject's visual performance at different distances and reflect the subjects' experience in most day-to-day tasks. (Fig. 1) Binocular contrast sensitivity was measured under photopic (target luminance, 85 candelas/square meter [cd/m<sup>2</sup>]) and mesopic (target luminance, 3 cd/m<sup>2</sup>) conditions. The functional acuity contrast test (E.A.C.T) was performed using the Optec 6500 view-in test system (Stereo Optical Co., Inc., Chicago, IL, USA). Ocular surface disease index (OSDI) scores were obtained before and after 1 week of multifocal soft contact lens wear. All participants were asked to complete a questionnaire regarding visual symptoms, dry eye symptoms, overall satisfaction with the lenses, satisfaction with glasses, and recommendation to others. The questionnaire, created by combining and modifying previously reported questionnaire to evaluate postoperative symptoms, spectacle independence, and satisfaction, is shown in Table 1<sup>5,15</sup>. Slit lamp examination was performed to evaluate corneal erosion, edema, and peripheral injection for safety analysis.

#### Multifocal soft contact lens

OptaCon ZOOM contact lens (INTEROJO Co. Ltd.. Pyeongtaek-si, Gyeonggi-do, Korea) is a disposable multifocal soft contact lens for presbyopia with an aspheric center-near design. (Fig. 2) The OptaCon ZOOM is designed with a 1 mm central radius optimized for near vision and exhibits gradual changes in power from the relative central plus to negative spherical aberration in the periphery. The OptaCon ZOOM is manufactured by cast molding with silicone hydrogel which is currently under a naming process. The OptaCon ZOOM has a base curve of 8.6 mm and a diameter of 14.2 mm. (Table 2) Lenses were prescribed for refractive error with a target of emmetropia, according to the manufacturer's trial set guidelines. The same lens used for distance vision was placed in a trial frame for near vision assessment. Near visual acuity was assessed using binocular examination. Overcorrection was performed to achieve the optimal best-corrected visual acuity for both far and near distances. The final lens prescription was adjusted based on subjective refraction, as the OptaCon ZOOM is available in



**Figure 1.** Defocus curve for the mean CDVA values before and after 1 week of multifocal soft contact lens wear. The curve had a U shape with a peak at 0.00D. The mean CDVA at the peak was not significantly different before and after 1 week of lens wear (*p*-value=0.317). The result of direct comparison of defocus curve revealed that multifocal soft contact lens wear showed significant improvement of visual outcome for both near and far. It resulted in a significantly improved outcome for near vision particularly at – 2.00D to – 5.00D, corresponding with a 20–50-cm range (*p*-value < 0.05 in all cases). The result also showed significant improvement in outcome for far vision at + 2.00, + 2.50 and + 3.00D (*p*-value < 0.05 in all cases).

various add powers. Participants with add power lower than 1.0 D received a low addition lens of 1.0 D, those with add power between 1.0 D and 1.75 D received a moderate addition lens of +1.50 D to +2.0 D, and those with add power of +2.0 D or higher received a high addition lens of 2.50 D. Participants were required to wear each daily disposable lens for at least 4 h daily for 1 week, including the follow-up visit. The lenses were discarded after a single day of wear. Lens removal was recommended in cases of severe adverse reactions including pain, abnormal discharge, and injection.

#### Statistical analysis

A non-parametric test was used for statistical analysis because of the small sample size. Wilcoxon's signed-rank test was used to compare the measurements obtained before and after OptaCon ZOOM wear. The statistical software package SPSS (version 22.0, IBM Corp., Armonk, NY, USA) was used for statistical analysis. A *p*-value of < 0.05 was considered statistically significant.

#### Results

A total of 22 patients with presbyopia were included in the study. The majority were female (16 [72.72%]), and the mean age was  $49.86 \pm 4.23$  years. The mean distance CDVA (DCDVA) was  $0.003 \pm 0.015$  logMAR for the right eye and  $0.00 \pm 0.00$  logMAR for the left eye. The mean DCNVA was  $0.19 \pm 0.15$  logMAR for the right eye and  $0.17 \pm 0.15$  for the left eye. The mean add power was  $1.51 \pm 0.49$  D for both eyes. The mean spherical equivalent (SE) was  $-1.98 \pm 1.80$  D for the right eye and  $-1.85 \pm 1.66$  D for the left eye. Flat curvature power (K1), steep curvature power (K2), and steeper axis were  $42.5 \pm 2.08$  D,  $43.23 \pm 2.29$  D, and  $89.86 \pm 21.90^{\circ}$  for the right eye and  $42.32 \pm 2.08$  D,  $43.03 \pm 2.20$  D, and  $87.45 \pm 25.60^{\circ}$  for the left eye, respectively. (Table 3).

The DCDVA was not significantly different before and after 1 week of contact lens use  $(0.00 \pm 0.11 \text{ vs.} 0.01 \pm 0.03, p = 0.127)$ . DCNVA showed significant improvement from  $0.19 \pm 0.15 \log$ MAR to  $0.04 \pm 0.08 \log$ -MAR after 1 week of contact lens wear (p < 0.001). (Table 4).

There was no significant difference in the results of the binocular contrast sensitivity test for any frequency under photopic and mesopic conditions (Fig. 3). Thus, multifocal soft contact lens wear was not significantly associated with a decrease in binocular contrast sensitivity under mesopic or photopic conditions.

The OSDI score was significantly increased after 1 week of multifocal soft contact lens use (p = 0.02; Table 4). No significant adverse events were observed on slit lamp examination. Dryness was the most common discomfort symptom with the mean score of  $1.38 \pm 1.24$  out of 4, followed by foreign sensation with the mean score of  $1.29 \pm 1.27$  out of 4. For visual symptom, diplopia and fluctuation of visual acuity were the most common symptoms with the mean score of  $1.81 \pm 1.24$  and  $1.57 \pm 1.21$  out of 4, respectively. The mean patient-reported satisfactory scores converted to a score of 100 obtained by the questionnaire were 96.2, 91.9, and 85.0 for distance, intermediate, and near vision, respectively. The lens satisfaction score comparing glasses and the score

Question	Never (0)	Sometimes (1)	Half of the time (2)	Almost always (3)	Always (4)	
1. Do you feel a sense of discomfort after wearing lenses?						
Photosensiti	vity					
Ocular pain						
Dryness						
Discomfort	or Foreign sensation					
2. Do you ex	perience any of the follo	wing visual sympto	ms after wearing lenses	?		
Glare						
Halo						
Diplopia						
Fluctuation	of visual acuity					
3. Do you ha	ve difficulty performing	the following tasks	;			
(a) Near dista	ance tasks					
Reading me	dication instruction					
Reading nev	rspaper					
Cleaning na	ls					
Using cellph	one					
Reading mag	gazine or books					
(b) Intermed	iate distance tasks					
Finding iten	is in the kitchen					
Computer w	Computer work					
Looking at the bathroom mirror						
Looking at the calendar						
Viewing a picture or a painting in a picture frame						
(c) Far distance tasks						
Reading road signs						
Identifying a person across the room						
Estimating the distance between cars						
Reading the house number or address						
Checking the time on the wall clock						
Question	Very dissatisfied (0)	Dissatisfied (1)	Neutral (2)	Satisfied (3)	Very satisfied(4)	
4. Do you fin	4. Do you find wearing OptaCon ZOOM contact lens more satisfactory compared to wearing glasses?					
5. Would you recommend OptaCon ZOOM contact lens to other people?						

 Table 1. Questionnaire regarding discomfort, visual symptoms, and ability to perform tasks.



**Figure 2.** Diagram illustrating the lens design of OptaCon ZOOM. OptaCon ZOOM is manufactured from silicone hydrogel by cast molding method, with base curve of 8.6 mm and lens diameter of 14.2 mm. OptaCon ZOOM is designed with near vision optimized central 1 mm radius and gradual changes to negative spherical aberration in the periphery, with max vertex refractive power in 4.5 mm zone.

for recommending the product were slightly above average, with scores of 2.10 and 2.05 out of 4, respectively. (Table 5).

Parameter	Details
Base curve	8.6 mm
Diameter	14.2 mm
Power	+ 0.00D ~ - 6.00D (0.25D steps)
ADD power	Low(+1.00D), Med(+1.50D, +2.00D), High(+2.50D)
Material	Silicone Hydrogel
Water content	45%
Manufacturing method	Cast molding
Oxygen permeability (Dk)	$70 \times 10^{-11} (\text{cm}^2/\text{s})[\text{ml O}_2/(\text{ml mmHg})$
Oxygen transmissibility (Dk/t) @-3.00D	$100 \times 10^{-9} (\text{cm/s}) [\text{ml O}_2/(\text{ml mmHg})]$

Table 2.	The detailed	parameter of O	ptaCon 1-da	ay ZOOM Multifocal contact lens.
----------	--------------	----------------	-------------	----------------------------------

Baseline characteristics	OD	OS		
Visual acuity				
DCDVA (log MAR)	$0.00 \pm 0.015$	$0.00 \pm 0$		
DCNVA (log MAR)	$0.19 \pm 0.15$	$0.17 \pm 0.15$		
Add power (D)	$1.51\pm0.49$	$1.51\pm0.49$		
Refractive errors				
Sph (D)	$-1.78 \pm 1.80$	$-1.68 \pm 1.69$		
Cyl (D)	$-0.40 \pm 0.35$	$-0.34 \pm 0.26$		
Axis (°)	116±63	$101 \pm 59$		
K values				
K1 (D)	$42.50\pm2.08$	$42.32\pm2.08$		
K2 (D)	$43.23 \pm 2.29$	$43.03 \pm 2.20$		
Axis (°)	90±22	87±26		

**Table 3.** Patient demographics. Values are presented as mean  $\pm$  standard deviation. OD—right eye, OS—lefteye, DCDVA—distance corrected distance visual acuity, DCNVA—distance corrected near visual acuity,OSDI—ocular surface disease index, Sph—spherical power of subjective refraction, CYL—cylindrical power of subjective refraction.

	Before wearing lens	After wearing lens	P value
DCDVA (logMAR)	$0.00 \pm 0.11$	$0.01 \pm 0.03$	0.13
DCNVA (logMAR)	$0.19 \pm 0.15$	$0.04 \pm 0.08$	< 0.001
OSDI	16.15	28.45	0.02

**Table 4.** Visual acuity of enrolled eye with and without wearing multifocal contact lenses. Values are presented as mean ± standard deviation. *DCDVA*—distance corrected distance visual acuity, *DCNVA*—distance corrected near visual acuity, *CNVA*—corrected near visual acuity, *OSDI*—ocular surface disease index.

.....

## Discussion

Presbyopia is expected to increase sharply with the rapid growth of the aging population; therefore, effective treatment modalities and strategies for the correction of presbyopia are becoming increasingly important (2). The present study analyzed the clinical performance of a refractive multifocal soft contact lens for the correction of presbyopia and found that the lens improved near vision while correcting distance vision. Previous studies on similar age groups that assessed the visual performance of presbyopic contact lens showed similar results<sup>16,17</sup>.

Monovision, which is a system for fitting one eye for distance vision and the other for near vision, is another common option for contact lens correction of presbyopia. In a previous randomized controlled study, both multi-focal and monovision contact lens wearers maintained at least 20/20 binocular vision, whereas multifocal wearers showed significantly lower vision at near distances under low-contrast conditions. However, self-reported quality of life and lens preference were significantly higher with multifocal contact lenses than in those with the baseline and monovision, possibly because of comparable visual acuity without compromising stereoacuity to the same degree as in monovision<sup>12</sup>. Monovision performs better for distance and near visual acuity, whereas a center-near aspheric multifocal contact lens provide a better balance of real-world function<sup>18</sup>.



**Figure 3.** Contrast sensitivity values before and after 1 week of OptaCon ZOOM wear for the correction of presbyopia. The images show binocular contrast sensitivity values under photopic (**A**) and mesopic (**B**) conditions.

Another previous study showed that both near and distance visual acuity were significantly improved in patients with presbyopia with multifocal contact lenses compared to uncorrected baseline visual acuity, with minor adverse effects<sup>7</sup>. Our study result showed significant improvement in DCNVA to fair value of  $0.04 \pm 0.08 \log$ MAR. For the result of other MCLs, near VA was reported to be  $0.05 \pm 0.06$  for PureVision multifocal,  $0.49 \pm 0.21$  in mild presbyopic group and  $0.48 \pm 0.20$  in moderate to severe presbyopia group for Air Optix<sup>\*</sup> Aqua Multifocal at 40 cm, and  $0.102 \pm 0.055$  for Morning QI-plus<sup>1,19,20</sup>.

The contrast sensitivity test is associated with the performance of spatial vision in both size and contrast, and it correlates better with identifying real-world objects, providing a better measure of visual sensitivity and quality than visual acuity alone<sup>18,21</sup>. Our study showed no significant difference in contrast sensitivity before and after multifocal contact lens wear at all spatial frequencies under both mesopic and photopic conditions. The results of the current study correlate with those of previous studies comparing different types of lenses<sup>18,22</sup>. However, some studies have reported conflicting results regarding the impact of multifocal contact lenses on contrast sensitivity, with a greater reduction in contrast sensitivity shown with multifocal lens wear than with single vision lenses under low-light conditions<sup>10,23,24</sup>.

In our study, the OSDI was significantly increased after wearing contact lenses compared with the baseline evaluation. These results appear to be consistent with those of other studies, in which contact lenses seemed to disrupt the ocular surface and tear film<sup>25,26</sup>. Given the high prevalence of dry eye disease in older population, dry eye symptoms must be considered when prescribing contact lens for presbyopia correction. However, the higher proportion of females in the current study may have led to an increase in the OSDI score, since previous studies have reported that females are more likely to report symptoms of dryness regardless of age<sup>26,27</sup>. Moreover,

	Mean ± SD			
1. Do you feel a sense of discomfort after wearing lenses?				
Photosensitivity	$0.52\pm0.68$			
Ocular pain	$0.43\pm0.68$			
Dryness	$1.38 \pm 1.24$			
Discomfort or Foreign sensation	$1.29 \pm 1.27$			
2. Do you experience any of the following visual symptoms after wearing lenses?				
Glare	$0.67\pm0.86$			
Halo	$1.19 \pm 1.21$			
Diplopia	$1.81 \pm 1.24$			
Fluctuation of visual acuity	$1.57 \pm 1.21$			
3. Do you have difficulty performing the following tasks?				
(a) Near distance				
Reading medication instruction	$3.05 \pm 1.02$			
Reading newspaper	$3.43\pm0.93$			
Cleaning nails	$3.33 \pm 1.15$			
Using cellphone	$3.33 \pm 0.87$			
Reading magazine or books	$4.00\pm0.97$			
(b) Intermediate distance				
Finding items in the kitchen	$3.14\pm0.00$			
Computer work	$3.71 \pm 1.01$			
Looking at the bathroom mirror	$3.71\pm0.72$			
Looking at the calendar	$3.62\pm0.80$			
Viewing a picture or a painting in a picture frame	$3.90\pm0.44$			
(c) Far distance				
Finding items in the kitchen	$3.90\pm0.44$			
Computer work	$3.90\pm0.44$			
Looking at the bathroom mirror	$3.90\pm0.45$			
Looking at the calendar	$3.71\pm0.72$			
Viewing a picture or a painting in a picture frame	$3.81 \pm 0.60$			
4. Do you find wearing OptaCon ZOOM contact lens more satisfactory compared to wearing glasses?	$2.10 \pm 1.04$			
5. Would you recommend OptaCon ZOOM contact lens to other people?	$2.05\pm0.80$			

 Table 5. Questionnaire regarding discomfort, visual symptoms, and ability to perform tasks.

the apparent lack of association between symptoms of dry eye and available clinical tests has been reported in previous studies<sup>26,28,29</sup>.

Although the multifocal refractive contact lens significantly improved near vision, patient satisfaction was limited by the lack of near-full correction expected by the patient. The OptaCon ZOOM has a near power zone of 1.0 mm, which may result in the symptoms of diplopia and fluctuation of vision and relatively narrow depth of focus for intermediate distances because of the sharp dioptric transition from near center zone to peripheral optic zone. Further studies are required to address this issue.

The current study had some limitations. First, the study had a small sample size and a short follow-up period of 1 week after lens wear. However, we expected relatively short adaptation period to be required since applying contact lenses provides an immediate improvement in VA. We expected a follow-up period of 1 week would be sufficient as it was in past studies. Still, the long-term efficacy of multifocal soft contact lenses should be further investigated in future studies with a larger sample size and a longer follow-up period. Furthermore, all patients included in this study had either no prior experience with contact lenses or had not used them for an extended period. The difficulty that the participants faced in adapting to the lenses may have influenced the results of the trial, including visual symptoms, ocular discomfort, and overall satisfaction. Taking this aspect into account during the prescription process may lead to higher levels of satisfaction.

In conclusion, refractive multifocal contact lenses improved both far and near vision in patients with presbyopia without compromising contrast sensitivity. The results of the current study suggest that an aspheric centernear designed multifocal contact lens can be considered a minimally invasive measure for presbyopia correction.

#### Data availability

All data generated or analyzed during this study are included in this published article. The datasets used during the current study are available from the corresponding author on reasonable request.

Received: 6 June 2023; Accepted: 28 February 2024 Published online: 17 April 2024

#### References

- Llorente-Guillemot, A., García-Lazaro, S., Ferrer-Blasco, T., Perez-Cambrodi, R. J. & Cerviño, A. Visual performance with simultaneous vision multifocal contact lenses. *Clin. Exp. Optom.* 95(1), 54–59 (2012).
- Park, W. C., Kim, C. H. & Jin, S. W. Clinical performance of presbyopia correction with a multifocal corneoscleral lens. Int. J. Ophthalmol. 14(4), 529 (2021).
- Pérez-Prados, R., Piñero, D. P., Pérez-Cambrodí, R. J. & Madrid-Costa, D. Soft multifocal simultaneous image contact lenses: A review. Clin. Exp. Optom. 100(2), 107–127 (2017).
- 4. Katz, J. A. et al. Presbyopia-a review of current treatment options and emerging therapies. Clin. Ophthalmol. 15, 2167–2178 (2021).
- 5. Park, S. Y. et al. Clinical Efficacy of Pinhole Soft Contact Lenses for the Correction of Presbyopia (Taylor & Francis, 2019).
- Fricke, T. R. *et al.* Global prevalence of presbyopia and vision impairment from uncorrected presbyopia: Systematic review, metaanalysis, and modelling. *Ophthalmology* 125(10), 1492–1499 (2018).
- Kim, H. Clinical results of multifocal soft contact lens wear in presbyopia. J. Korean Ophthalmic Opt. Soc. 25(3), 265–272 (2020).
   Wolffsohn, J. S., Davies, L. N. & Sheppard, A. L. New insights in presbyopia: Impact of correction strategies. BMJ Open Ophthalmol.
  - **8**(1), e001122 (2023).
- 9. Zanen, A. The advantages and disadvantages of bifocal lenses. Bull. Soc. Belge Dophtalmologie 264, 71-78 (1997).
- Remon, L., Pérez-Merino, P., Macedo-de-Araujo, R. J., Amorim-de-Sousa, A. I. & González-Méijome, J. M. Bifocal and multifocal contact lenses for presbyopia and myopia control. J. Ophthalmol. 2020, 1–18 (2020).
- Gispets, J., Arjona, M., Pujol, J., Vilaseca, M. & Cardona, G. Task oriented visual satisfaction and wearing success with two different simultaneous vision multifocal soft contact lenses. J. Optom. 4(3), 76–84 (2011).
- Richdale, K., Mitchell, G. L. & Zadnik, K. Comparison of multifocal and monovision soft contact lens corrections in patients with low-astigmatic presbyopia. Optom. Vis. Sci. 83(5), 266–273 (2006).
- Bakaraju, R. C., Ehrmann, K., Ho, A. & Papas, E. Inherent ocular spherical aberration and multifocal contact lens optical performance. Optom. Vis. Sci. 87(12), 1009–1022 (2010).
- Bron, A. J. et al. Methodologies to diagnose and monitor dry eye disease: Report of the diagnostic methodology subcommittee of the international dry eye workshop (2007). Ocular Surf. 5(2), 108–152 (2007).
- Kim, J. S., Jung, J. W., Lee, J. M., Seo, K. Y. & Kim, E. K. Clinical outcomes following implantation of diffractive multifocal intraocular lenses with varying add powers. Am. J. Ophthalmol. 160(4), 702–709 (2015).
- Kim, E., Bakaraju, R. C. & Ehrmann, K. Power profiles of commercial multifocal soft contact lenses. Optom. Vis. Sci. 94(2), 183 (2017).
- Sha, J. et al. Visual performance of daily-disposable multifocal soft contact lenses: A randomized, double-blind clinical trial. Optom. Vis. Sci. 95(12), 1096–1104 (2018).
- 18. Gupta, N., Naroo, S. A. & Wolffsohn, J. S. Visual comparison of multifocal contact lens to monovision. *Optom. Vis. Sci.* 86(2), E98-E105 (2009).
- 19. Sha, J. *et al.* Short-term visual performance of soft multifocal contact lenses for presbyopia. *Arquivos Brasileiros de Oftalmologia* **79**, 73–77 (2016).
- Kim, H. K., Kim, H. M. & Chung, S. K. Effectiveness of multifocal soft contact lens for presbyopia. J. Korean Ophthalmol. Soc. 49(5), 727–731 (2008).
- Owsley, C. & Sloane, M. E. Contrast sensitivity, acuity, and the perception of 'real-world'targets. Br. J. Ophthalmol. 71(10), 791–796 (1987).
- Collins, M. J., Brown, B. & Bowman, K. J. Contrast sensitivity with contact lens corrections for presbyopia. Ophthalmic Physiol. Opt. 9(2), 133–138 (1989).
- Nti, A. N., Gregory, H. R., Ritchey, E. R., Wolffsohn, J. S. & Berntsen, D. A. Contrast sensitivity with center-distance multifocal soft contact lenses. Optom. Vis. Sci. 99(4), 342–349 (2022).
- 24. Fisher, K. Presbyopic visual performance with modified monovision using multifocal soft contact lenses. *Int. Contact Lens Clin.* 24(3), 91–100 (1997).
- Lemp, A. Report of the national eye institute/industry workshop on clinical trials in dry eyes. Eye Contact Lens 21(4), 221–232 (1995).
- 26. du Toit, R., Situ, P., Simpson, T. & Fonn, D. The effects of six months of contact lens wear on the tear film, ocular surfaces, and symptoms of presbyopes. *Optom. Vis. Sci.* 78(6), 455–462 (2001).
- Doughty, M. J. et al. A patient questionnaire approach to estimating the prevalence of dry eye symptoms in patients presenting to optometric practices across Canada. Optom. Vis. Sci. Off. Publ. Am. Acad. Optom. 74(8), 624–631 (1997).
- Rueff, E. M., Jones-Jordan, L. A. & Bailey, M. D. A randomised clinical trial of multifocal contact lenses and contact lens discomfort. Ophthalmic Physiol. Opt. 41(1), 93–104 (2021).
- Pult, H., Purslow, C., Berry, M. & Murphy, P. J. Clinical tests for successful contact lens wear: relationship and predictive potential. Optom. Vis. Sci. 85(10), E924–E929 (2008).

#### Acknowledgements

This work was supported by the Korea Medical Device Development Fund grant funded by the Korea government (the Ministry of Science and ICT, the Ministry of Trade, Industry and Energy, the Ministry of Health & Welfare, the Ministry of Food and Drug Safety) (Project Number: 1711196706, RS-2023-00253749).

#### **Author contributions**

H.M.A., S.K.L., S.Y.K., I.H.J., K.Y.S. and T.I.K. participated in study design, hypothesis generation, data collection and screening. D.Y.K., H.M.A., S.K.L., and T.I.K. participated in data analysis, literature research and evidence synthesis. D.Y.K., H.M.A., I.H.J., K.Y.S. and T.I.K. wrote the manuscript with contributions from all co-authors. All authors have reviewed the manuscript.

#### **Competing interests**

The authors declare no competing interests.

#### Additional information

Correspondence and requests for materials should be addressed to T.-I.K.

Reprints and permissions information is available at www.nature.com/reprints.

**Publisher's note** Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

**Open Access** This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons licence, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons licence, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons licence and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this licence, visit http://creativecommons.org/licenses/by/4.0/.

© The Author(s) 2024