

OPEN

Received: 03 November 2014 Accepted: 03 July 2015 Published: 07 August 2015

Microfluidic Chip-LC/MS-based Glycomic Analysis Revealed Distinct *N*-glycan Profile of Rat Serum

Wei-Na Gao¹, Lee-Fong Yau¹, Liang Liu¹, Xing Zeng², Da-Can Chen², Min Jiang³, Ju Liu³, Jing-Rong Wang¹ & Zhi-Hong Jiang¹

The rat is an important alternative for studying human pathology owing to certain similarities to humans. Glycomic studies on rat serum have revealed that variations in the *N*-glycans of glycoproteins correlated with disease progression, which is consistent with the findings in human serum. Therefore, we comprehensively characterized the rat serum *N*-glycome using microfluidic chip-LC-ESI-QTOF MS and MS/MS techniques. In total, 282 *N*-glycans, including isomers, were identified. This study is the first to present comprehensive profiling of *N*-glycans containing *O*-acetylated sialic acid, among which 27 *N*-glycans are novel. In addition, the co-existence of *N*-acetylneuraminic acid (NeuAc) and *N*-glycolylneuraminic acid (NeuGc) in a single *N*-glycan ('mixed' *N*-glycan) was detected and represents a new type of *N*-glycan in rat serum. The existence of *O*-acetylated sialic acid is the characteristic feature of rat serum that distinguishes it from mouse and human sera. Comparisons between the rat, mouse, and human serum glycomes revealed that the rat glycome is more similar to that of human sera than to that of mouse sera. Our findings highlight the similarities between the glycomic profile of rat and human sera and provided important selection criteria for choosing an appropriate animal model for pathological and pharmacological studies.

Owing to the inherent similarity to humans, rats have been used as subjects for scientific research in many experimental studies that have contributed significantly to our understanding of a variety of scientific areas including genetics, human diseases, and pharmacological studies. Among rodents, rats are physiologically most similar to humans. For instance, the number of detoxifying P450 gene in rats is similar to that in humans than to that in mice, which explains superiority of rats over mice in pharmacologic investigations¹. Besides, ribosomal protein genes, which play important role to build up ribosomes, have high sequence similarity between rat and human². A rat model is therefore superior to a mouse model for assaying the pharmacodynamics and toxicities of drugs¹ and for studying various human diseases, such as cardiovascular disease, arthritis, diabetes mellitus, and autoimmune disorders³, and numerous rat models have been established for such types of *in vivo* research.

In the post-genomic era, glycosylation has been brought to public attention as a frontier⁴. The *N*-glycans of glycoproteins are directly involved in almost every biological process and play a crucial role in human diseases due to their unusual biological selectivities⁵. Consequently, glycomics, which is one of the latest members in the "-omics" family of the life sciences, has started to gain its respectable role

²State Key Laboratory of Quality Research in Chinese Medicines, Macau Institute for Applied Research in Medicine and Health, Macau University of Science and Technology, Macau, China. ²Guangdong Provincial Hospital of Chinese Medicine, Second Clinical College of Guangzhou University of Chinese Medicine, Guangdong Provincial Academy of Chinese Medical Sciences, Guangzhou, China. ³Division of Rheumatology, Jiujiang First People's Hospital, Jiujiang, China. Correspondence and requests for materials should be addressed to J.-R.W. (email: jrwang@must.edu.mo) or Z.-H.J. (email: zhjiang@must.edu.mo)

based on its substantially impact on research towards diagnosing, preventing, and treating diseases^{6,7}. For example, differences in the glycosylation of serum proteins between cancer patients and healthy people have been distinguished for potential biomarker identification⁸. These non-invasive *N*-glycan biomarkers are promising for disease prognosis. As it is very important to develop reliable strategies to monitor progression of disease in the animal model, the glycomic investigation by using animal model is a necessity. Therefore, glycomics-based studies using animal models have been increasingly reported. For example, glycomic studies of liver cancer in rats revealed that variations in fucosylation of *N*-glycans in serum glycoproteins was closely related to the progression of hepatocellular carcinomas. This was consistent with observations for human liver cancer⁹⁻¹¹. However, in comparison to the vast array of research into the human glycome, limited studies have been performed on the rat glycome. Less than 50 *N*-glycan compositions of rat serum have been reported to date^{9,10}. Therefore, a comprehensive glycomic study in rats is much needed. The findings from which may aid our understanding of the physiology and pathology associated with glycosylation for various diseases.

Although direct measurements of free N-glycans may seem to be limited in that the information associated with the structure of the glycoprotein is lost owing to the deglycosylation step, such measurements are still significative because: (1) N-glycans are often the crucial functional elements in cellular and biomolecular interactions; (2) glycomic techniques are methodologically easier than the glycoproteomic approaches ⁵. To carry out a comprehensive glycomic study, the primary analytical challenge is the detection and identification of rat-specific N-glycans, especially those in low abundance and unknown N-glycans. The analytical method must therefore display high detection sensitivity. Moreover, because of the inherent complexity of N-glycans, numerous N-glycans with similar structures might be co-eluted, leading to significant ionization suppression and a compromised MS/MS spectrum. Therefore, highly efficient separation of N-glycans is also imperative. In this study, we developed a microfluidic porous graphitized carbon (PGC) chip-LC/MS-based approach for the profiling of rat serum N-glycans. As a type of nano-LC, this microfluidic-chip has many advantages over routine HPLC and UHPLC^{12,13}. First, the microfluidic chip-LC shows favorable mass spectrometric sensitivity because of reduced flow rates, minimized dead volume, and lower peak dispersion¹⁴. Second, the PGC-based stationary phase enables isomer-specific separation of N-glycans, and thus facilitates maximum resolution of the overlapping ions¹⁵. Therefore, employment of this microfluidic PGC chip-LC/MS technique in glycomic studies would provide high sensitivity and high resolution. Notably, in this study, higher pH and high ionic strength mobile phase instead of routine mobile phase was applied to PGC-Chip for the first time, chromatographic and MS conditions were then carefully optimized, leading to remarkable improvement in the detection of sialylated N-glycans, especially multi-sialylated N-glycans. This report describes the comprehensive glycome profiling of rat serum using this newly developed technique and then reports glycomic comparisons of rat, mouse, and human sera.

Results

Preparation of *N***-glycan sample for glycomic profiling.** The initial step for *N*-glycan analysis is the release of *N*-glycans from glycoproteins. Enzymatic hydrolysis using peptide *N*-glycosidase F (PNGase F) is the most frequently used method for this step. A variety of protocols with varied digestion buffers, denaturation conditions, and digestion times have been employed in previous studies^{16,17}. To select the conditions that allow for highly efficient release of *N*-glycans without affecting their stability, the three most frequently used protocols were selected for comparison, these included: (1) digestion with thermal denaturation in SDS-containing buffer, with the addition of DTT¹⁸; (2) digestion in ammonium bicarbonate buffer^{16,17}; (3) digestion with thermal denaturation in ammonium bicarbonate buffer¹⁹. The results indicated that the intensities of acidic *N*-glycans varied greatly under different digestion conditions, while those of neutral *N*-glycans obtained by different methods were relatively stable (Supplementary Figure S1). Detailed comparison of individual *N*-glycans suggested that acidic *N*-glycans with one sialic acid were more sensitive to the digestion conditions than those with two or more sialic acids. Hence, serum samples were enzymatically digested using Method (1) in the following experiments.

The next step in the process is purification of the released *N*-glycans from the digestion mixture. In the digestion mixture, de-glycosylated proteins and various salts need to be removed prior to LC-MS analysis so as to avoid contamination of the chip and ion suppression of the signal²⁰. These procedures are especially crucial for the nano-LC system, because the proteins and salts may precipitate to block the micro-column in the chip and the spray needle²¹. Purification of *N*-glycans was performed using four different strategies: a PGC cartridge^{19,22,23}, a reverse-phase cartridge^{24,25}, ultrafiltration²⁶, and precipitation²⁷. The results were evaluated based on the recovery rate of both neutral and acidic *N*-glycans. The highest recovery rate was achieved using the PGC method, and neutral and acidic *N*-glycans were both sensitive to the enrichment methods (Supplementary Figure S2). Moreover, the PGC method is more efficient than the other methods because salt, *N*-glycans, and de-glycosylated proteins can be sequentially eluted from the PGC cartridge²³. Given the labile property of sialylated *N*-glycan during the sample preparation, chromatographic separation and ionization in MS^{28,29}, the stability of di-sialylated *N*-glycan was determined in our experiment with varied concentration. The results suggested that native sialylated *N*-glycans are quite stable when using the method described. Thus the determined *N*-glycans represent a "genuine" and "primary" profile of rat serum (Data not show).

Separation of N-qlycans on the microfluidic PGC chip. When using routine mobile phases, i.e., 0.1% formic acid in water and acetonitrile (pH 3), highly sialylated N-glycans (containing three or four sialic acids) were out of the scope of analysis due to their poor peak shape³⁰, leading to the missed detection of important N-glycans. To improve the detection of these highly sialylated N-glycans, several mobile phases with varied pH values and ionic strengths were compared. As a result, 0.5% formic acid of pH 3.0 (adjusted by adding ammonium hydroxide) in water and acetonitrile was shown to be the optimal mobile phase for the separation of acidic N-glycans. A significantly improved peak shape for all acidic N-glycans, especially the highly sialylated N-glycans, was obtained using this mobile phase. In addition, the pH of this mobile phase allowed separation of acidic N-glycans from neutral N-glycans, thus avoiding ion suppression, which in turn enhanced the sensitivity of detection of acidic N-glycans. Furthermore, the improved peak shape further facilitated resolution of several groups of sialylated N-glycan isomers on the PGC nano-column. For example, the two peaks corresponding to isomers of Hex₅HexNAc₄NeuAc₂ were observed under the developed chromatographic conditions (Supplementary Figure S3). With the good separation of N-glycans, several consistent trends in the retention times of N-glycans on the PGC stationary phase were observed and are summarized in Supplementary Figures S4 and S5. These retention behaviors provided additional information for the assignment of N-glycans. In addition, the repeatability of N-glycan profiles as analyzed by our platform was high (Supplementary Figure S6). This study indicated that chip-based-LC/MS technology is a stable technology suitable for the comparative study of the N-glycan profiles of rat, mouse, and human sera.

Methods for the comprehensive characterization of rat serum *N***-glycans.** With the improved separation and detection, the raw data were analyzed on the basis of retention time, accurate mass, isotope abundance, and MS/MS fragmentation. In total, 282 *N*-glycans including 192 acidic *N*-glycans and 90 neutral *N*-glycans were identified, among which 172 *N*-glycans were reported from rat serum for the first time, and 27 were novel *N*-glycans (Fig. 1a and Supplementary Table S1). The analytical procedure is described step-by-step below.

Characterization of N-glycans based on accurate mass. For the characterization of N-glycans, a MassHunter Personal Compound Database and Library (PCDL) was created. This theoretical N-glycan database was established based on the knowledge of mammalian serum N-glycan biosynthesis and the possible modifications^{9,10,19,31}. As a result, a theoretical library of rat serum N-glycans consisting of more than 4000 N-glycan compositions was created for identification purposes.

For the determination of N-glycan composition, accurate mass and isotopic clusters for each compound were matched with the theoretical mass values of N-glycans in the PCDL. The ultra-high accurate mass value obtained by Q-TOF MS facilitated assignment of N-glycans at the MS level. In particular, isobaric N-glycans can be well differentiated based on the high-resolution MS data. For example, the N-glycan at m/z 1141.4015 was identified as disialylated N-glycan Hex₅HexNAc₄NeuAc₁NeuGc₁ + OAc based on its measured mass value ($C_{8c}H_{140}N_6O_{64}$, theoretical mass 2280.7885), while the N-glycan at m/z 1141.4197 was characterized as Hex₅HexNAc₅dHex₁NeuAc₁ on the basis of its measured mass value ($C_{87}H_{144}N_6O_{63}$, theoretical mass 2280.8247). Both N-glycans could be unambiguously distinguished according to their accurate mass values (Fig. 1b).

MS/MS screening for NeuAc- and NeuGc-containing N-glycans. After each N-glycan was assigned based on its accurate mass, the monosaccharide compositions were confirmed using targeted MS/MS experiments. N-glycolylneuraminic acid (NeuAc) and N-glycolylneuraminic acid (NeuGc) represent two of the most common sialic acids in nature. N-glycans containing NeuAc and NeuGc were widespread on all mammalian cell surfaces in a species-specific manner^{19,32-35}, and were associated with pathological conditions³⁶⁻³⁸. This indicates the importance of identification of NeuAc and NeuGc. For organisms expressing both NeuAc and NeuGc, such as rat, it is necessary to differentiate these two species on the basis of chemical evidence.

The only difference between the two sialic acids, NeuAc and NeuGc, is an additional oxygen atom in the *N*-glycolyl group of NeuGc (acetyl amino group)³⁹. Such subtle differences in the molecular composition led to the formation of a series of isomers of these two species. For example, the molecular formula of Hex₄HexNAc₃NeuAc₁ (C₅₉H₉₈N₄O₄₄) is exactly the same as that of Hex₃HexNAc₃dHex₁NeuGc₁ (C₅₉H₉₈N₄O₄₄). For differentiating these types of isomers, targeted MS/MS experiments were carried out. Fragment ions derived from respective sialic acids, *i.e.*, product ions at *m/z* 274.09 (NeuAc-H₂O), *m/z* 292.10 (NeuAc), and *m/z* 657.23 (Hex₁HexNAc₁NeuAc₁) generated from NeuAc-containing *N*-glycan, as well as fragment ions at *m/z* 290.08 (NeuGc-H₂O), *m/z* 308.09 (NeuGc), and *m/z* 673.22 (Hex₁HexNAc₁NeuGc₁) derived from NeuGc-containing *N*-glycan, were employed as diagnostic ions for differentiation of these two species (Fig. 2a). This approach facilitated the identification of a large number of sialylated *N*-glycans, among which 129 were NeuAc-containing *N*-glycans, 59 were NeuGc-containing *N*-glycans, and four were 'mixed' *N*-glycans containing both NeuAc and NeuGc (Fig. 2b). Among the four mixed *N*-glycans, *i.e.*, Hex₅HexNAc₄NeuAc₁NeuGc₁, Hex₅HexNAc₄NeuAc₁NeuGc₁ + OAc, Hex₅HexNAc₄NeuAc₁NeuGc₁, and Hex₅HexNAc₄dHex₁N euAc₁NeuGc₁ was found on bovine plasminogen and porcine




Figure 1. Overlaid extracted compound chromatograms (ECC) of rat serum N-glycans identified by PGC nano-LC/MS (a) and differentiation of isobaric N-glycans based on the accurate mass value (b).

vitronectin 40,41 and mice serum 42 . This is the first time that this type of sialylated N-glycan has been identified in rats.

Screening for modified N-glycans. Generally, additional modifications of N-glycans occur on the hydroxyl group of the monosaccharide, these may include O-acetylation, methylation, lactylation, sulfation, or phosphorylation⁴³⁻⁴⁵. These modifications result from the actions of different transferases and are monosaccharide specific. For example, mannose is often modified by phosphorylation^{44,46}, while galactose and N-acetylhexosamine are generally modified by sulfation^{45,47}. Notably, modification of sialic acid is not only the most abundant but also the most structurally diverse because all of the aforementioned modifications can occur on sialic acid⁴³.

To survey the modifications of rat serum N-glycans, the aforementioned five types of modifications were all included in the PCDL for characterization. On the basis of accurate mass, isotopic abundance, and targeted MS/MS, O-acetylation was found to be the dominant modification of rat serum N-glycans. O-acetylation occurred on the terminal sialic acid, which has only be reported in the serum of Atlantic salmon⁴⁸ and several glycoproteins of rats^{31,49}. Taking di-sialylated N-glycan $Hex_5HexNAc_4NeuAc_2$ and its mono-O-acetylation as an example, the identification procedures were performed as follows. First,




Figure 2. Targeted MS/MS of sialylated *N*-glycans in rat serum. (a) Characteristic fragment ions derived from *N*-glycans containing NeuAc (Left), NeuGc (Middle), and both NeuAc and NeuGc (Right). (b) Full MS/MS spectra of rat serum *N*-glycan Hex₅HexNAc₄NeuAc₁NeuGc₁ (doubly protonated *m/z* 1120.40).

N-glycans were identified based on PCDL. Because the accurate mass and isotopic abundance of compounds in the MS scan were highly similar to the theoretical data on PCDL, they were preliminarily identified as $\text{Hex}_5\text{HexNAc}_4\text{NeuAc}_2$ and $\text{Hex}_5\text{HexNAc}_4\text{NeuAc}_2 + \text{OAc}$. Then, O-acetylation was verified by analysis of an associated fragment produced by targeted MS/MS. Regular variations were found by comparing the MS/MS spectra of the two N-glycans (Fig. 3a). A series of fragments, such as m/z 274.09, 292.10, and 657.23, shown in the $\text{Hex}_5\text{HexNAc}_4\text{NeuAc}_2$ MS/MS spectrum (Fig. 3b) indicated the existence of $\text{HexNAc}_4\text{Hex} + \text{NeuAc}$, that is sequentially attached to the core structure of the N-glycan. For $\text{Hex}_5\text{HexNAc}_4\text{NeuAc}_2 + \text{OAc}$ (Fig. 3c), the series of fragments were 42 Da more than the corresponding fragments in $\text{Hex}_5\text{HexNAc}_4\text{NeuAc}_2$, *i.e.*, m/z 316.09 (274.09 + 42), 334.10 (292.10 + 42), and 699.24 (657.23 + 42). These characteristic ions facilitated identification of O-acetylated sialic acid in N-glycans. Using this approach, N-glycans $\text{Hex}_5\text{HexNAc}_4\text{NeuAc}_2$ with between one and three acetyl groups were identified (Fig. 3d,e).

It is worth mentioning that uncommon disialylated N-glycans with three O-acetyl groups, $Hex_5HexNAc_4NeuAc_2 + 3OAc$ and $Hex_5HexNAc_4dHex_1NeuAc_2 + 3OAc$, were identified in rat serum using this approach. In these N-glycans, one O-acetyl group was located at one terminal NeuAc, while the other two O-acetyl groups were located at another terminal NeuAc. This approach also led to the identification of a novel O-acetylated 'mixed' N-glycan possessing both NeuAc and NeuGc in the molecule, *i.e.*, $Hex_5HexNAc_4NeuAc_1NeuGc_1 + OAc$ (m/z 1141.40) (Fig. 3f). Acetylation-associated fragments are summarized in Table 1 and the characterized O-acetylated N-glycans of rat serum are listed in Table 2. It was noted that, although both NeuAc and NeuGc exist in rat serum N-glycans, O-acetylation exclusively occurs on NeuAc.

Comparison of the *N***-glycomic profiles of rat, mouse, and human sera.** According to different biosynthetic procedures, serum *N*-glycans were divided into five classes: high mannose (high Man) *N*-glycans; undecorated complex/hybrid (C/H) *N*-glycans; fucosylated complex/hybrid (C/H-F)




Figure 3. Targeted MS/MS of O-acetylated N-glycans in rat serum. (a) Diagnostic fragment ions generated from N-glycans with O-acetylated NeuAc (Left) and O-acetylated NeuGc (Right). (\mathbf{b} - \mathbf{e}) Full MS/MS spectra of a series of N-glycans with varied numbers (0–3) of acetylation sites on NeuAc and (\mathbf{f}) N-glycan Hex₅HexNAc₄NeuAc₁ NeuGc₁ + OAc.

N-glycans; sialylated complex/hybrid (C/H-S) N-glycans; and fucosylated-sialylated complex/hybrid (C/H-FS) N-glycans. Using the established glycomic approach, N-glycans of mouse and human sera were also comprehensively profiled and compared with N-glycans of rat serum (Fig. 4a). Relative abundances were assigned to each N-glycan class on the basis of the total ion counts for the individual N-glycan signals (Fig. 4b). The results showed that sialylated N-glycans were much more abundant than neutral N-glycans in all three species, accounting for about 70–80% of the total N-glycans in the serum of each species. Di-sialylated N-glycans were always the most dominant species, accounting for more than 50% of the total N-glycans, followed by mono-, tri-, and tetra-sialylated N-glycans (Fig. 5a). In addition, fuco-sylated N-glycans (25–40%) were significantly less abundant than their counterparts (60–75%) in the sera of all three species (Fig. 5b), and more than 80% of the N-glycans in rat, mouse, and human sera possessed di-antennary structures, whereas mono- and tri-antennary structures were less abundant (Fig. 5c).

Despite similarities in the total levels of *N*-glycans in sera among the three species, the differences in the composition of serum *N*-glycans were remarkable among rats, mice, and humans. First, the most notable difference in *N*-glycans among rat, mouse, and human sera was the type of sialic acid present (Fig. 6a). In humans, sialic acids in the serum *N*-glycans were exclusively in the form of NeuAc, whereas NeuGc was the sole sialic acid moiety in mouse serum. However, both NeuAc and NeuGc were identified in the *N*-glycans of rat serum. It is noted that the content of NeuAc in rat serum is almost 20 times that of NeuGc, which is quit similar to human glycome in which NeuAc is the solely acidic species. Second, a large number of sialylated *N*-glycans with *O*-acetylation were detected in rat serum, which, by contrast, were of low abundance in human serum and extremely low abundance in mouse serum (Fig. 6b). The distinct *O*-acetylation of rat glycome was structurally a kind of modification, and this modification exclusively occurred on NeuAc. Our results indicated a species-specific composition of sialic acid in the *N*-glycans of human, mouse and rat sera. However, the glycomic profile of rat serum was more similar to that of human serum based on both the high abundance of NeuAc and similar glycomic distributions (Fig. 4a).

Discussion

A comprehensive glycomic study of rat serum is not available in the literature. To date, less than 50 N-glycan compositions have been reported, significantly less than those identified in human sera^{50,51} and mouse sera¹⁹ (about 100 compositions for each species). Using optimized chromatographic conditions,

			Charge state	O-acetylation		
Composition	Symbol	m/z		Non-	Mono-	Di-
NeuAc - H ₂ O	→ -H ₂ O	274.09		+	+	+
NeuAc	•	292.10	1	1 +		+
NeuAc+OAc-H ₂ O	→ +OAc -H ₂ O	316.09	316.09 1		+	+
NeuAc + OAc	◆ +OAc	334.10	1		+	+
NeuAc + 2OAc - H ₂ O	+2OAc -H₂O	358.11	1			+
NeuAc+2OAc	→ +2OAc	376.12	1			+
Hex ₁ NeuAc ₁ - H ₂ O	◆-0 -H ₂ O	454.15	1 +		+	+
$Hex_1NeuAc_1 + OAc - H_2O$	♦•• +OAc-H ₂ O	496.17	1		+	+
$\text{Hex}_1\text{NeuAc}_1 + 2\text{OAc} - \text{H}_2\text{O}$	◆-○ +2OAc-H ₂ O	538.18	1			+
$Hex_1HexNAc_1NeuAc_1$	◆	657.23	1	+	+	+
$Hex_1HexNAc_1NeuAc_1 + OAc$	◆•••+OAc-H ₂ O	699.24	1		+	+
$Hex_1HexNAc_1NeuAc_1 + 2OAc$	◆	741.26	1			+
Hex ₂ HexNAc ₁ NeuAc ₁	◆	819.28	8 1		+	+
$Hex_2HexNAc_1NeuAc_1 + OAc$	◆•••+OAc-H ₂ O	861.30	1		+	+
$\text{Hex}_2\text{HexNAc}_1\text{NeuAc}_1 + 2\text{OAc}$	◆ • • +2OAc-H ₂ O	903.31	1			+
Hex ₅ HexNAc ₄ NeuAc ₁	-H₂O	966.85	2 +		+	+
Hex ₅ HexNAc ₄ NeuAc ₁ + OAc	+OAc-H ₂ O	987.85	2		+	+
Hex ₅ HexNAc ₄ NeuAc ₁ + 2OAc	+2OAc-H ₂ O	1008.86	2			+
Hex ₄ HexNAc ₂ NeuAc ₁	••••	1346.47	1 +		+	+
Hex ₄ HexNAc ₂ NeuAc ₁ + OAc	+OAc-H ₂ O	1388.48	1		+	+
Hex ₄ HexNAc ₂ NeuAc ₁ + 2OAc	+2OAc-H ₂ O	1430.49	1			+
Hex ₃ HexNAc ₃ NeuAc ₁	••••-H ₂ O	1405.50	1	+	+	+
Hex ₃ HexNAc ₃ NeuAc ₁ +OAc	+OAc-H ₂ O	1447.52	1		+	+
$Hex_3HexNAc_3NeuAc_1 + 2OAc$	◆◆■◆ **+2OAc-H ₂ O	1489.52	1			+
Hex ₄ HexNAc ₃ NeuAc ₁	-H ₂ O	1567.56	1 +		+	+
Hex ₄ HexNAc ₃ NeuAc ₁ + OAc	+OAc-H ₂ O	1609.57	1		+	+
Hex ₄ HexNAc ₃ NeuAc ₁ + 2OAc	+2OAc-H ₂ O	1651.58	1			+
		-1	1	1	I .	1

Table 1. Typical O-acetylated NeuAc-associated fragments. *the fragment ion was detected on MS/MS spectrum.

overlapping N-glycans were resolved as far as possible on the PGC-chip, and ionization suppression was significantly reduced, thereby allowing for the detection of low abundance N-glycans. Another advantage of the optimized chromatographic conditions was the improved separation of highly sialylated N-glycans. Fourteen tri- and two tetra-sialylated N-glycans that could not be detected using routine methods were fully characterized in this study. Collectively, our microfluidic chip-LC-MS-based approach enabled improved detection of both neutral and acidic N-glycans with highly enhanced sensitivity of the low abundance N-glycans. This approach resulted in the detection of 282 N-glycans, including 27 novel N-glycans in rat serum.

In addition to the overall enhanced sensitivity, employment of microfluidic chips also allowed the separation of isomeric/isobaric species. Efficient separation of the isomeric/isobaric N-glycans is crucial for the acquisition of the target MS/MS spectrum, which is in turn essential for the differentiation of these isomeric/isobaric species. In the current study, isomeric species of NeuAc and NeuGc were successfully discriminated in rat serum. This was a significant finding as it is the first report that both types of sialic acid exist in single N-glycans in this animal model.

Overall enhanced sensitivity, separation-facilitated reliable MS/MS experiments, together with annotation on the basis of well-established PCDL, further led to the discovery of 102 N-glycans containing O-acetylated sialic acid in rat serum. Of these, 27 O-acetylated N-glycans from 12 compositions were reported for the first time in rat serum. The extensive O-acetylation of N-glycans is clearly a distinguishing feature of rat serum, as the O-acetylation level of sialylated N-glycans in mouse and human sera is extremely low. O-acetylation, the most common modification of sialic acid⁵², has been suggested to be involved in a variety of biological processes in mammals, such as lectin recognition, virus binding, tumor antigenicity, tissue morphogenesis, and cell-cell interactions^{43,53}. The implications of the extensive O-acetylation of N-glycans detected in rat serum need to be explored in future studies.

The glycomic analysis carried out in this study provided a comprehensive *N*-glycan profile of rat serum. The results revealed a high abundance of acidic *N*-glycans containing NeuAc and/or NeuGc. The degree of sialylation varied from one to four, with di-sialylated *N*-glycans as the most dominant species. Significantly, structural analysis of the rat serum *N*-glycans revealed extensive *O*-acetylation of NeuAc, which was more diverse than any previous report^{30,48}. The degree of *O*-acetylation at distinct sites varied from one to three, with mono-acetylation being the most dominant species. To the best of our knowledge, this is the first report showing the extensive *O*-acetylation of the *N*-glycome of rat serum.

Animal models have been used increasingly for investigations in the field of glycobiology. Such investigations are based on the assumption that high similarity exists between the glycomes of mice or rats and humans, which has not yet been extensively studied. The aim of the current study was therefore to carry out a comprehensive comparison of the glycome of rat, mouse, and human sera. Similarities in the serum glycome of these three species were revealed based on the patterns of sialylation, fucosylation, and the degree of branching of *N*-glycans. This is ascribable to the fact that the biosyntheses of *N*-glycans are primarily regulated by hundreds of transferases encoded by the same genes in all three mammals^{54,55}.

However, there were some notable differences in the types of sialic acid modifications on the N-glycans of these three species, i.e., NeuAc, NeuGc, and O-acetylated NeuAc. These subtle differences between species may be explained by evolution and epigenetic regulation. First, because the CMAH gene mutation in the human genome occurred by a lethal microbial pathogen that required cell-surface NeuGc for effective infection, humans are now deficient in NeuGc biosynthesis⁵⁶, whereas such a mutation is rare in other vertebrates⁵⁷. Thus, species-specific sialic acids existed: NeuAc for human sera, NeuGc for mouse sera, and both NeuAc and NeuGc for rat sera. Second, glycosylation is not a simple addition to a protein, but an intricately and carefully regulated physiological process, which is affected by epigenetic regulation, which includes competing enzyme activities and the availability of activated monosaccharide donors⁵⁸⁻⁶⁰. Epigenetic regulation can modify biological structures and add/remove covalent groups, such as O-acetyl and methyl groups, which play important roles in the adaptation to specific environmental conditions^{61,62}. Among these modifications, O-acetylation plays a crucial role in the regulation of ligand function of terminal sialic acids on N-glycans⁶³. It is well established that the attachment of influenza A and B viruses can be inhibited by O-acetylation of sialylated ligands owing to reducing viral sialidase activity^{48,63}. Thus, it could be speculated that the unique and extensive O-acetylation of N-glycans in rat serum may be partially associated with the strong immune resistance of rats against various ailments⁶⁴. Despite the distinct features of rat serum glycome, it was more similar to human serum glycome than the mouse serum glycome. In conclusion, our results provide important glycomic data that will aid selection of the most appropriate animal model for pathological and pharmacological studies.

Methods

This study was approved by the Ethics Committee of the Guangdong Provincial Chinese Medicine Hospital. The methods were carried out in accordance with the approved guidelines.

N-glycan nomenclature and abbreviations. Nomenclature and cartoons of *N*-glycan species are presented in accordance with the Consortium for Functional Glycomics.

Materials and chemicals. Peptide *N*-glycosidase (PNGase F) was purchased from New England Biolabs (Ipswich, MA, USA). Dithiothreitol (DTT) was purchased from GE Healthcare (Uppsala, Sweden), HyperSep Hypercarb PGC was from Thermo Scientific (Runcorn WA7 1TA, UK). Sep-Pak Vac (1cc) was obtained from Waters (Milford, MA, USA). Centrifuge filters for ultrafiltration were a product of Millipore (Carrigtowohill, Ireland). Acetonitrile (LC/MS grade) and water were from Avantor (Center Valley, PA, USA), formic acid and ammonium bicarbonate (LC/MS grade) were products of Sigma-Aldrich (St. Louis, MO, USA). Reagent grade ammonia and acetone were obtained from Scharlab S. L (Sentmenat, Spain).

2° Hex,HenNac,NeuAc ₁ +2OAc 4.3.0.1° C63 H102 N4 O46 1650.5765 1650.5662 -6.27 3° Hex,HenNac,NeuAc ₂ + OAc 4.4.0.1° C69 H113 N5 O59 1811.6443 1811.6445 -0.48 7 4 Hex,HenNac,NeuAc ₂ + OAc 5.4.0.1° C75 H122 N5 O55 1973.6982 1973.6993 0.60 7.68.7 5° Hex,HenNac,NeuAc ₂ + OAc 5.3.0.2° C78 H127 N5 O58 2061.7142 2061.7202 -2.53 7° Hex,HenNac,Relack,NeuAc ₂ + OAc 5.3.0.2° C78 H127 N5 O58 2061.7142 2061.7202 -2.53 8° Hex,HenNac,Giller,NeuAc ₂ + OAc 5.4.1.1° C81 H133 N5 O59 2119.7561 2119.7595 -0.01 8.75 9 Hex,HenNac,Giller,NeuAc ₂ + OAc 5.4.0.2° C86 H140 N6 O63 2264.7936 2264.7934 -0.08 7.93.8 10° Hex,HenNac,NeuAc ₂ NeuGc ₂ + OAc 5.4.0.2° C86 H140 N6 O64 2280.8802 2280.872 -0.00 7.93.8 10° Hex,HenNac,NeuAc ₂ NeuGc ₂ + OAc 5.5.1.1° C88 H142 N6 O64 2306.8022	о.	Name	Abbr. N-glycan Composition	Molecular Formula	Calculated Mass	Observed Mass	Error (ppm)	t _R (min)	
38	:	Hex ₄ HexNAc ₃ NeuAc ₁ + OAc	4_3_0_1*a	C61 H100 N4O45	1608.5660	1608.5637	-1.40	8.40, 9.87	
4 Hes, HerNAc, NeuAc, + OAC 5,4,0,1° C75 H123 N5 O55 1973,6982 1973,6993 0.60 7.68,7 5° Hes, HerNAc, NeuAc, + OAC 4,5,0,1° C77 H126 N6 O55 2014.7247 2014.726 0.43 66 Hes, HerNAc, NeuAc, + OAC 5,3,0,2° C78 H127 N5 O58 2061.7142 2061.7202 -2.53 7° Hes, HerNAc, NeuAc, + OAC 5,4,1,1° C81 H133 N5 O59 2119.7561 2119.7595 -0.01 8.75 8° Hes, HerNAc, HeuAc, + OAC 5,4,1,1° C83 H135 N5 O60 2161.7666 2161.7697 1.42 8.8 9 Hes, HerNAc, NeuAc, + OAC 5,4,0,2° C86 H140 N6 O63 2264.7936 2264.7934 -0.08 7.93,8 10° Hes, HerNAc, NeuAc, + OAC 5,4,0,2° C86 H140 N6 O63 2264.7936 2264.7934 -0.50 11° Hes, HerNAc, NeuAc, + OAC 5,4,0,2° C88 H140 N6 O64 2230.7885 2280.7874 -0.50 11° Hes, HerNAc, NeuAc, + OAC 5,5,1,1° C89 H146 N6 O64 2302.8845 2322.4427 3.11 13° Hes, HerNAc, NeuAc, + OAC 5,5,1,1° C99 H146 N6 O65 2348.8147 2348.8135 -0.50 919.9 15° Hes, HerNAc, HeuAc, + OAC 5,5,1,1° C91 H148 N6 O65 2348.8460 2246.8431 -1.25 9.94 16° Hes, HerNAc, HerNAc, NeuAc, + OAC 5,5,1,1° C92 H150 N6 O67 2410.8515 2410.8511 -0.14 8.39,8 17° Hes, HerNAc, HerNAc, NeuAc, + OAC 5,4,1,2° C92 H150 N6 O67 2410.8515 2410.8511 -0.14 8.39,8 17° Hes, HerNAc, HerNAc, NeuAc, + OAC 5,4,1,2° C92 H150 N6 O67 2410.8515 2410.8511 -0.14 8.39,8 17° Hes, HerNAc, HerNAc, HerNAc, + OAC 5,4,1,2° C94 H152 N6 O68 2426.8464 2426.8349 -4.75 0.8 17° Hes, HerNAc, HerNAc, + OAC 5,4,1,2° C94 H152 N6 O68 2426.8464 2426.8349 -4.75 0.8 17° Hes, HerNAc, HerNAc, + OAC 5,4,1,2° C94 H152 N6 O68 2426.8464 2426.8349 -4.75 0.8 18° Hes, HerNAc, HerNAc, + OAC 5,4,1,2° C94 H152 N6 O68 2426.8464 2426.8349 -4.75 0.8 19° Hes, HerNAc, HerNAc, + OAC 5,4,0,3° C94 H152 N7 O75 2597.896 2597.901 0.57 0.57 0.57 0.57 0.57 0.57 0.57 0.57 0.57 0.57 0.57 0.57 0.57 0.57 0.57 0.57 0.57 0.57 0.5		Hex ₄ HexNAc ₃ NeuAc ₁ +2OAc	4_3_0_1*b	C63 H102 N4 O46	1650.5765	1650.5662	-6.27	8.14	
5° Hex,HexNAc,NeuAc, + OAc 4_5_0_1^{12} C77 H126 N6 O55 20147247 2014726 0.43 6 Hex,HexNAc,NeuAc_+ OAc 5_3_0_2^{24} C78 H127 N5 O58 2061.7142 2061.7202 -2.53 7° Hex,HexNAc,Gletex,NeuAc_+ OAc 5_4_1_1^{-16} C81 H133 N5 O59 2119.7561 2119.7595 -0.01 8.75 8° Hex,HexNAc,Gletax,NeuAc_+ 2OAc 5_4_0_2^{-16} C86 H140 N6 O63 2264.7936 2264.7934 -0.08 7.93.8 10° Hex,HexNAc_NeuAc_+ COAc 5_4_0_2^{-16} C86 H140 N6 O64 22280.7885 2280.7874 -0.50 11 Hex,HexNAc_NeuAc_+ COAc 5_4_0_2^{-16} C88 H142 N6 O64 2306.8042 2306.8022 -0.31 12° Hex,HexNAc_dHex,NeuAc_+ COAc 5_5_1_1^{-16} C88 H144 N6 O64 2322.8354 2322.8427 3.11 13 Hex,HexNAc_dHex,NeuAc_+ COAc 5_5_1_1^{-16} C89 H144 N6 O65 2348.8147 2348.8135 -0.50 9.19.9 14 Hex,HexNAc_dHex,NeuAc_+ COAc 5_5_1_1^{-2} C92 H150 N6 O67 2410.8515 2410.8511 </td <td></td> <td>$Hex_4HexNAc_4NeuAc_1 + OAc$</td> <td>4_4_0_1*a</td> <td>C69 H113 N5 O50</td> <td>1811.6453</td> <td>1811.6445</td> <td>-0.48</td> <td colspan="2">7.68, 8.42</td>		$Hex_4HexNAc_4NeuAc_1 + OAc$	4_4_0_1*a	C69 H113 N5 O50	1811.6453	1811.6445	-0.48	7.68, 8.42	
6 Hex,HexNAc,NeuAc; +OAc 5_4_1_1** C81 H133 N5 O59 2119.7561 2119.7595 -0.01 8.75 8* Hex,HexNAc,dHex,NeuAc; +OAc 5_4_1_1** C81 H133 N5 O59 2119.7561 2119.7595 -0.01 8.75 8* Hex,HexNAc,dHex,NeuAc; +OAc 5_4_1_1** C81 H133 N5 O59 2119.7561 2119.7595 -0.01 8.75 8* Hex,HexNAc,dHex,NeuAc; +OAc 5_4_0_2** C86 H140 N6 O63 2264.7936 2264.7934 -0.08 7.93, 8 10* Hex,HexNAc,NeuAc; +OAc 5_4_0_2** C86 H140 N6 O63 2264.7936 2264.7934 -0.08 7.93, 8 10* Hex,HexNAc,NeuAc; +OAc 5_4_0_2** C86 H140 N6 O64 2280.7885 2280.7874 -0.50 11* Hex,HexNAc,NeuAc; +OAc 5_4_0_2** C86 H140 N6 O64 2306.8042 2306.8022 -0.83 8.34, 8 12* Hex,HexNAc,NeuAc; +OAc 5_5_1_1** C89 H146 N6 O64 2322.8354 2322.8427 3.11 13* Hex,HexNAc,NeuAc; +OAc 5_5_1_1** C89 H146 N6 O64 2322.8354 2322.8427 3.11 14* Hex,HexNAc,NeuAc; +OAc 5_5_1_1** C91 H148 N6 O65 2348.8147 2348.8135 -0.50 9.19, 9 10.04 14* Hex,HexNAc,dHex,NeuAc; +OAc 5_5_1_1** C91 H148 N6 O65 2364.8460 2364.8431 -1.25 9.2 15* Hex,HexNAc,dHex,NeuAc; +OAc 5_4_1_2** C92 H150 N6 O67 2410.8515 2410.8511 -0.14 8.839, 8 17* Hex,HexNAc,dHex,NeuAc; +OAc 5_4_1_2** C92 H150 N6 O68 2426.8464 2426.8349 -4.75 8.8 18* Hex,HexNAc,dHex,NeuAc; +OAc 5_4_1_2** C92 H150 N6 O68 2467.8729 2467.8555 -7.06 19* Hex,HexNAc,dHex,NeuAc; +OAc 5_4_1_2** C94 H153 N7 O68 2467.8729 2467.8555 -7.06 19* Hex,HexNAc,dHex,NeuAc; +OAc 5_4_1_2** C94 H153 N7 O68 2467.8729 2467.8555 -7.06 19* Hex,HexNAc,dHex,NeuAc; +OAc 5_4_0.3** C97 H157 N7 O71 2555.8890 2555.8877 -0.51 8.85, 9. 21 Hex,HexNAc,NeuAc; +OAc 5_4_0.3** C97 H157 N7 O71 2555.8890 2555.8877 -0.51 8.85, 9. 22 Hex,HexNAc,dHex,NeuAc; +OAc 5_4_1_3** C103 H167 N7 O75 2701.9469 2701.9497 1.02 9.85, 10* Hex,HexNAc,dHex,NeuAc; +OAc 5_4_1_3** C103 H167 N7 O75 2701.9469 2701.9497 1.02 9.85, 10* Hex,HexNAc,dHex,NeuAc; +OAc 5_4_1_3** C103 H167 N7 O75 2701.9469 2701.9477 1.02 9.85, 10* Hex,HexNAc,dHex,NeuAc; +OAc 6_5_0_3** C105 H169 N7 O76 2743.9575 2743.9537 -1.39 10.07, 10* Hex,HexNAc,NeuAc; +OAc 6_5_0_3** C105 H169 N7 O76 2743.9575 2743.9537 -1.39 10.07,		$Hex_5HexNAc_4NeuAc_1 + OAc$	5_4_0_1*a	C75 H123 N5 O55	1973.6982	1973.6993	0.60	7.68, 7.98, 8.72, 9.40	
7% Hex,HexNAc,dHex,NeuAc, + OAc 5_4_1_1^3 C81 H133 N5 O59 2119.7561 2119.7595 -0.01 8.75 8% Hex,HexNAc,dHex,NeuAc, + 2OAc 5_4_1_1^* C83 H135 N5 O60 2161.7666 2161.7697 1.42 8 9 Hex,HexNAc,NeuAc,+OAC 5_4_0_2^* C86 H140 N6 O63 2264.7936 2264.7934 -0.08 7.93.8 108 Hex,HexNAc,NeuAc,+OAC 5_4_0_2^* C86 H140 N6 O64 2280.7885 2280.7874 -0.50 11 Hex,HexNAc,NeuAc,+OAC 5_4_0_2^* C88 H142 N6 O64 2306.8042 2306.8022 -0.83 8.34.8 128 Hex,HexNAc,NeuAc,+OAC 5_5_1_1^* C89 H146 N6 O64 2308.8042 2306.8022 -0.83 8.34.8 128 Hex,HexNAc,NeuAc,+OAC 5_5_1_1^* C89 H146 N6 O64 2322.8354 2322.8427 3.11 13 Hex,HexNAc,NeuAc,+OAC 5_5_1_1^* C91 H148 N6 O65 2348.8147 2348.8135 -0.50 919.9 14 Hex,HexNAc,NeuAc,+OAC 5_5_1_2^* C92 H150 N6 O67 2410.8515 2410.8		$Hex_4HexNAc_5NeuAc_1 + OAc$	4_5_0_1*a	C77 H126 N6 O55	2014.7247	2014.726	0.43	9.75	
8		$Hex_5HexNAc_3NeuAc_2 + OAc$	5_3_0_2*a	C78 H127 N5 O58	2061.7142	2061.7202	-2.53	6.65	
9 Hes,HexNAc,NeuAc, + OAC 5.4.0.2° C86 H140 N6 O63 2264.7936 2264.7934 -0.08 7.93, 8 10 ⁸ Hex,HexNAc,NeuAc,NeuAc, + OAC 5.4.0.2° C86 H140 N6 O64 2280.7885 2280.7874 -0.50 11 Hex,HexNAc,NeuAc, + OAC 5.4.0.2° C88 H142 N6 O64 2306.8042 2306.8042 2306.8022 -0.83 8.34, 8. 12 ⁸ Hex,HexNAc,dHex,NeuAc, + OAC 5.5.1.1° C89 H146 N6 O64 2322.8354 2322.8427 3.11 13 Hex,HexNAc,dHex,NeuAc, + 3OAC 5.4.0.2° C90 H144 N6 O65 2348.8147 2348.8135 -0.50 910.04 14 Hex,HexNAc,dHex,NeuAc, + 2OAC 5.5.1.1° C91 H148 N6 O65 2364.8460 2364.8431 -1.25 9. 15 Hex,HexNAc,dHex,NeuAc, + OAC 5.4.1.2° C92 H150 N6 O67 2410.8515 2410.8511 -0.14 89.70.09 16 ⁸ Hex,HexNAc,dHex,NeuAc, + OAC 5.4.1.2° C92 H150 N6 O68 2426.8464 2426.8349 -4.75 88 17 ⁸ Hex,HexNAc,dHex,NeuAc, + OAC 5.4.1.2° C94 H153 N7 O68 2467.8729 2467.8555 -7.06 19 ⁸ Hex,HexNAc,dHex,NeuAc, + OAC 5.4.1.2° C96 H154 N6 O69 2494.8926 2494.8699 -1.10 9.82, 10. 20 Hex,HexNAc,NeuAc, + OAC 5.4.0.3° C99 H159 N7 O72 2597.8996 2597.9010 0.57 9.85, 9. 21 Hex,HexNAc,NeuAc, + OAC 5.4.0.3° C99 H159 N7 O72 2597.8996 2597.9010 0.57 9.05, 9.7 11 Hex,HexNAc,dHex,NeuAc, + OAC 5.4.0.3° C101 H161 N7 O73 2639.9101 2639.9032 -2.61 9.94, 10. 22 Hex,HexNAc,dHex,NeuAc, + OAC 5.4.0.3° C101 H161 N7 O73 2639.9101 2639.9032 -2.61 9.94, 10. 23 Hex,HexNAc,dHex,NeuAc, + OAC 5.4.1.3° C103 H167 N7 O75 2701.9469 2701.9497 1.02 9.85, 11.00 1.00 1.00 1.00 1.00 1.00 1.00 1.		$Hex_5 HexNAc_4 dHex_1 NeuAc_1 + OAc \\$	5_4_1_1*a	C81 H133 N5 O59	2119.7561	2119.7595	-0.01	8.75, 9.35, 9.87	
108		$\rm Hex_5 HexNAc_4 dHex_1 NeuAc_1 + 2OAc$	5_4_1_1*b	C83 H135 N5 O60	2161.7666	2161.7697	1.42	8.60, 9.20	
11 Hes,HexNAc,NeuAc, + 2OAc 5,4,0,2° C88 H142 N6 O64 2306.8042 2306.8042 -0.83 8.34, 8. 12* Hes,HexNAc,GelHes,NeuAc, + OAc 5,5,1,1° C89 H146 N6 O64 2322.8354 2322.8427 3.11 13 Hes,HexNAc,QeuAc, + 3OAc 5,4,0,2° C90 H144 N6 O65 2348.8147 2348.8135 -0.50 919.9 910.04 14 Hes,HexNAc,GelHes,NeuAc, + 2OAc 5,5,1,1° C91 H148 N6 O65 2364.8460 2364.8431 -1.25 9. 15 Hes,HexNAc,GelHes,NeuAc, + OAc 5,4,1,2° C92 H150 N6 O67 2410.8515 2410.8511 -0.14 8.39,8 8.39,8 8.39,8 9.70,9,3 16* Hes,HexNAc,GelHes,NeuAc,PeuAc, + OAc 5,4,1,2° C92 H150 N6 O68 2426.8464 2426.8349 -4.75 8.8 17* Hes,HexNAc,GelHes,NeuAc, + 2OAc 5,4,1,2° C94 H152 N6 O68 2452.8620 2452.8623 0.12 9.35,9 10 18* Hes,HexNAc,GelHes,NeuAc, + OAc 5,5,0,2° C94 H152 N6 O68 2467.8729 2467.8555 -7.06 19* Hex,HexNAc,GelHex,NeuAc, + 3OAc 5,4,1,2° C94 H153 N7 O68 2467.8729 2467.8555 -7.06 19* Hex,HexNAc,GelHex,NeuAc, + 3OAc 5,4,1,2° C96 H154 N6 O69 2494.8926 2494.8699 -1.10 9.82,100 20 Hex,HexNAc,GelHex,NeuAc, + 3OAc 5,4,0,3° C97 H157 N7 O71 2555.8890 2555.8877 -0.51 8.85, 9. 21 Hex,HexNAc,NeuAc, + 2OAc 5,4,0,3° C99 H159 N7 O72 2597.8996 2597.9010 0.57 9.05, 97.11 11 11 11 11 11 11 11 11 11 11 11 11		$Hex_5 Hex NAc_4 Neu Ac_2 + OAc$	5_4_0_2*a	C86 H140 N6 O63	2264.7936	2264.7934	-0.08	7.93, 8.47, 8.80, 9.12, 9.40	
128 Hex,HexNAc,dHex,NeuAc,+ OAc 5_5_1_1^2 C89 H146 N6 O64 2322.8354 2322.8427 3.11 Hex,HexNAc,dNeuAc,+ 3OAc 5_4_0_2^c C90 H144 N6 O65 2348.8147 2348.8135 -0.50 9,19,9 10.04 Hex,HexNAc,dHex,NeuAc,+ 2OAc 5_4_1_2^a C91 H148 N6 O65 2348.8147 2348.8135 -0.50 9,19,9 10.04 Hex,HexNAc,dHex,NeuAc,+ 2OAc 5_4_1_2^a C92 H150 N6 O67 2410.8515 2410.8511 -0.14 8.39,8 9.30,9 9.30,	&	$Hex_5HexNAc_4NeuAc_1NeuGc_1 + OAc$	5_4_0_2*#a	C86 H140 N6 O64	2280.7885	2280.7874	-0.50	9.02	
13		Hex ₅ HexNAc ₄ NeuAc ₂ +2OAc	5_4_0_2*b	C88 H142 N6 O64	2306.8042	2306.8022	-0.83	8.34, 8.94, 9.52, 10.12	
13 Hex ₃ HexNAc ₄ NeuAc ₁ + 2OAc 5_4_1_2* C91 H148 N6 O65 2364.8460 2364.8431 -1.25 9. 15 Hex ₃ HexNAc ₄ dHex ₁ NeuAc ₂ + OAc 5_4_1_2* C92 H150 N6 O67 2410.8515 2410.8511 -0.14 8.39, 8.70, 9.1 16* Hex ₃ HexNAc ₄ dHex ₁ NeuAc ₁ NeuGc ₁ + OAc 5_4_1_2** C92 H150 N6 O68 2426.8464 2426.8349 -4.75 8.70, 9.1 17* Hex ₃ HexNAc ₄ dHex ₁ NeuAc ₂ + 2OAc 5_4_1_2** C92 H150 N6 O68 2426.8464 2426.8349 -4.75 8.70, 9.1 18* Hex ₃ HexNAc ₄ dHex ₁ NeuAc ₂ + OAc 5_4_1_2** C94 H152 N6 O68 2452.8620 2452.8623 0.12 9.35, 9.1 18* Hex ₃ HexNAc ₅ NeuAc ₂ + OAc 5_5_0_2** C94 H153 N7 O68 2467.8729 2467.8755 -7.06 19* Hex ₃ HexNAc ₄ dHex ₁ NeuAc ₂ + 3OAc 5_4_1_2** C96 H154 N6 O69 2494.8926 2494.8699 -1.10 9.82, 10.1 20 Hex ₃ HexNAc ₄ NeuAc ₃ + OAc 5_4_0_3** C97 H157 N7 O71 2555.8890 2555.8877 -0.51 8.85, 9. 21 Hex ₃ HexNAc ₄ NeuAc ₃ + 2OAc 5_4_0_3** C99 H159 N7 O72 2597.8996 2597.9010 0.57 9.05, 9.7 11* Hex ₃ HexNAc ₄ NeuAc ₃ + 3OAc 5_4_0_3** C101 H161 N7 O73 2639.9101 2639.9032 -2.61 9.94, 10. 23 Hex ₃ HexNAc ₄ dHex1NeuAc ₃ + OAc 5_4_1_3** C103 H167 N7 O75 2701.9469 2701.9497 1.02 9.85, 10.07 24 Hex ₃ HexNAc ₄ dHex1NeuAc ₃ + 2OAc 5_4_1_3** C105 H169 N7 O76 2743.9575 2743.9537 -1.39 10.07, 10.07 25 Hex ₃ HexNAc ₄ dHex1NeuAc ₃ + 3OAc 5_4_1_3** C105 H169 N7 O76 2785.9680 2785.9679 -0.04 10.32, 11.1 26 Hex ₅ HexNAc ₄ dHex1NeuAc ₃ + OAc 6_5_0_3** C111 H180 N8 O81 2921.0212 2921.0160 -1.77 10.01.0.	&	$Hex_5 HexNAc_5 dHex_1 NeuAc_1 + OAc \\$	5_5_1_1*a	C89 H146 N6 O64	2322.8354	2322.8427	3.11	9.15	
15		$Hex_5 Hex NAc_4 Neu Ac_2 + 3 OAc$	5_4_0_2*c	C90 H144 N6 O65	2348.8147	2348.8135	-0.50	9.19, 9.39, 9.64, 9.80, 10.04, 10.22, 10.41	
16* Hex ₃ HexNAc ₄ dHex ₁ NeuAc ₂ +OAc 5_4_1_2^* C92 H150 N6 O68 2426.8464 2426.8349 -4.75 88 17* Hex ₃ HexNAc ₄ dHex1NeuAc ₂ +2OAc 5_4_1_2^* C92 H150 N6 O68 2426.8464 2426.8349 -4.75 88 17* Hex ₃ HexNAc ₄ dHex1NeuAc ₂ +2OAc 5_4_1_2^* C94 H152 N6 O68 2452.8620 2452.8623 0.12 9.35, 9.10 18* Hex ₃ HexNAc ₄ NeuAc ₂ +OAc 5_5_0_2^* C94 H153 N7 O68 2467.8729 2467.8555 -7.06 19* Hex ₃ HexNAc ₄ dHex ₁ NeuAc ₂ +3OAc 5_4_1_2^* C96 H154 N6 O69 2494.8926 2494.8699 -1.10 9.82, 10. 20 Hex ₃ HexNAc ₄ NeuAc ₃ +OAc 5_4_0_3^* C97 H157 N7 O71 2555.8890 2555.8877 -0.51 8.85, 9. 21 Hex ₃ HexNAc ₄ NeuAc ₃ +2OAc 5_4_0_3^* C99 H159 N7 O72 2597.8996 2597.9010 0.57 9.05, 9.7 22 Hex ₃ HexNAc ₄ NeuAc ₃ +3OAc 5_4_0_3^* C101 H161 N7 O73 2639.9101 2639.9032 -2.61 9.94, 10. 23 Hex ₃ HexNAc ₄ dHex1NeuAc ₃ +OAc 5_4_1_3^* C103 H167 N7 O75 2701.9469 2701.9497 1.02 9.85, 24 Hex ₃ HexNAc ₄ dHex1NeuAc ₃ +2OAc 5_4_1_3^* C105 H169 N7 O76 2743.9575 2743.9537 -1.39 10.07, 1		$Hex_5 HexNAc_5 dHex_1 NeuAc_1 + 2OAc \\$	5_5_1_1*b	C91 H148 N6 O65	2364.8460	2364.8431	-1.25	9.48, 10.10,	
17% Hex ₃ HexNAc ₄ dHex1NeuAc ₂ +2OAc 5_4_1_2'* C94 H152 N6 O68 2452.8620 2452.8623 0.12 9.35, 9. 10 18% Hex ₃ HexNAc ₅ NeuAc ₂ + OAc 5_5_0_2'* C94 H153 N7 O68 2467.8729 2467.8555 -7.06 19% Hex ₃ HexNAc ₄ dHex ₁ NeuAc ₂ +3OAc 5_4_1_2'* C96 H154 N6 O69 2494.8926 2494.8699 -1.10 9.82, 10. 20 Hex ₃ HexNAc ₄ NeuAc ₃ + OAc 5_4_0_3'* C97 H157 N7 O71 2555.8890 2555.8877 -0.51 8.85, 9. 21 Hex ₃ HexNAc ₄ NeuAc ₃ +2OAc 5_4_0_3'* C99 H159 N7 O72 2597.8996 2597.9010 0.57 9.05, 9.7 11 22 Hex ₃ HexNAc ₄ NeuAc ₃ +3OAc 5_4_0_3'* C101 H161 N7 O73 2639.9101 2639.9032 -2.61 9.94, 10. 23 Hex ₃ HexNAc ₄ dHex1NeuAc ₃ +OAc 5_4_1_3'* C103 H167 N7 O75 2701.9469 2701.9497 1.02 9.85, 10.07, 1		$Hex_5 Hex NAc_4 dHex_1 Neu Ac_2 + OAc \\$	5_4_1_2*a	C92 H150 N6 O67	2410.8515	2410.8511	-0.14	8.39, 8.69, 8.97, 9.20, 9.70, 9.85, 10.37, 10.64	
17° Hex ₅ HexNAc ₄ drexINeuAc ₂ +2OAc 5_4_1_2° C94 H153 N7 O68 2467.8729 2467.8555 -7.06 19 ⁸ Hex ₅ HexNAc ₄ dHex ₁ NeuAc ₂ +3OAc 5_4_1_2° C96 H154 N6 O69 2494.8926 2494.8699 -1.10 9.82, 10.0 20 Hex ₅ HexNAc ₄ NeuAc ₃ +OAc 5_4_0_3° C97 H157 N7 O71 2555.8890 2555.8877 -0.51 8.85, 9. 21 Hex ₅ HexNAc ₄ NeuAc ₃ +2OAc 5_4_0_3° C99 H159 N7 O72 2597.8996 2597.9010 0.57 9.05, 9.7 11 22 Hex ₅ HexNAc ₄ NeuAc ₃ +3OAc 5_4_0_3° C101 H161 N7 O73 2639.9101 2639.9032 -2.61 9.94, 10. 23 Hex ₅ HexNAc ₄ dHex1NeuAc ₃ +OAc 5_4_1_3° C103 H167 N7 O75 2701.9469 2701.9497 1.02 9.85, 24 Hex ₅ HexNAc ₄ dHex1NeuAc ₃ +2OAc 5_4_1_3° C105 H169 N7 O76 2743.9575 2743.9537 -1.39 10.07, 10 25 Hex ₅ HexNAc ₄ dHex1NeuAc ₃ +3OAc 5_4_1_3° C107 H171 N7 O77 2785.9680 2785.9679 -0.04 10.32, 26 Hex ₆ HexNAc ₅ NeuAc ₅ +OAc 6_5_0_3° C111 H180 N8 O81 2921.0212 2921.0160 -1.77 10.0, 10.	&	$Hex_5 HexNAc_4 dHex_1 NeuAc_1 NeuGc_1 + OAc \\$	5_4_1_2*#a	C92 H150 N6 O68	2426.8464	2426.8349	-4.75	8.67, 8.89	
19 ^{&} Hex ₃ HexNAc ₄ Hex ₁ NeuAc ₂ +3OAc 5 ₋ 4 ₋ 1 ₋ 2' ⁻ C96 H154 N6 O69 2494.8926 2494.8699 -1.10 9.82, 10.0 20 Hex ₃ HexNAc ₄ NeuAc ₃ +OAc 5 ₋ 4 ₋ 0 ₋ 3' ⁻ C97 H157 N7 O71 2555.8890 2555.8877 -0.51 8.85, 9. 21 Hex ₃ HexNAc ₄ NeuAc ₃ +2OAc 5 ₋ 4 ₋ 0 ₋ 3' ⁻ C99 H159 N7 O72 2597.8996 2597.9010 0.57 9.05, 9.7 11 12 Hex ₃ HexNAc ₄ NeuAc ₃ +3OAc 5 ₋ 4 ₋ 0 ₋ 3' ⁻ C101 H161 N7 O73 2639.9101 2639.9032 -2.61 9.94, 10. 22 Hex ₃ HexNAc ₄ NeuAc ₃ +OAc 5 ₋ 4 ₋ 1 ₋ 3' ⁻ C103 H167 N7 O75 2701.9469 2701.9497 1.02 9.85, 23 Hex ₃ HexNAc ₄ dHex1NeuAc ₃ +OAc 5 ₋ 4 ₋ 1 ₋ 3' ⁻ C105 H169 N7 O76 2743.9575 2743.9537 -1.39 10.07, 24 Hex ₃ HexNAc ₄ dHex1NeuAc ₃ +2OAc 5 ₋ 4 ₋ 1 ₋ 3' ⁻ C107 H171 N7 O77 2785.9680 2785.9679 -0.04 10.32, 25 Hex ₃ HexNAc ₄ dHex1NeuAc ₃ +OAc 6 ₋ 5 ₋ 0 ₋ 3' ⁻ C111 H180 N8 O81 2921.0212 2921.0160 -1.77 10.0, 10.	&	Hex ₅ HexNAc ₄ dHex1NeuAc ₂ +2OAc	5_4_1_2*b	C94 H152 N6 O68	2452.8620	2452.8623	0.12	9.35, 9.57, 9.94, 10.42, 10.84, 11.59	
19 Hex ₅ HexNAc ₄ NeuAc ₃ + OAc 5_4_0_3'a C97 H157 N7 O71 2555.8890 2555.8877 -0.51 8.85, 9. 20 Hex ₅ HexNAc ₄ NeuAc ₃ + OAc 5_4_0_3'a C97 H157 N7 O71 2555.8890 2555.8877 -0.51 8.85, 9. 21 Hex ₅ HexNAc ₄ NeuAc ₃ + 2OAc 5_4_0_3'b C99 H159 N7 O72 2597.8996 2597.9010 0.57 9.05, 9.7 11 22 Hex ₅ HexNAc ₄ NeuAc ₃ + 3OAc 5_4_0_3'c C101 H161 N7 O73 2639.9101 2639.9032 -2.61 9.94, 10. 23 Hex ₅ HexNAc ₄ dHex1NeuAc ₃ + OAc 5_4_1_3'a C103 H167 N7 O75 2701.9469 2701.9497 1.02 9.85, 24 Hex ₅ HexNAc ₄ dHex1NeuAc ₃ + 2OAc 5_4_1_3'b C105 H169 N7 O76 2743.9575 2743.9537 -1.39 10.07, 10 25 Hex ₅ HexNAc ₄ dHex1NeuAc ₃ + 3OAc 5_4_1_3'c C107 H171 N7 O77 2785.9680 2785.9679 -0.04 10.32, 11 26 Hex ₆ HexNAc ₅ NeuAc ₅ + 2OAc 6_5_0_3'a C111 H180 N8 O81 2921.0212 2921.0160 -1.77 10.0, 10.	&	Hex ₅ HexNAc ₅ NeuAc ₂ + OAc	5_5_0_2*a	C94 H153 N7 O68	2467.8729	2467.8555	-7.06	11.02	
21 Hex ₅ HexNAc ₄ NeuAc ₃ + 2OAc 5_4_0_3'b C99 H159 N7 O72 2597.8996 2597.9010 0.57 9.05, 9.7 11 22 Hex ₅ HexNAc ₄ NeuAc ₃ + 3OAc 5_4_0_3'c C101 H161 N7 O73 2639.9101 2639.9032 -2.61 9.94, 10. 23 Hex ₅ HexNAc ₄ dHex1NeuAc ₃ + OAc 5_4_1_3'a C103 H167 N7 O75 2701.9469 2701.9497 1.02 9.85, 24 Hex ₅ HexNAc ₄ dHex1NeuAc ₃ + 2OAc 5_4_1_3'b C105 H169 N7 O76 2743.9575 2743.9537 -1.39 10.07, 10 25 Hex ₅ HexNAc ₄ dHex1NeuAc ₃ + 3OAc 5_4_1_3'c C107 H171 N7 O77 2785.9680 2785.9679 -0.04 10.32, 11 26 Hex ₆ HexNAc ₅ NeuAc ₃ + OAc 6_5_0_3'a C111 H180 N8 O81 2921.0212 2921.0160 -1.77 10.0, 10.	&	Hex ₅ HexNAc ₄ dHex ₁ NeuAc ₂ +3OAc	5_4_1_2*c	C96 H154 N6 O69	2494.8926	2494.8699	-1.10	9.82, 10.00, 10.31, 10.49, 10.67	
21 Hex ₅ HexNAc ₄ NeuAc ₃ + 2OAc 5_4_0_3* C99 H159 N/ O/2 2597.8996 2597.9010 0.57 11 22 Hex ₅ HexNAc ₄ NeuAc ₃ + 3OAc 5_4_0_3*c C101 H161 N7 O73 2639.9101 2639.9032 -2.61 9.94, 10. 23 Hex ₅ HexNAc ₄ dHex1NeuAc ₃ + OAc 5_4_1_3*a C103 H167 N7 O75 2701.9469 2701.9497 1.02 9.85, 24 Hex ₅ HexNAc ₄ dHex1NeuAc ₃ + 2OAc 5_4_1_3*b C105 H169 N7 O76 2743.9575 2743.9537 -1.39 10.07, 10 25 Hex ₅ HexNAc ₄ dHex1NeuAc ₃ + 3OAc 5_4_1_3*c C107 H171 N7 O77 2785.9680 2785.9679 -0.04 10.32, 11 26 Hex ₆ HexNAc ₅ NeuAc ₃ + OAc 6_5_0_3*a C111 H180 N8 O81 2921.0212 2921.0160 -1.77 10.0, 10.		$Hex_5HexNAc_4NeuAc_3 + OAc$	5_4_0_3*a	C97 H157 N7 O71	2555.8890	2555.8877	-0.51	8.85, 9.25, 9.75, 10.00, 10.46	
23 Hex ₅ HexNAc ₄ dHex1NeuAc ₃ +OAc 5_4_1_3'a C103 H167 N7 O75 2701.9469 2701.9497 1.02 9.85, 24 Hex ₅ HexNAc ₄ dHex1NeuAc ₃ +2OAc 5_4_1_3'b C105 H169 N7 O76 2743.9575 2743.9537 -1.39 10.07, 10 10.07, 10 10.32, 11 Hex ₅ HexNAc ₄ dHex1NeuAc ₃ +3OAc 5_4_1_3'c C107 H171 N7 O77 2785.9680 2785.9679 -0.04 10.32, 11 10 10.07, 11 10		Hex ₅ HexNAc ₄ NeuAc ₃ + 2OAc	5_4_0_3*b	C99 H159 N7 O72	2597.8996	2597.9010	0.57	9.05, 9.77, 10.19, 10.57, 11.19, 11.77	
24 Hex ₅ HexNAc ₄ dHex1NeuAc ₃ +2OAc 5_4_1_3*b C105 H169 N7 O76 2743.9575 2743.9537 -1.39 10.07, 10 25 Hex ₅ HexNAc ₄ dHex1NeuAc ₃ +3OAc 5_4_1_3*c C107 H171 N7 O77 2785.9680 2785.9679 -0.04 10.32, 11 26 Hex ₆ HexNAc ₅ NeuAc ₃ +OAc 6_5_0_3*a C111 H180 N8 O81 2921.0212 2921.0160 -1.77 10.0, 10.		Hex ₅ HexNAc ₄ NeuAc ₃ + 3OAc	5_4_0_3*c	C101 H161 N7 O73	2639.9101	2639.9032	-2.61	9.94, 10.27, 10.79, 11.31, 12.01	
25 Hex ₅ HexNAc ₄ dHex1NeuAc ₃ +2OAc 5_4_1_3' C105 H169 N/ O/6 2/43.95/5 2/43.95/5 -1.39 10 25 Hex ₅ HexNAc ₄ dHex1NeuAc ₃ +3OAc 5_4_1_3' C107 H171 N7 O77 2785.9680 2785.9679 -0.04 10.32, 11 26 Hex ₆ HexNAc ₅ NeuAc ₃ +OAc 6_5_0_3' C111 H180 N8 O81 2921.0212 2921.0160 -1.77 10.0, 10.		$Hex_5 Hex NAc_4 dHex \\ 1 Neu Ac_3 + OAc$	5_4_1_3*a	C103 H167 N7 O75	2701.9469	2701.9497	1.02	9.85, 10.22, 10.84	
26 Hex ₆ HexNAc ₅ NeuAc ₃ + OAc 6_5_0_3*a C111 H180 N8 O81 2921.0212 2921.0160 -1.77 10.0, 10. 27 Hex Hex NAc NeuAc ₃ + 2OAc 6_5_0_3*b C113 H182 N8 O82 2963.0318 2963.0250 -2.28 10.32,		$Hex_5 Hex NAc_4 dHex 1 Neu Ac_3 + 2 OAc \\$	5_4_1_3*b	C105 H169 N7 O76	2743.9575	2743.9537	-1.39	10.07, 10.32, 10.51, 10.96, 11.47	
26 Hex ₆ HexNAc ₅ NetIAc ₃ + OAc 6_5_0_3 ** C111 H180 N8 O81 2921.0212 2921.0160 -1.//		Hex ₅ HexNAc ₄ dHex1NeuAc ₃ + 3OAc	5_4_1_3*c	C107 H171 N7 O77	2785.9680	2785.9679	-0.04	10.32, 10.79, 11.11, 11.59, 12.46	
		Hex ₆ HexNAc ₃ NeuAc ₃ + OAc	6_5_0_3*a	C111 H180 N8 O81	2921.0212	2921.0160	-1.77	10.0, 10.37, 10.81, 11.27, 11.46	
		Hex ₆ HexNAc ₅ NeuAc ₃ + 2OAc	6_5_0_3*b	C113 H182 N8 O82	2963.0318	2963.0250	-2.28	10.32, 10.72, 11.02, 11.44, 12.07	
28 Hex ₆ HexNAc ₅ NeuAc ₃ +3OAc 6_5_0_3'c C115 H184 N8 O83 3005.0423 3005.0390 -1.10		$Hex_6HexNAc_5NeuAc_3 + 3OAc$	6_5_0_3*c	C115 H184 N8 O83	3005.0423	3005.0390	-1.10	11.42	
29 Hex ₆ HexNAc ₅ NeuAc ₄ +2OAc 6_5_0_4'b C124 H199 N9 O90 3254.1272 3254.1221 -1.57 11		Hex ₆ HexNAc ₅ NeuAc ₄ + 2OAc	6_5_0_4*b	C124 H199 N9 O90	3254.1272	3254.1221	-1.57	11.29, 11.77	
30 Hex ₆ HexNAc ₅ NeuAc ₄ +3OAc 6_5_0_4'c C126 H201 N9 O91 3296.1378 3296.1379 0.06 11		Hex ₆ HexNAc ₅ NeuAc ₄ + 3OAc	6_5_0_4*c	C126 H201 N9 O91	3296.1378	3296.1379	0.06	11.49, 11.96	

Table 2. Identified *N***-glycans containing** *O***-acetylated sialic acid.** $^{\&}$ novel *N*-glycan; * NeuAc; $^{\#}$ NeuGc; a mono-O-acetylation; b bi-O-acetylation; c tri-O-acetylation.

Collection of serum samples. Sera of male Sprague-Dawley rats weighing 200–250 g and ICR male mice weighing 20–22 g were collected via orbital eye bleeding. Sera from volunteers were obtained from the Guangdong Provincial Chinese Medicine Hospital.

Sample preparation. *N*-glycan release and purification were performed according to the PNGase F protocol provided by New England Biolabs (Ipswich, MA, USA). Briefly, 5μ L of serum were thermally denatured in denaturing buffer containing 0.5% SDS and 40 mM DTT, 10 min prior to digestion by PNGase F (New England Biolabs). Digestion was performed in 50 mM sodium phosphate and 1% NP-40 overnight at 37 °C. Released *N*-glycans were purified by Hypercarb solid-phase extraction (SPE) cartridges, washed with water, then eluted with 40% acetonitrile and 0.05% trifluoroacetic acid (v/v) in water^{19,22,23}. The samples were dried *in vacuo* before MS analysis.




Figure 4. The overall relative abundances of 97 distinct *N*-glycan compositions in rat, mouse, and human sera (a) and relative abundances of each *N*-glycan class. (b) High mannose (HM), undecorated complex/hybrid (C/H), fucosylated complex/hybrid (C/H-Fuc), sialylated complex/hybrid (C/H-Sia), and fucosylated/sialylated complex/hybrid (C/H-Fuc&Sia). For the full name of each *N*-glycan refer to Supplementary Table S1.




Figure 5. Species-specific variations of N-glycans with different numbers of sialic acids, (a) with or without fucose, (b) and with different numbers of antennary, (c) in rat, mouse, and human sera.

Chromatographic separations and MS analysis of the serum *N*-glycome. Samples were analyzed using a microfluidic chip-LC coupled with the Agilent 6550 iFunnel Accurate Mass Quadrupole Time-of-Flight Mass Spectrometer System (Agilent Technologies, Santa Clara, CA, USA) equipped with an auto-sampler (maintained at 5 °C), capillary pumps, nano pumps, and a chip/MS interface. The microfluidic chip consisted of a 40-nL enrichment column and a 43 × 0.075 mm i.d. analytical column, both




Figure 6. Variant relative abundance of *N*-glycans with NeuAc or NeuGc, or with both NeuAc and NeuGc in rat, mouse, and human sera (a) and the retention behavior of O-acetylated N-glycans on the PGC chip (b) Insert table is the relative abundance of O-acetylated *N*-glycans in three species sera.

packed with $5\,\mu m$ graphitized carbon as the stationary phase, with integrated nano-ESI spray tips. For each sample, $1.0\,\mu L$ of sample solution was loaded onto the enrichment column and washed with a solution of 0.1% formic acid (v/v) in water. A rapid N-glycan elution gradient was delivered at $0.5\,\mu L/m$ min using solutions of (A) 65 mM formic acid buffered to pH 3 in water and (B) 100% acetonitrile, at the following proportions and time points: 5-60%, $0-12\,m$ in; Remaining non-glycan compounds were flushed out with 80% B at $0.5\,\mu L/m$ in for 3 min, while the enrichment column was re-equilibrated with 0.1% formic acid at $3\,\mu L/m$ in for $10\,m$ in. The drying gas temperature was set at $225\,^{\circ}C$ with a flow rate of $11\,L/m$ in (filtered nitrogen gas). MS spectra were acquired over a mass range of m/z 500–3000 with an acquisition time of $1.0\,s$ per spectrum in positive ionization mode. The instrument was operated using the target MS/MS mode, with the m/z range from m/z 50 to 3000 with an acquisition time of $1.5\,s$ per spectrum.

Mass correction was enabled using a reference mass of m/z 922.0098 as the internal standard (G1969-85001; Agilent Technologies). The collision energy was set at 5–20 V. The full-width half maximum of the quadrupole mass bandpass used during MS/MS precursor isolation was set to medium ($\sim 4 m/z$).

LC/MS data processing and *N*-glycan identification using accurate mass data. LC/MS raw data were processed using the Molecular Feature Extractor (MFE) algorithm (Version B.06.00; Agilent Technologies). MFE is able to perform chemical relationship testing and chromatographic covariance testing, identify charge carriers (such as sodium) and multimers and group them, and reconstruct spectra by including isotope information. The critical parameters setting in the MFE includes the *N*-glycan model in the isotope model item, three as the maximum charge states, and 20 ppm as the accurate mass

criteria. MS peaks were filtered with a signal-to-noise ratio of 5.0 and parsed into individual ion species. Using the expected isotopic distribution, charge state information and retention time, all ion species associated with a single compound (e.g., the doubly protonated ion, the triply protonated ion, and all associated isotopologues) were totaled, and the neutral monoisotopic mass of the compound was calculated. Using this information, a list of all peaks in the sample was generated with abundances represented by chromatographic peak areas. Computerized algorithms were used to identify N-glycan compositions by accurate mass. By combining these empirical findings with previous research into the N-glycans of mammals, a virtual personal compound database and library (PCDL) was established resulting in over 4000 N-glycan compositions, which contained all biologically plausible rat serum N-glycan compositions with modifications such as O-acetylation, methylation, lactylation, sulfation, phosphorylation, and glycosylation. Deconvoluted experimental masses were compared against theoretical N-glycan masses using a mass error tolerance of 5 ppm. The scoring of the generated formulas was based on three factors: first, the measured mass (or m/z) was compared with the value predicted from the proposed formula; second, the abundance pattern of the measured isotope cluster was compared with values predicted from the proposed formula; third, the m/z spacing between the lowest m/z ion and the A + 1 and A + 2 ions were compared with the values predicted from the proposed formula. These individual factors were computed as match probabilities. Combining the individual match probabilities into an overall score was done as a weighted average rather than as a product. On the basis of known rat serum N-glycosylation patterns, N-glycan compositions containing hexose (Hex), N-acetylhexosamine (HexNAc), deoxyhexose (dHex), N-acetylneuraminic acid (NeuAc), N-glycolylneuraminic acid (NeuGc), and O-acetylation (OAc) were considered.

References

- 1. Lindblad-Toh, K. Genome sequencing: three's company. Nature 428, 475-476 (2004).
- 2. Ohmachi, T., Fukuoka, R., Kimura, Y., Asada, Y. & Ennis, H. L. The characterization of two Dictyostelium discoideum genes encoding ribosomal proteins with sequence similarity to rat L27a and L37a. *Biosci. Biotechnol. Biochem.* **62**, 2008–2015 (1998).
- 3. Abbott, A. Laboratory animals: the Renaissance rat. Nature 428, 464-466 (2004).
- Nakahara, T. et al. Glycoconjugate Data Bank: Structures—an annotated glycan structure database and N-glycan primary structure verification service. Nucleic. Acids Res. 36, D368–371 (2008).
- Alley, W. R., Jr., Mann, B. F. & Novotny, M. V. High-sensitivity analytical approaches for the structural characterization of glycoproteins. Chem. Rev. 113, 2668–2732 (2013).
- 6. Hart, G. W. & Copeland, R. J. Glycomics hits the big time. Cell 143, 672-676 (2010).
- 7. Saldova, R. et al. Increase in Sialylation and Branching in the Mouse Serum N-glycome Correlates with Inflammation and Ovarian Tumour Progression. PLoS One 8, e71159 (2013).
- 8. Adamczyk, B., Tharmalingam, T. & Rudd, P. M. Glycans as cancer biomarkers. Biochim. Biophys. Acta 1820, 1347-1353 (2012).
- Nakagawa, T. et al. Glycomic analyses of glycoproteins in bile and serum during rat hepatocarcinogenesis. J. Proteome. Res. 9, 4888–4896 (2010).
- 10. Desmyter, L. et al. Rating of CCl(4)-induced rat liver fibrosis by blood serum glycomics. J. Gastroenterol. Hepatol. 22, 1148–1154 (2007)
- 11. Fang, M. et al. Serum N-glycome biomarker for monitoring development of DENA-induced hepatocellular carcinoma in rat. Mol. Cancer 9, 215 (2010).
- Bereman, M. S., Williams, T. I. & Muddiman, D. C. Development of a nanoLC LTQ orbitrap mass spectrometric method for profiling glycans derived from plasma from healthy, benign tumor control, and epithelial ovarian cancer patients. *Anal. Chem.* 81, 1130–1136 (2009).
- Bereman, M. S., Young, D. D., Deiters, A. & Muddiman, D. C. Development of a robust and high throughput method for profiling N-linked glycans derived from plasma glycoproteins by NanoLC-FTICR mass spectrometry. J. Proteome. Res. 8, 3764–3770 (2009).
- 14. Houbart, V. et al. Development of a nano-liquid chromatography on chip tandem mass spectrometry method for high-sensitivity hepcidin quantitation. J. Chromatogr. A 1218, 9046–9054 (2011).
- Palmisano, G., Larsen, M. R., Packer, N. H. & Thaysen-Andersen, M. Structural analysis of glycoprotein sialylation part II: LC-MS based detection. RSC Advances 3, 22706–22726 (2013).
- 16. Shah, P. et al. Mass spectrometric analysis of sialylated glycans with use of solid-phase labeling of sialic acids. Anal. Chem. 85, 3606–3613 (2013).
- Morris, H. R. et al. Gender-specific glycosylation of human glycodelin affects its contraceptive activity. J Biol. Chem. 271, 32159–32167 (1996).
- 18. Maley, F., Trimble, R. B., Tarentino, A. L. & Plummer Jr, T. H. Characterization of glycoproteins and their associated oligosaccharides through the use of endoglycosidases. *Anal. Biochem.* **180**, 195–204 (1989).
- 19. Hua, S. et al. Isomer-specific LC/MS and LC/MS/MS profiling of the mouse serum N-glycome revealing a number of novel sialylated N-glycans. Anal. Chem. 85, 4636–4643 (2013).
- Carlier, J. et al. Integrated microfabricated systems including a purification module and an on-chip nano electrospray ionization interface for biological analysis. J. Chromatogr. A 1071, 213–222 (2005).
- 21. Noga, M., Sucharski, F., Suder, P. & Silberring, J. A practical guide to nano-LC troubleshooting. J. Sep. Sci. 30, 2179-2189 (2007).
- 22. Chen, X. & Flynn, G. C. Analysis of *N*-glycans from recombinant immunoglobulin G by on-line reversed-phase high-performance liquid chromatography/mass spectrometry. *Anal. Biochem.* **370**, 147–161 (2007).
- Packer, N. H., Lawson, M. A., Jardine, D. R. & Redmond, J. W. A general approach to desalting oligosaccharides released from glycoproteins. Glycoconj. J. 15, 737–747 (1998).
- 24. Harvey, D. J. et al. Differentiation between isomeric triantennary N-linked glycans by negative ion tandem mass spectrometry and confirmation of glycans containing galactose attached to the bisecting (beta1-4-GlcNAc) residue in N-glycans from IgG. Rapid Commun. Mass Spectrom. 22, 1047–1052 (2008).
- 25. Yuen, C. T., Gee, C. K. & Jones, C. High-performance liquid chromatographic profiling of fluorescent labelled *N*-glycans on glycoproteins. *Biomed. Chromatogr.* **16**, 247–254 (2002).
- Abdul Rahman, S. et al. Filter-Aided N-Glycan Separation (FANGS): A Convenient Sample Preparation Method for Mass Spectrometric N-Glycan Profiling. J. Proteome. Res. 13, 1167–1176 (2014).

- 27. Ruhaak, L. R. et al. Enrichment strategies in glycomics-based lung cancer biomarker development. Proteomics Clin. Appl. 7. 664-676 (2013).
- 28. Sekiya, S., Wada, Y. & Tanaka, K. Derivatization for stabilizing sialic acids in MALDI-MS. Anal. Chem. 77, 4962-4968 (2005).
- 29. Yang, S. & Zhang, H. Glycomic analysis of glycans released from glycoproteins using chemical immobilization and mass spectrometry. Curr. Protoc. Chem. Biol. 6, 191–208 (2014).
- 30. Hua, S. et al. Isomer-specific chromatographic profiling yields highly sensitive and specific potential N-glycan biomarkers for epithelial ovarian cancer. J. Chromatogr. A 1279, 58–67 (2013).
- 31. Sano, K. et al. Survival signals of hepatic stellate cells in liver regeneration are regulated by glycosylation changes in rat vitronectin, especially decreased sialylation. J. Biol. Chem. 285, 17301–17309 (2010).
- 32. Hedlund, M. et al. N-glycolylneuraminic acid deficiency in mice: implications for human biology and evolution. Mol. Cell. Biol. 27, 4340-4346 (2007)
- 33. Raju, T. S., Briggs, J. B., Borge, S. M. & Jones, A. J. Species-specific variation in glycosylation of IgG: evidence for the species-specific sialylation and branch-specific galactosylation and importance for engineering recombinant glycoprotein therapeutics. *Glycobiology* 10, 477–486 (2000).
- 34. Hua, S. et al. Comprehensive native glycan profiling with isomer separation and quantitation for the discovery of cancer biomarkers. Analyst 136, 3663–3671 (2011).
- 35. Zhu, J. et al. A simple integrated system for rapid analysis of sialic-acid-containing N-glycopeptides from human serum. Proteomics 13, 1306–1313 (2013).
- Tangvoranuntakul, P. et al. Human uptake and incorporation of an immunogenic nonhuman dietary sialic acid. Proc. Natl. Acad. Sci. USA 100, 12045–12050 (2003).
- 37. Lin, S. Y. et al. Precise mapping of increased sialylation pattern and the expression of acute phase proteins accompanying murine tumor progression in BALB/c mouse by integrated sera proteomics and glycomics. J. Proteome Res. 7, 3293–3303 (2008).
- 38. Varki, A. N-glycolylneuraminic acid deficiency in humans. Biochimie 83, 615-622 (2001).
- 39. Varki, A. Glycan-based interactions involving vertebrate sialic-acid-recognizing proteins. Nature 446, 1023–1029 (2007).
- 40. Marti, T. et al. The N- and O-linked carbohydrate chains of human, bovine and porcine plasminogen. Species specificity in relation to sialylation and fucosylation patterns. Eur. J. Biochem. 173, 57–63 (1988).
- 41. Yoneda, A. et al. Structures of the N-linked oligosaccharides on porcine plasma vitronectin. Eur. J. Biochem. 218, 797-806 (1993).
- 42. Chaze, T. et al. Alteration of the serum N-glycome of mice locally exposed to high doses of ionizing radiation. Mol. Cell. Proteomics 12, 283-301 (2013).
- 43. Schauer, R. Achievements and challenges of sialic acid research. Glycoconj. J. 17, 485-499 (2000).
- 44. Bones, J. et al. Identification of N-glycans displaying mannose-6-phosphate and their site of attachment on therapeutic enzymes for lysosomal storage disorder treatment. Anal. Chem. 83, 5344–5352 (2011).
- 45. Pabst, M. & Altmann, F. Influence of electrosorption, solvent, temperature, and ion polarity on the performance of LC-ESI-MS using graphitic carbon for acidic oligosaccharides. *Anal. Chem.* **80**, 7534–7542 (2008).
- 46. Takashiba, M., Chiba, Y. & Jigami, Y. Identification of phosphorylation sites in N-linked glycans by matrix-assisted laser desorption/ionization time-of-flight mass spectrometry. Anal. Chem. 78, 5208–5213 (2006).
- 47. Wheeler, S. F. & Harvey, D. J. Extension of the in-gel release method for structural analysis of neutral and sialylated *N*-linked glycans to the analysis of sulfated glycans: application to the glycans from bovine thyroid-stimulating hormone. *Anal. Biochem.* **296.** 92–100 (2001).
- 48. Liu, X. et al. O-acetylation of sialic acids in N-glycans of Atlantic salmon (Salmo salar) serum is altered by handling stress. Proteomics 8, 2849–2857 (2008).
- 49. Nakano, M., Kakehi, K., Tsai, M. H. & Lee, Y. C. Detailed structural features of glycan chains derived from alpha1-acid glycoproteins of several different animals: the presence of hypersialylated, *O*-acetylated sialic acids but not disialyl residues. *Glycobiology* 14, 431–441 (2004).
- 50. Ruhaak, L. R., Miyamoto, S., Kelly, K. & Lebrilla, C. B. N-Glycan Profiling of Dried Blood Spots. Anal. Chem. 84, 396-402 (2011).
- 51. Aldredge, D., An, H. J., Tang, N., Waddell, K. & Lebrilla, C. B. Annotation of a Serum *N*-Glycan Library for Rapid Identification of Structures. *J. Proteome. Res.* 11, 1958–1968 (2012).
- 52. Klein, A. & Roussel, P. O-acetylation of sialic acids. *Biochimie* **80**, 49–57 (1998).
- 53. Varki, A. Diversity in the sialic acids. *Glycobiology* **2,** 25–40 (1992).
- 54. Varki, A. Evolutionary forces shaping the Golgi glycosylation machinery: why cell surface glycans are universal to living cells. Cold Spring Harb Perspect Biol. 3, a005462 (2011).
- 55. Batzoglou, S., Pachter, L., Mesirov, J. P., Berger, B. & Lander, E. S. Human and mouse gene structure: comparative analysis and application to exon prediction. *Genome. Res.* **10**, 950–958 (2000).
- Varki, A. Loss of N-glycolylneuraminic acid in humans: Mechanisms, consequences, and implications for hominid evolution. Am. J. Phys. Anthropol. Suppl 33, 54–69 (2001).
- 57. Schauer, R. Sialic acids: fascinating sugars in higher animals and man. *Zoology (Jena)* **107**, 49–64 (2004).
- 58. Horvat, T., Zoldos, V. & Lauc, G. Evolutional and clinical implications of the epigenetic regulation of protein glycosylation. *Clin. Epigenetics.* 2, 425–432 (2011).
- 59. Horvat, T. et al. Reversibility of membrane N-glycome of HeLa cells upon treatment with epigenetic inhibitors. PLoS One 8, e54672 (2013).
- 60. Horvat, T., Muzinic, A., Barisic, D., Bosnar, M. H. & Zoldos, V. Epigenetic modulation of the HeLa cell membrane N-glycome. *Biochim. Biophys. Acta* 1820, 1412–1419 (2012).
- 61. Lauc, G., Vojta, A. & Zoldoš, V. Epigenetic regulation of glycosylation is the quantum mechanics of biology. *Biochi. Biophysi. Acta* **1840**, 65–70 (2014).
- 62. Zoldos, V. et al. Epigenetic silencing of HNF1A associates with changes in the composition of the human plasma N-glycome. Epigenetics 7, 164–172 (2012).
- Higa, H. H., Rogers, G. N. & Paulson, J. C. Influenza virus hemagglutinins differentiate between receptor determinants bearing N-acetyl-, N-glycollyl-, and N,O-diacetylneuraminic acids. Virology 144, 279–282 (1985).
- 64. Gibbs, R. A. et al. Genome sequence of the Brown Norway rat yields insights into mammalian evolution. *Nature* 428, 493–521 (2004).

Acknowledgments

We are especially grateful to Qiong Meng for technical guidance and assistance. This project was sponsored by the Macao Science and Technology Development Fund (020/2013/A1 to Z.-H. Jiang) and Tertiary Education Services Office, Macau Special Administrative Region.

Author Contributions

Z.H.J. and J.R.W. conceived the study. L.L., J.R.W., D.C.C. and M.J. designed the experiments. W.N.G. and L.F.Y. carried out the experiments. X.Z. and J.L. analyzed the experimental data. Z.H.J., J.R.W. and W.N.G. wrote the manuscript. All authors have read and approved the final manuscript.

Additional Information

Supplementary information accompanies this paper at http://www.nature.com/srep

Competing financial interests: The authors declare no competing financial interests.

How to cite this article: Gao, W.-N. et al. Microfluidic Chip-LC/MS-based Glycomic Analysis Revealed Distinct N-glycan Profile of Rat Serum. Sci. Rep. 5, 12844; doi: 10.1038/srep12844 (2015).

This work is licensed under a Creative Commons Attribution 4.0 International License. The images or other third party material in this article are included in the article's Creative Commons license, unless indicated otherwise in the credit line; if the material is not included under the Creative Commons license, users will need to obtain permission from the license holder to reproduce the material. To view a copy of this license, visit http://creativecommons.org/licenses/by/4.0/