

NEW PRODUCTS AND FOCUS ON IMPRESSION MATERIALS AND SMALL EQUIPMENT/CONSUMABLES

Please send product news information and images to Kate Quinlan at the *BDJ*, Nature Publishing Group, The Macmillan Building, 4-6 Crinan Street, London, N1 9XW. Product news is provided as a service to readers using text and images from the manufacturer, supplier or distributor and does not imply endorsement by the *BDJ*. Normal and prudent research should be exercised before purchase or use of any product mentioned.

BUY A TOOTHBRUSH FOR A CHILD IN NEED

The Dental Wellness Trust has teamed up with Henry Schein, CTS and Whitewash Laboratories to launch the Good Deed brand.

The competitively priced Good Deed brushes allow dental professionals and consumers to do a good deed every time they sell or buy a brush. For each Good Deed brush sold, a brush will be given to a child somewhere who has no access to a toothbrush normally.

The collaboration aim to give 3 million toothbrushes to children over the course of the next year.

The brushes will be available to dental professionals to sell to patients in their practices. Visit www.gooddeed.info for pricing and to find out how to order your Good Deed brushes.

COLOUR CODED DENTAL DISINFECTION

Mistakes in dental disinfection procedures can be extremely dangerous. For this reason, Dürr Dental has colour coded its preparations with a special four-colour system to provide users with greater clarity. It conveniently classifies over 30 products into their area of application: green for surfaces, blue for instruments, yellow for special areas and suction units, and pink for skin and hands.

Dürr Dental concentrate for instrument disinfection ID 213 or drill disinfection ID 220, for example, are coded with the colour blue. Both preparations are bactericidal, tuberculocidal, fungicidal and fully virucidal against enveloped and non-enveloped viruses and meet all requirements of the Robert Koch Institute.

Pink is for hand disinfection and features products such as the rub-in preparation HD 410. Application time for hygienic disinfection is cut in half from 30 to 15 seconds, and surgical hand disinfection requires only 1.5 minutes as opposed to the normal five minutes. It is also moisturising and gentle on skin.

Green is the Dürr System Hygiene colour for disinfection of contaminated surfaces. Wipe disinfectant



FD 322 boasts a 15-second application time and fully virucidal quick disinfectant FD 333 an application time of one minute.

Disinfection and cleaning products in the fourth, yellow-coded speciality area, remove pathogens from suction units. Foam-free liquid concentrate Orotol plus with long-term effect is suited to simultaneous disinfection, cleaning, deodorising and care of suction units and amalgam separators. Find out more at www.duerrdental.de.

EXACTLY THE RIGHT X-RAY DOSE



The launch of Dürr Dental's new X-ray system, VistaIntra, works in perfect partnership with their image plate scanner, VistaScan.

VistaIntra offers the quality you would expect from Dürr Dental. The unit is modern and slim in design and is easy to operate and position. It is pre-programmed with exactly the right X-ray dose when using Dürr Dental image plates and sensors. This ensures it is ready for use and always delivers clear, perfectly exposed images. Moreover, the

patient radiation dose is reduced by 25% compared to conventional AC generators, ensuring greater safety for both patient and practitioner.

It makes sense to use equipment that works in synergy with each other. Dürr VistaIntra and VistaScan do just that, offering you a complete system from a single provider. Optimally matched components ensure you have the best image results.

For more information visit www.duerrdental.de.