## **HITGEN LTD**

www.hitgen.com

## HitGen - OpenDEL

## Making encoded library technology a viable alternative to HTS

HitGen's new OpenDEL encoded library of over 145 million small molecules, all with drug-like properties, can be screened using only 2 mg of target protein, at costs comparable to HTS on much smaller decks. The OpenDEL library is based on 70 different templates and is designed for high diversity within drug-like space.

OpenDEL is just one of the encoded library collections available for screening at HitGen; the other libraries, including the Macrocycle Lead-Finder library, have demonstrated the ability to generate exciting hits even for challenging targets such as PPIs.

## **Contact details**

Jin Li, CEO | F7-10, Building B3, Tianfu Life Science Park, 88 South Keyuan Road, Chengdu 610041, Sichuan, China Tel: +86 288 519 7385 | Email: leadgen.info@hitgen.com



Harnessing the power of Encoded Library Technology for your new targets