PRESERVE AND PROTECT

The cleaning and disinfection of your treatment centre will preserve its life as well as protecting both patients and staff against infection. That is why Belmont has developed a range of care products to ensure their new equipment retains its aesthetics.

The products can be divided into two categories: chair maintenance (B100 & B300) and dental unit waterline protection (B700 & B900).

Chairs need to be wiped down between patients and for such regular use you need a fast-acting upholstery cleaner that will not damage artificial leather; B300 does just that. It has reduced levels of alcohol so that over time it will prevent brittle, cracked surfaces. However, it is bactericidal (even killing TB), yeasticidal and has limited virucidal impact, being effective against non-enveloped viruses, including norovirus.

The perfect adjunct to this is B100, an intensive cleaner designed for occasional use, to remove stains of all kinds as well as discolouration. Despite its rigorous cleaning action, it will not damage your upholstery.

The daily care of your waterlines is another routine hygiene task that must be carried out. B700 is a tried and tested maintenance solution for waterlines. It maintains the water flowing through your unit and will minimise the potential for biofilm

formation. B700 is supplied in a single-dose sachet. B700 is also non-effervescent, which means it can be used instantly.

Before you start using B700 it is essential that you test and 'shock' your system using the B900 kit. This will remove any existing biofilm. Thereafter, you're advised to perform this quarterly to prevent the build-up of future biofilm. The B900 kit contains a box of five shock bottles, tester strips and ten dipstick slides for 'before' and 'after' evaluation.

For your convenience Belmont have also put together a starter kit, called B500, so that you can test and treat your water lines with minimal investment.

Whether you have a new or older Belmont unit, you can contact your dealer to order products from the range and ensure your unit remains looking and functioning optimally.



A CLEANER AND DISINFECTANT WITH A

OPT*IM* 1 is the market-leading wipe in North America, now available in Europe. OPT*IM* 1 wipes:

1-MINUTE CONTACT TIME

- Clean and disinfect, with a 1-minute contact time
- Are gentle and non-toxic and will not irritate the eyes, skin, or respiratory system
- Are compatible with a range of hard, nonporous materials and surfaces
- Have hydrogen peroxide as the active ingredient, which breaks down into water and oxygen, making them kind to the environment
- Have been tested and passed for all major dental chair upholstery for added protection.

OPTIM 1 wipes are effective against

difficult-topathogens, includi COVID-19 (coronavirus).

Available from SciCan, recently acquired by COLTENE.

For information about OPTIM 1 visit: https://www.scican.com/eu/products/cleaners-disinfectants/optim-1/.

www.coltene.com, email info.uk@coltene.com or call 01444 235486.

PROTECT AGAINST INFECTIOUS DISEASES

Don't let infectious diseases threaten the safety of patients or staff. Ensure your dental instruments are thoroughly cleaned and sterilised with reliable, proven solutions from W&H. The new Teon+thermal washer disinfector offers dedicated programmes to ensure instruments are gently and effectively processed. With its user-friendly display and high-capacity chamber, the Teon+ can be combined with the new generation of W&H sterilisers – including the Lara and Lisa type B vacuum steriliser – to further streamline your decontamination workflow.

Lisa has been innovatively designed with patented Eco Dry+ technology, which reduces the cycle time and optimises energy consumption for truly efficient sterilisation procedures.

Explore the decontamination range by visiting the W&H website today.

To find out more visit www.wh.com/en_uk, call 01727 874990 or email office. uk@wh.com.

HAND CARE THAT MAKES A DIFFERENCE

Hand hygiene is incredibly important to follow in order to prevent the spread of disease in your dental practice.

So why not ensure pathogens have no chance of transferring to other surfaces by choosing the Steri-7 Xtra range of disinfectants from Initial Medical?

The Personal Care Range includes a hypoallergenic surgical level hand rub and wash. These products can be used in the surgery but also in washrooms and other spaces where germs may be spread. They kill 99.9999% of pathogens and also have Reactive Barrier Technology – an innovative feature that prevents pathogens recolonising treated surfaces for up to 72

For further information visit www. initial.co.uk/medical or telephone 0870 850 4045.

hours.

