

New products and focus on implants and periodontology

Please send product news information and images to Kate Maynard at the *BDJ*, Nature Publishing Group, The Macmillan Building, 4–6 Crinan Street, London, N1 9XW.

Product news is provided as a service to readers using text and images from the manufacturer, supplier or distributor and does not imply endorsement by the *BDJ*. Normal and prudent research should be exercised before purchase or use of any product mentioned.

A FUN WAY TO KEEP DECAY AT BAY

To help give children the best start possible where good oral hygiene habits are concerned, recommend Spry Toothgel with Xylitol for the ideal way to keep cavities away.

Xylitol is a naturally occurring sweetener that is proven to be a powerful tool in the fight against tooth decay. Thanks to its unique chemical composition, xylitol is unfavourable to the metabolism of plaque bacteria, reducing their

ability to proliferate and adhere.

Available in a handy 60 ml tube, Spry Toothgel is a great way to introduce children to oral hygiene and comes in a number of child-friendly flavours to make teeth cleaning even more fun.

Based on all-natural ingredients, Spry Toothgel is completely safe for infants – it is safe to swallow and fluoride-free. For the very young, Spry Toothgel can be used to gently massage into gums and teeth or can be used in conjunction with the Spry Baby Banana Brush as the ideal introduction to regular brushing.

For further information contact Anyone 4 Tea Ltd on 01428 652131 or visit www.anyone4tea.com.



TABLET OPERATED DRILL MOTOR

Nobel Biocare has launched OsseoCare Pro – the first drill motor to be operated by an iPad.

The free OsseoCare Pro app enables clinicians to plan the treatment sequence prior to surgery, is customisable and comes packed with safety features such as guided drilling protocols. With multiple-user log-ins, treatment data can be shared with colleagues quickly and conveniently.

The OsseoCare Pro LED handpiece can effortlessly access areas of limited space and has one of the lightest

and most powerful motors on the market. The device also has ample and stable lighting during surgery as well as optimal cooling.

Speed, torque, irrigation flow and light intensity can all be controlled and modified through the app, which also offers a built-in recording and exporting function.

New app features and functions will be added and updated regularly.

For more information contact Nobel Biocare on 020 8756 3300 or visit www.nobelbiocare.com.

GET CREDIT FOR YOUR BLEACHING LAMPS

In a bid to help dental practices meet the new tooth whitening directive, which came into force at the beginning of this month, Enlighten is offering to buy back your bleaching lights.

European law dictates that UK dentists can no longer supply tooth whitening products which release or contain more than 6% hydrogen peroxide (18% carbamide peroxide). As a result there may be numerous practices that still have bleaching lights on the premises, which cannot be legally used.

You can send your floor-standing bleaching lamps to Enlighten in exchange for £200 credit towards Enlighten Whitening Kits. Alternatively, if you prefer, Enlighten will make a £75 donation to the Bridge2Aid charity.

To take advantage of this offer, pack your lamp in a suitable box, with suitable protective packaging, and Enlighten will pick it up from your practice for free (one lamp per practice and the lamp must be in good working order on delivery; practices must be in the UK and GDC-registered).

Enlighten reserves the right to refuse to take any individual lamp and may withdraw the promotion without prior notice.

Contact Enlighten Smiles on 020 7424 3270 or visit www.enlightensmiles.com/buy-back.

A REVOLUTIONARY PHYSIODISPENSER

The iChiropro by Bien-Air is a revolutionary physiodispenser controlled from an iPad. Simple to use, the intuitive, ergonomic application ensures significant time savings during operations.

The multi-user iChiropro incorporates the complete sequences of the leading implant manufacturers. It also offers the option of saving and exporting operation data attached to the patient file built into the control interface.

Telephone 01293 550200 or email ba-uk@bienair.com for more information.